
Molecular Imaging Documentation

Release 0.1

Molecular Imagers

April 28, 2016

1	General Protocols	1
1.1	Common Buffers	1
2	Cell Culture Media	5
2.1	Dulbecco's Modified Eagle Media	5
2.2	Freezing Medium	5
3	Indices and tables	7

General Protocols

0.1

1.1 Common Buffers

1.1.1 Cell Culture

1x PBS

To 800ml ddH₂O add:

8g NaCl 0.2g KCl 1.44g Na₂HPO₄ (sodium phosphate dibasic) 0.24g KH₂PO₄ (monopotassium phosphate)

pH to 7.4, and bring up to 1L with ddH₂O. Autoclave to sterilize.

10x PBS

To 800ml ddH₂O add:

- 80g NaCl
- 2.0g KCl
- 14.4g Na₂HPO₄ (sodium phosphate dibasic)
- 2.4g KH₂PO₄ (monopotassium phosphate)

pH to 7.4, and bring up to 1L with ddH₂O. Autoclave to sterilize.

50x TAE

For each litre of solution:

242g Tris Base 57.1ml Glacial Acetic Acid 100ml 0.5M EDTA

Mix Tris with stir bar to dissolve in about 600ml of ddH₂O. Add the EDTA and Acetic Acid. Bring final volume to 1L with ddH₂O. Store at room temperature.

Note: Final (1x) working concentration:

0.04M Tris-Acetate 0.001M EDTA

1.1.2 Western Blotting

Lysis Buffer

Reagent	Stock Solution	Volume
ddH ₂ O	–	7.68 ml
Tris pH 8.0	20 mM	200 ul
NaCl	0.15 M	750 ul
EDTA	2 mM	40 ul
NP40	–	100 ul
Glycerol	–	1000 ul
Na ₃ VO ₄	1 mM	100 ul

The above recipe is for 10ml. This buffer can be made in bulk and used as a wash buffer, or allotted and frozen for future use. Immediately before use add protease inhibitors at the manufacturers recommended concentration, and phosphatase inhibitors if required.

Laemmli's Buffer, 4x

2.4 ml 1 M Tris pH 6.8 (Same as upper gel buffer) 0.8 g SDS stock 4 ml 100% glycerol 0.01% bromophenol blue. Final Concentration is .02% 2.8 ml ddH₂O

Before use add 1/10th volume of β -mercaptoethanol

Laemmli's Buffer, 6x

1.2g SDS (sodium dodecyl sulfate) 0.01% bromophenol blue 4.7ml glycerol 1.2ml Tris 0.5M pH6.8 2.1ml ddH₂O

Before use add 1/8th volume of β -mercaptoethanol

4x Lower Gel Buffer

182g Tris-Base (1.5M) 1L dH₂O

pH to 8.8; do not overshoot. Use glass pipettes to pH.

4x Upper Gel Buffer

30.28g Tris-HCl (0.5M) 500ml dH₂O

pH to 6.8, do not overshoot. Use glass pipettes.

10x Running Buffer

30.4g Tris-HCl 144.2g Glycine 10g SDS

Dissolve in 1L dH₂O. Dilute 1:10 in ddH₂O for use as running buffer

10X Transfer Buffer (Towbin Buffer)

30.3g Tris-Base 144.15g Glycine 100ml 10% SDS (0.1%)

Bring upto 1L in ddH₂O.

1X Transfer Buffer (Towbin Buffer)

100ml 10X buffer 200ml methanol 700ml ddH₂O

10x TBS

302.5g Tris-Base 400g NaCl 18g KCl

Dilute to 5L in dH₂O, pH to 7.5 with saturated HCl.

Wash Buffer (TBST)

to 1X TBS ADD:

Standard 0.1% Tween-20

Strong 0.1% Tween-20 + 0.1% NP-40

Extra-Strong 0.1% Tween-20 + 1.5% NP-40

Note: Standard buffer is for most application. Strong/extra strong are only for antibodies with a high degree of non-specificity.

DNA Buffers

6X DNA Loading Buffer

6.7 ml of ddH₂O 10 mg of bromophenol blue 10 mg of xylene cyanol FF (optional) 3.3ml glycerol

Notes:

Bromophenol Blue runs at ~300bp (should be added to avoid over-running of gels) Xylene cyanol FF runs at ~4000bp (optional dye) 60 mM EDTA can be added to inhibit DNA-altering enzymes (600ul of 1M EDTA; reduce water appropriately)

Cell Culture Media

2.1 Dulbecco's Modified Eagle Media

DMEM	500 ml	
Glutamine 100x	5 ml	2 mM
FBS	50 ml	10 %
Pen-Strep 100x	5 ml	5 ug/l ??

2.2 Freezing Medium

Complete medium	950 µl	90 %
DMSO	50 µl	10 %

Indices and tables

- `genindex`
- `search`