

---

# **pynetdicom3 Documentation**

*Release 1.0*

**scaramallion**

**Apr 28, 2018**



---

## Contents:

---

<b>1</b>	<b>ApplicationEntity</b>	<b>1</b>
1.1	Initialisation . . . . .	1
1.2	Service Class User . . . . .	1
1.3	Service Class Provider . . . . .	1
1.4	Association . . . . .	2
<b>2</b>	<b>Introduction to Association</b>	<b>3</b>
<b>3</b>	<b>Examples</b>	<b>5</b>
3.1	Service Class User Examples . . . . .	5
3.2	SCP examples . . . . .	7
<b>4</b>	<b>Applications</b>	<b>9</b>
4.1	echoscu . . . . .	9
4.2	echoscp . . . . .	11
<b>5</b>	<b>Package documentation</b>	<b>13</b>
5.1	ApplicationEntity . . . . .	13
5.2	Association . . . . .	13
5.3	SOP Class . . . . .	13
5.4	ACSE . . . . .	13
5.5	DIMSE . . . . .	13
5.6	DUL . . . . .	13
<b>6</b>	<b>Indices and tables</b>	<b>15</b>



ApplicationEntity (or AE) is the main class for constructing a DICOM Application Entity.

### 1.1 Initialisation

```
from pydicom.uid import ImplicitVRLittleEndian
from pynetdicom3 import AE, VerificationSOPClass

ae = AE(ae_title='PYNETDICOM',
        port=11112,
        scu_sop_class=[VerificationSOPClass],
        scp_sop_class=[VerificationSOPClass],
        transfer_syntax=[ImplicitVRLittleEndian])
```

### 1.2 Service Class User

#### 1.2.1 Maximum PDU Size

Each association request can specify its own maximum PDU receive size. A value of None indicates that there is no maximum size limit. >>> ae = AE(scu\_sop\_class=[VerificationSOPClass]) >>> assoc1 = ae.associate('127.0.0.1', 11112, max\_pdu=16382) >>> assoc2 = ae.associate('127.0.0.1', 11112, max\_pdu=None)

### 1.3 Service Class Provider

The following parameters should be set prior to calling *start()*.

### 1.3.1 AE Title Matching

The called and/or calling AE title can be set to be required to match against an expected value, with the association being rejected if they fail. A value of '' means that no matching will be performed. >>> ae = AE(port=11112, scp\_sop\_class=[VerificationSOPClass]) >>> ae.require\_calling\_aet = 'CALLING\_AET' >>> ae.require\_called\_aet = 'CALLED\_AET' >>> ae.start()

### 1.3.2 Maximum Number of Associations

The maximum number of simultaneous associations that the AE will support. Any additional association requests be rejected with a reason of 'Local limit exceeded'. >>> ae = AE(port=11112, scp\_sop\_class=[VerificationSOPClass]) >>> ae.maximum\_associations = 3 >>> ae.start()

### 1.3.3 Timeouts

Timeouts for ACSE, DIMSE and network messages >>> ae = AE(scu\_sop\_class=[VerificationSOPClass]) >>> ae.dimse\_timeout = 30 >>> ae.acse\_timeout = 60 >>> ae.network\_timeout = 60 >>> ae.start()

### 1.3.4 Maximum PDU Size

A value of None indicates that there is no maximum PDU size limit. >>> ae = AE(port=11112, scp\_sop\_class=[VerificationSOPClass]) >>> ae.maximum\_pdu\_size = 16382 >>> ae.start()

## 1.4 Association

## CHAPTER 2

---

### Introduction to Association

---



## 3.1 Service Class User Examples

### 3.1.1 Verification SCU

The most basic way of verifying and troubleshooting a connection with a peer Application Entity (AE) is to send a DICOM C-ECHO, which utilises the *Verification SOP Class*

```
from pynetdicom3 import AE

# The Verification SOP Class has a UID of 1.2.840.10008.1.1
ae = AE(scu_sop_class=['1.2.840.10008.1.1'])

# Try and associate with the peer AE
# Returns the Association thread
print('Requesting Association with the peer')
assoc = ae.associate(addr, port)

if assoc.is_established:
    print('Association accepted by the peer')
    # Send a DIMSE C-ECHO request to the peer
    # status is a pydicom Dataset object with (at a minimum) a
    # (0000, 0900) Status element
    status = assoc.send_c_echo()

    # Output the response from the peer
    if status:
        print('C-ECHO Response: 0x{0:04x}'.format(status.Status))

    # Release the association
    assoc.release()
elif assoc.is_rejected:
    print('Association was rejected by the peer')
```

```
elif assoc.is_aborted:
    print('Received an A-ABORT from the peer during Association')
```

### 3.1.2 Storage SCU

A common use of a DICOM SCU is to use a DICOM C-STORE message to send DICOM datasets to peer AEs.

```
from pydicom import read_file
from pynetdicom3 import AE
from pynetdicom3 import StorageSOPClassList

# StorageSOPClassList contains all the Standard SOP Classes supported
# by the Storage Service Class (see PS3.4 Annex B.5)
ae = AE(scu_sop_class=StorageSOPClassList)

# Try and associate with the peer AE
# Returns the Association thread
print('Requesting Association with the peer')
assoc = ae.associate(addr, port)

if assoc.is_established:
    print('Association accepted by the peer')

    # Read the DICOM dataset from file 'dcmfile'
    dataset = read_file('dcmfile')

    # Send a DIMSE C-STORE request to the peer
    status = assoc.send_c_echo(dataset)
    print('C-STORE status: %s' %status)

    # Release the association
    assoc.release()
```

### 3.1.3 Query/Retrieve - Find SCU

Query the peer AE to see if it contains any Instances with attributes matching those specified by the user-created *dataset*.

```
from pydicom.dataset import Dataset

from pynetdicom3 import AE
from pynetdicom3 import QueryRetrieveSOPClassList

# QueryRetrieveSOPClassList contains the SOP Classes supported
# by the Query/Retrieve Service Class (see PS3.4 Annex C.6)
ae = AE(scu_sop_class=QueryRetrieveSOPClassList)

# Try and associate with the peer AE
# Returns the Association thread
print('Requesting Association with the peer')
assoc = ae.associate(addr, port)

if assoc.is_established:
    print('Association accepted by the peer')
```

```
# Create a new DICOM dataset with the attributes to match against
# In this case match any patient's name at the PATIENT query
# level. See PS3.4 Annex C.6 for the complete list of possible
# attributes and query levels.
dataset = Dataset()
dataset.PatientName = '*'
dataset.QueryRetrieveLevel = "PATIENT"

# Send a DIMSE C-FIND request to the peer
# query_model is the Query/Retrieve Information Model to use
# and is one of 'W', 'P', 'S', 'O'
# 'W' - Modality Worklist (1.2.840.10008.5.1.4.31)
# 'P' - Patient Root (1.2.840.10008.5.1.4.1.2.1.1)
# 'S' - Study Root (1.2.840.10008.5.1.4.1.2.2.1)
# 'O' - Patient/Study Only (1.2.840.10008.5.1.4.1.2.3.1)
responses = assoc.send_c_find(dataset, query_model='P')

for (status, dataset) in responses:
    # While status is pending we should get the matching datasets
    if status == 'Pending':
        print(dataset)
    elif status == 'Success':
        print('C-FIND finished, releasing the association')
    elif status == 'Cancel':
        print('C-FIND cancelled, releasing the association')
    elif status == 'Failure':
        print('C-FIND failed, releasing the association')

# Release the association
assoc.release()
```

## 3.2 SCP examples



## 4.1 echoscu

```
echoscu [options] peer port
```

### 4.1.1 Description

The `echoscu` application implements a Service Class User (SCU) for the *Verification SOP Class* (UID 1.2.840.10008.1.1)<sup>1</sup>. It establishes an Association with a peer Application Entity (AE) which it then sends a DICOM C-ECHO-RQ message<sup>2</sup> and waits for a response. The application can be used to verify basic DICOM connectivity.

The following simple example shows what happens when it is successfully run on an SCP that supports the *Verification SOP Class*:

```
user@host: echoscu 192.168.2.1 11112
user@host:
```

When attempting to send a C-ECHO to an SCP that doesn't support the *Verification SOP Class*:

```
user@host: echoscu 192.168.2.1 11112
E: No Acceptable Presentation Contexts
user@host:
```

When attempting to associate with a non-DICOM peer

```
user@host: echoscu 192.168.2.1 11112
E: Association Request Failed: Failed to establish association
E: Peer aborted Association (or never connected)
E: TCP Initialisation Error: (Connection refused)
user@host:
```

<sup>1</sup> See DICOM Standard 2015b PS3.6 Table A-1

<sup>2</sup> See DICOM Standard 2015b PS3.7 Sections 9.1.5 and 9.3.5

Using the `--propose-pc [n]` option, the echoscu application can also propose  $n$  Presentation Contexts<sup>3</sup> (all with an Abstract Syntax of *Verification SOP Class*) in order to provide debugging assistance for association negotiation

## 4.1.2 Options

### Logging

- `-q --quiet` quiet mode, prints no warnings or errors
- `-v --verbose` verbose mode, prints processing details
- `-d --debug` debug mode, prints debugging information
- `-ll --log-level [l]evel (str)` One of ['critical', 'error', 'warning', 'info', 'debug'], prints logging messages with corresponding level or higher
- `-lc --log-config [f]ilename (str)` use python logging config<sup>4</sup> file `f` for the logger

### Application Entity Titles

- `-aet --calling-aet [a]etitle (str)` set the local AE title (default: ECHOSCU)
- `-aec --called-aet [a]etitle (str)` set the called AE title for the peer AE (default: ANY-SCP)

### Association Negotiation Debugging

- `-pts --propose-ts [n]umber (int)` propose  $n$  transfer syntaxes (1-3)

### Miscellaneous DICOM

- `-to --timeout [s]econds (int)` timeout for connection requests (default: unlimited)
- `-ta --acse-timeout [s]econds (int)` timeout for ACSE messages (default: 30)
- `-td --dimse-timeout [s]econdsr (int)` timeout for DIMSE messages (default: unlimited)
- `-pdu --max-pdu [n]umber of bytes (int)` set maximum receive PDU bytes to  $n$  bytes (default: 16384)
- `--repeat [n]umber (int)` repeat echo  $n$  times
- `--abort` abort association instead of releasing it

## 4.1.3 DICOM Conformance

The echoscu application supports the following SOP Class as an SCU:

Verification SOP Class	1.2.840.10008.1.1
------------------------	-------------------

Unless the `--propose-ts` option is used, the echoscu application will only propose the *Little Endian Implicit VR Transfer Syntax* (UID 1.2.840.10008.1.2). The supported Transfer Syntaxes<sup>5</sup> are:

---

<sup>3</sup> See DICOM Standard 2015b PS3.8 Sections 7.1.1.13 and 9.3.2.2

<sup>4</sup> See the Python documentation

<sup>5</sup> See DICOM Standard 2015b PS3.5 Section 10 and Annex A

Little Endian Implicit VR	1.2.840.10008.1.2
Little Endian Explicit VR	1.2.840.10008.1.2.1
Big Endian Explicit VR	1.2.840.10008.1.2.2

## 4.2 echoscp

```
echoscp [options] port
```

### 4.2.1 Description

The `echoscp` application implements a Service Class Provider (SCP) for the *Verification SOP Class* (UID 1.2.840.10008.1.1)<sup>1</sup>. It establishes an Association with peer Application Entities (AEs) and receives DICOM C-ECHO-RQ<sup>2</sup> message to which it responds with a DICOM C-ECHO-RSP message. The application can be used to verify basic DICOM connectivity.

The following example shows what happens when it is started and receives a C-ECHO from a peer:

```
user@host: echoscp 11112
```

More information is available when a connection is received while running with the `-v` option:

```
user@host: echoscp 11112 -v
I: Association Received
I: Association Acknowledged
I: Received Echo Request (MsgID 1)
I: Association Released
```

When a peer AE attempts to send non C-ECHO message:

```
user@host: echoscu 192.168.2.1 11112 -v
I: Association Received
I: Association Acknowledged
I: Association Aborted
```

Much more information is available when a connection is received while running with the `-d` option:

```
user@host: echoscp 11112 -d
D: $echosco.py v0.2.0 2016-03-15 $
D:
D: Starting DICOM UL service "Thread-1"
D: PDU Type: Associate Request, PDU Length: 215 + 6 bytes PDU header
D: 01 00 00 00 00 d1 00 01 00 00 41 4e 59 2d 53 43
...
D: Request Parameters:
D: ===== BEGIN A-ASSOCIATE-RQ =====
D: Their Implementation Class UID: 1.2.826.0.1.3680043.9.381.0.9.0
...
I: Received Echo Request (MsgID 1)
...
I: Association Released
D: DICOM UL service "Thread-1" stopped
```

<sup>1</sup> See DICOM Standard 2015b PS3.6 Table A-1

<sup>2</sup> See DICOM Standard 2015b PS3.7 Sections 9.1.5 and 9.3.5

## 4.2.2 Options

### Logging

- `-q --quiet` quiet mode, prints no warnings or errors
- `-v --verbose` verbose mode, prints processing details
- `-d --debug` debug mode, prints debugging information
- `-ll --log-level [l]evel (str)` One of ['critical', 'error', 'warning', 'info', 'debug'], prints logging messages with corresponding level l or higher
- `-lc --log-config [f]ilename (str)` use python logging config<sup>3</sup> file f for the logger

### Application Entity Titles

- `-aet --aetitle [a]etitle (str)` set my AE title (default: ECHOSCP)

### Miscellaneous DICOM

- `-to --timeout [s]econds (int)` timeout for connection requests (default: unlimited)
- `-ta --acse-timeout [s]econds (int)` timeout for ACSE messages (default: 30)
- `-td --dimse-timeout [s]econdsr (int)` timeout for DIMSE messages (default: unlimited)
- `-pdu --max-pdu [n]umber of bytes (int)` set maximum receive PDU bytes to n bytes (default: 16384)

### Preferred Transfer Syntaxes

- `-x= --prefer-uncompr` prefer explicit VR local byte order (default)
- `-xe --prefer-little` prefer explicit VR little endian transfer syntax
- `-xb --prefer-big` prefer explicit VR big endian transfer syntax
- `-xi --implicit` accept implicit VR little endian transfer syntax only

## 4.2.3 DICOM Conformance

The echoscp application supports the following SOP Class as an SCP:

Verification SOP Class	1.2.840.10008.1.1
------------------------	-------------------

The supported Transfer Syntaxes<sup>4</sup> are:

Little Endian Implicit VR	1.2.840.10008.1.2
Little Endian Explicit VR	1.2.840.10008.1.2.1
Big Endian Explicit VR	1.2.840.10008.1.2.2

---

<sup>3</sup> See DICOM Standard 2015b PS3.8 Sections 7.1.1.13 and 9.3.2.2

<sup>4</sup> See the Python documentation

## **5.1 ApplicationEntity**

## **5.2 Association**

## **5.3 SOP Class**

## **5.4 ACSE**

## **5.5 DIMSE**

### **5.5.1 DIMSEServiceProvider**

### **5.5.2 DIMSE primitives**

### **5.5.3 DIMSEMessage**

## **5.6 DUL**

### **5.6.1 DULServiceProvider**

### **5.6.2 PDU**

### **5.6.3 PDU primitives**



## CHAPTER 6

---

### Indices and tables

---

- `genindex`
- `modindex`
- `search`