

---

# xrdtools Documentation

*Release 0.1.1*

**Benedikt Ziegler**

**Apr 12, 2018**



---

## Contents

---

<b>1</b>	<b>Table of Content</b>	<b>3</b>
1.1	Getting xrddtools . . . . .	3
1.2	Quick start tutorial . . . . .	3
1.3	Command line tool . . . . .	4
1.4	xrddtools package . . . . .	5
<b>2</b>	<b>Indices and tables</b>	<b>7</b>
	<b>Python Module Index</b>	<b>9</b>



The **xrdtools** is a package intended to load *.xrdml* files created by Panalytical XPert XRD machines. No full feature support is guaranteed at the current state.



# CHAPTER 1

---

## Table of Content

---

### 1.1 Getting xrdtools

The simplest way to install the packages is via pip:

```
$ pip install xrdtools
```

The most reasoned version can be cloned from git by running following commands in your command line:

```
$ git clone https://github.com/paruch-group/xrdtools.git
$ cd xrdtools
$ python setup.py install
```

In case you do not have git installed on your system, go to

<https://github.com/paruch-group/xrdtools.git>

and download the compressed package. Unpack the file in a directory of your choosing and run the following commands in the our terminal:

```
$ cd /path/to/the/extracted/zip/file
$ python setup.py install
```

### 1.2 Quick start tutorial

The main functionality of the **xrdtools** package is to read \*.xrdml files. This can be easily achieved by running the following code, e.g. in a ipython prompt:

```
import xrdtools
data = xrdtools.read_xrdml('foo.xrdml')
```

The data returned from `xrdtools.read_xrdml()` is stored in a `dict`. In case of a simple line scan (e.g. *2theta-omega* scan) we can get the *xy* data as simple as:

```
x = data['x']
y = data['data']
```

And plot it for example with `matplotlib.pyplot`:

```
from matplotlib import pyplot as plt

plt.plot(x, y)
plt.show()
```

## 1.3 Command line tool

Together with the xrdtools package a command line tool `xrdml` is installed. It allows to extract the recorded data from `*.xrdml` files into text files or the command prompt.

Export data into a text file:

```
$ xrdml my_xrdml_file.xrdml #another_file.xrdml ...
```

In case `my_xrdml_file.xrdml` is a simple *2Theta-Omega* scan, this will create a text file with two columns, one for the *2Theta* angles and one for the *Intensity*.

Export data to the prompt:

```
$ xrdml my_xrdml_file.xrdml -o stdout

# 2Theta-Omega  Intensity
1.5000000000000000e+01      9.000000000000000222e-01
1.501869158878504606e+01    5.99999999999999778e-01
...
```

The type of delimiter can be changed with the `--delimiter` keyword argument:

```
$ xrdml my_xrdml_file.xrdml -o stdout --delimiter=','
# 2Theta-Omega, Intensity
1.5000000000000000e+01, 9.000000000000000222e-01
1.501869158878504606e+01, 5.99999999999999778e-01
...
```

The output format can be changed with the `--fmt` keyword argument:

```
$ xrdml my_xrdml_file.xrdml -o stdout --fmt='%.2f'

# 2Theta-Omega  Intensity
15.00  0.90
15.02  0.60
...
```

## 1.4 xrddtools package

### 1.4.1 xrddtools module

`xrddtools.read_xrdml(filename)`

Load a Panalytical XRDML file.

**Parameters** `filename` (`str`) – The filename of the xrdml file to be loaded.

**Returns** A dictionary with all relevant data of the measurement.

**Return type** `dict`

### 1.4.2 xrddtools.io module

`xrddtools.io.read_xrdml(filename)`

Load a Panalytical XRDML file.

**Parameters** `filename` (`str`) – The filename of the xrdml file to be loaded.

**Returns** A dictionary with all relevant data of the measurement.

**Return type** `dict`

`xrddtools.io.validate_xrdml_schema(filename)`

Validate the xml schema of a given file.

**Parameters** `filename` (`str`) – The Filename of the .xrdml file to test.

**Returns** Returns the version number as float or None if the file was not matching any provided xml schema.

**Return type** `float` or `None`

### 1.4.3 xrddtools.utils module

`xrddtools.utils.angle2qvector(tt, om, lam=1.54)`

Convert angles to q vector.

Calculate the q-vector from the 2theta `tt` and omega `om` angle and the x-ray wavelength lambda `lam`.

#### Parameters

- `tt` (`array-like`) – Array containing the 2Theta values.
- `om` (`array-like`) – Array containing the Omega values.
- `lam` (`float`) – The wavelength lambda in Angstrom [Default: 1.54].

#### Returns

- `kpar` (`ndarray`)
- `kperp` (`ndarray`)

`xrddtools.utils.angles(hkl, lam=1.54, lattice_param=(3.905, 3.905, 3.905))`

Compute the angle for a given hkl position.

Compute the 2Theta, Omega and Delta angle for a given hkl point, wavelength lambda and unit cell lattice parameters.

#### Parameters

- **hkl** (*dict*) – A dictionary containing the hkl values.
- **lam** (*float*) – The wavelength lambda in Angstrom. Defaults to 1.54.
- **lattice\_param** (*tuple*) – A tuple of three floats for the lattice parameter.

**Returns**

- **tt** (*ndarray*)
- **omega** (*ndarray*)
- **delta** (*ndarray*)

`xrdtools.utils.get_qmap(data, omega_offset=0)`

Function to calculate kpar, kperp.

**Parameters**

- **data** (*dict*) – A xrdml data dictionary.
- **omega\_offset** (*float*) – Offset for the omega angle.

**Returns**

- **kpar** (*ndarray*)
- **kperp** (*ndarray*)

`xrdtools.utils.q2hkl_map(x, y, lattice_params=(3.905, 3.905, 3.905), hkl=None)`

Compute the hk coordinates for a given q vector.

**Parameters**

- **x** (*ndarray*) –
- **y** (*ndarray*) –
- **lattice\_params** (*tuple*) – A tuple of three floats for the lattice parameter.
- **hkl** (*dict*) – A dictionary containing the hkl values. Defaults to 001 if not given [Default: None].

**Returns**

- **x** (*ndarray*)
- **y** (*ndarray*)

## 1.4.4 Subpackages

### xrdtools.tools package

#### xrdtools.tools.clt module

`xrdtools.tools.clt.xrdml()`

Command line tool to export measurement data from xrdml files.

- o, --output** [str] Choices: ‘stdout’, ‘txt’ [default: ‘txt’]
- delimiter** [str] Default: ‘ ’
- fmt** [str] Default: ‘%.18e’

# CHAPTER 2

---

## Indices and tables

---

- genindex
- modindex
- search



---

## Python Module Index

---

### X

`xrdtools`, 5  
`xrdtools.io`, 5  
`xrdtools.tools.clt`, 6  
`xrdtools.utils`, 5



---

## Index

---

### A

`angle2qvector()` (in module `xrdtools.utils`), 5  
`angles()` (in module `xrdtools.utils`), 5

### G

`get_qmap()` (in module `xrdtools.utils`), 6

### Q

`q2hkl_map()` (in module `xrdtools.utils`), 6

### R

`read_xrdml()` (in module `xrdtools`), 5  
`read_xrdml()` (in module `xrdtools.io`), 5

### V

`validate_xrdml_schema()` (in module `xrdtools.io`), 5

### X

`xrdml()` (in module `xrdtools.tools.clt`), 6  
`xrdtools` (module), 5  
`xrdtools.io` (module), 5  
`xrdtools.tools.clt` (module), 6  
`xrdtools.utils` (module), 5