
BITalino*islDocumentation*

Release 0.0.3

Fernando Suárez Jiménez

Nov 03, 2018

Contents:

1	Introduction	1
1.1	Getting started	1
2	API documentation	3

CHAPTER 1

Introduction

A python module called *bitalino-lsl* to stream BITalino data through the Lab Streaming Layer (LSL). This module gets data from the BITalino device through the *bitalino python api* and uses the *Lab Stream Layer* to stream the data.

The module should work with python versions ≥ 2.7 although it has only been tested for:

- Python 2.7.15
- Python 3.6.5

1.1 Getting started

You can install the module with pip:

```
$ pip install bitalino_lsl
```

This examples streams BITalino data through the LSL and reads it from the stream for 5 seconds:

```
import bitalino_lsl
from pylsl import StreamInlet, resolve_stream
import time

# MAC address of the BITalino device
MAC_ADDRESS_BITALINO_DEVICE = "20:17:11:20:51:60"

# List with channels of the BITalino device to be streamed to the LSL
# This channels can be specified as a list or as a dictionary with their
# position in the 10-20 system. BITalino uses bipolar electrodes so the
# position will be defined by two points
# CHANNELS = {0: 'Fp1-Fp2', 1: 'P3-T5'}
CHANNELS = [0,1]

# Connect with the BITalino device
device = bitalino_lsl.BitalinoLSL(MAC_ADDRESS_BITALINO_DEVICE)
```

(continues on next page)

(continued from previous page)

```
# Create the Stream with the channels information
device.create_lsl_EEG(CHANNELS)

# Start the stream getting data from the BITalino device
device.start()

# Get the Stream to read the data from
inlet = StreamInlet(resolve_stream('type', 'EEG')[0])

# Read the BITalino data for 5 seconds
t_end = time.time() + 5
while time.time() < t_end:
    sample, timestamp = inlet.pull_sample()
    print(sample)

# Stop the device
device.stop()

# Close the connection with the BITalino device
device.close()
```

CHAPTER 2

API documentation
