

Diffcalc User Guide

Release 2.1

Diamond Light Source

Dec 03, 2019

Contents

1	Introduction	3
2	Overview 2.1 Theory	5 5
3	Getting Help	7
4	Diffcalc's Scannables	9
5	Crystal orientation5.1Start a new UB calculation5.2Load a UB calculation5.3Generate a U matrix from two reflections5.4Generate a U matrix from one reflection5.5Edit reflection list5.6Generate a U matrix from two lattice directions5.7Calculate a UB matrix5.8Calculate a U matrix from crystal mismount5.9Manually specify U matrix5.10Refining UB matrix with one reflection5.12Set the reference vector	11 11 12 13 14 14 14 15 15 15 15 15 16 17
	5.13 Set the surface normal vector	17
6	Motion6.1Constraining solutions for moving in hkl space6.2Example constraint modes6.3Changing constrained values6.4Configuring limits and cuts6.5Moving in hkl space	19 19 20 21 21 22
7	Scanning in hkl space 7.1 Fixed hkl scans 7.2 Scanning hkl 7.3 Multidimension scans	25 25 25 26
8	Commands	27

	8.1 8.2	Motion commands	27 29
9	Refe	rences	31
Bi	bliogra	aphy	33

Diffcalc: A diffraction condition calculator for diffractometer control

Author Rob Walton

Contact rob.walton (at) diamond (dot) ac (dot) uk

Website https://github.com/DiamondLightSource/diffcalc

See also the quickstart guide at github

Introduction

This manual assumes that you are running Diffcalc within OpenGDA or have started it using IPython. It assumes that Diffcalc has been configured for the six circle diffractometer pictured here:

in the second second

Fig. 1: 4s + 2d six-circle diffractometer, from H.You (1999)

Your Diffcalc configuration may have been customised for the geometry of your diffractometer and possibly the types of experiment you perform. For example, a five-circle diffractometer might be missing the nu circle above.

The laboratory frame is shown above. With all settings at zero as shown the crystal cartesian frame aligns with the laboratory frame. Therefore a cubic crystal mounted squarely in a way that the U matrix (defined below) is unitary will have hllallx, kllblly & lllcllz, crystal and reciprocal-lattice coordinate frames are defined with respect to the beam and to gravity to be (for a cubic crystal):

Overview

The following assumes that the diffractometer has been properly levelled, aligned with the beam and zeroed. See the SPEC fourc manual.

Before moving in hkl space you must calculate a UB matrix by specifying the crystal's lattice parameters (which define the B matrix) and finding two reflections (from which the U matrix defining any mismount can be inferred); and, optionally for surface-diffraction experiments, determine how the surface of the crystal is oriented with respect to the phi axis.

Once a UB matrix has been calculated, the diffractometer may be driven in hkl coordinates. A valid diffractometer setting maps easily into a single hkl value. However for a diffractometer with more than three circles there are excess degrees of freedom when calculating a diffractometer setting from an hkl value. Diffcalc provides modes for using up the excess degrees of freedom.

Diffcalc does not perform scans directly. Instead, Scannables that use diffcalc to map between reciprocal lattice space and real diffractometer settings are scanned using the Gda's (or minigda's) generic scan mechanism.

2.1 Theory

Thanks to Elias Vlieg for sharing his DOS based DIF software that Diffcalc has borrowed heavily from. The version of Diffcalc described here is based on papers by pHH. You. [?] and Busing & Levy [?]. (See also the THANKS.txt file.)

Getting Help

There are few commands to remember. If a command is called without arguments in some cases Diffcalc will prompt for arguments and provide sensible defaults which can be chosen by pressing enter.

Orientation. The helpub command lists all commands related with crystal orientation and the reference vector (often used with surfaces). See the *Orientation Commands* section at the end of this manual:

>>> help ub

HKL movement. The help hkl list all commands related to moving in reciprocal-lattice space. See the *Motion Commands* section at the end of this manual:

>>> help hkl

Call help on any command. e.g.:

```
>>> help loadub
loadub (diffcalc command):
loadub 'name' | num -- load an existing ub calculation
```

Diffcalc's Scannables

To list and show the current positions of your beamline's scannables use pos with no arguments:

>>> pos

Results in:

Energy and wavelength scannables:

energy 12.3984 wl: 1.0000

Diffractometer scannables, as a group and in component axes (in the real GDA these have limits):

```
sixc: mu: 0.0000 delta: 0.0000 gamma: 0.0000 omega: 0.0000 chi: 0.0000 phi: 0.0000
chi: 0.0000
delta: 0.0000
gamma: 0.0000
omega: 0.0000
phi: 0.0000
```

Dummy counter, which in this example simply counts at 1hit/s:

ct: 0.0000

Hkl scannable, as a group and in component:

```
hkl: Error: No UB matrix
h: Error: No UB matrix
k: Error: No UB matrix
l: Error: No UB matrix
```

Parameter scannables, used in some modes, these provide a scannable alternative to the *Motion* section. Some constrain of these constrain virtual angles:

alpha: --beta: --naz: --psi: --qaz: ---

and some constrain physical angles:

phi_con: --chi_con: --delta_con: --gam_con: --mu_con: ---

Crystal orientation

Before moving in hkl space you must calculate a UB matrix by specifying the crystal's lattice parameters (which define the B matrix) and finding two reflections (from which the U matrix can be inferred); and, optionally for surfacediffraction experiments, determine how the surface of the crystal is oriented with respect to the phi axis.

5.1 Start a new UB calculation

A *UB calculation* contains the description of the crystal-under-test, any saved reflections, reference angle direction, and a B & UB matrix pair if they have been calculated or manually specified. Starting a new UB calculation will clear all of these.

Before starting a UB-calculation, the ub command used to summarise the state of the current UB-calculation, will reflect that no UB-calculation has been started:

```
>>> ub
<<< No UB calculation started >>>
```

A new UB-calculation calculation may be started and lattice specified explicitly:

```
>>> newub 'example'
>>> setlat '1Acube' 1 1 1 90 90 90
```

or interactively:

```
>>> newub
calculation name: example
crystal name: 1Acube
crystal system
1) Triclinic
2) Monoclinic
3) Orthorhombic
4) Tetragonal
5) Rhombohedral
```

(continues on next page)

(continued from previous page)

where a is unit cell basis vector in Angstroms for cubic crystal system.

The ub command will show the state of the current UB-calculation (and the current energy for reference):

```
>>> ub
UBCALC
  name:
             example
REFERNCE
         1.00000 0.00000 0.00000 <- set
  n_hkl:
SURFACE NORMAL
  n_phi: 0.00000 0.00000 1.00000 <- set
CRYSTAL
             1Acube
  name:
  a, b, c: 1.00000 1.00000 1.00000
            90.00000 90.00000 90.00000 Cubic
  B matrix: 6.28319 0.00000 0.00000
             0.00000 6.28319 0.00000
             0.00000 0.00000 6.28319
UB MATRIX
  <<< none calculated >>>
REFLECTIONS
  <<< none specified >>>
CRYSTAL ORIENTATIONS
  <<< none specified >>>
```

5.2 Load a UB calculation

To load the last used UB-calculation:

```
>>> lastub
Loading ub calculation: 'mono-Si'
```

To load a previous UB-calculation:

5.3 Generate a U matrix from two reflections

The normal way to calculate a U matrix is to find the position of **two** reflections with known hkl values. Diffcalc allows many reflections to be recorded. After adding first two reflections UB matrix will be calculated automatically. If there are multiple recorded reflections, the indices or tags can be provided to calcub command as arguments to calculate UB matrix from any two given reflections.

Find U matrix from two reflections:

```
>>> pos wl 1
wl: 1.0000
>>> c2th [0 0 1]
59.9999999999999
>>> pos sixc [0 60 0 30 90 0]
sixc: mu: 0.0000 delta: 60.0000 gam: 0.0000 eta: 30.0000 chi: 90.0000 phi: _
-0.0000
>>> addref [0 0 1]
>>> pos sixc [0 90 0 45 45 90]
sixc: mu: 0.0000 delta: 90.0000 gam: 0.0000 eta: 45.0000 chi: 45.0000 phi: _
-90.0000
>>> addref [0 1 1]
Calculating UB matrix.
```

Check that it looks good:

>>>	checkub								
	ENERGY	Н	K	L	H_COMP	K_COMP	L_COMP	TAG	
1	12.3984	0.00	0.00	1.00	0.0000	0.0000	1.0000		
2	12.3984	0.00	1.00	1.00	0.0000	1.0000	1.0000		

After adding another reflection we can use the first and the third reflections to recalculate UB matrix:

```
>>> addref [1 0 1]
>>> calcub 1 3
>>> checkub
    ENERGY
             Η
                  K
                       L
                           H_COMP
                                    K_COMP
                                            L_COMP
                                                      TAG
  12.3984 0.00 0.00 1.00 0.0000 0.0000
1
                                            1.0000
2 12.3984 0.00 1.00 1.00
                          0.0000 1.0000
                                            1.0000
3 12.3984 1.00 0.00 1.00
                          1.0000
                                   0.0000
                                            1.0000
```

5.4 Generate a U matrix from one reflection

To estimate based on first reflection only:

5.5 Edit reflection list

Use showref to show the reflection list:

owref									
NERGY	Н	K	L	MU	DELTA	GAM	ETA	CHI	PHI 👝
2.398	0.00	0.00	1.00	0.0000	60.0000	0.0000	30.0000	90.0000	0.0000
2.398	0.00	1.00	1.00	0.0000	90.0000	0.0000	45.0000	45.0000	90.0000
1	owref NERGY 2.398 2.398	owref NERGY H 2.398 0.00 2.398 0.00	owref NERGY H K 2.398 0.00 0.00 2.398 0.00 1.00	owref NERGY H K L 2.398 0.00 0.00 1.00 2.398 0.00 1.00 1.00	owref NERGY H K L MU 2.398 0.00 0.00 1.00 0.0000 2.398 0.00 1.00 1.00 0.0000	owref MU DELTA NERGY H K L MU DELTA 2.398 0.00 0.00 1.00 0.0000 60.0000 2.398 0.00 1.00 1.00 0.0000 90.0000	owref NERGY H K L MU DELTA GAM 2.398 0.00 0.00 1.00 0.0000 60.0000 0.0000 2.398 0.00 1.00 1.00 0.0000 90.0000 0.0000	owref NERGY H K L MU DELTA GAM ETA 2.398 0.00 0.00 1.00 0.0000 60.0000 0.0000 30.0000 2.398 0.00 1.00 1.00 0.0000 90.0000 0.0000 45.0000	owref NERGY H K L MU DELTA GAM ETA CHI 2.398 0.00 0.00 1.00 0.0000 60.0000 0.0000 30.0000 90.0000 2.398 0.00 1.00 0.0000 90.0000 0.0000 45.0000 45.0000

Use swapref to swap reflections:

Use delref to delete a reflection:

>>> delref 1

5.6 Generate a U matrix from two lattice directions

Another approach to calculate a U matrix is to provide orientation of **two** crystal lattice directions using addorient command after aligning sample in laboratory frame of reference. The first lattice direction should be aligned along the selected direction in the laboratory frame. For the purpose of finding azimuthal orientation in U matrix calculation it is sufficient for the projection of the second lattice direction to be aligned to the given orientation in the laboratory frame in the plane perpendicular to the first lattice orientation.

Find U matrix from two lattice directions:

```
>>> addorient [0 0 1] [0 0 1]
>>> addorient [1 0 0] [1 1 0]
Calculating UB matrix.
```

5.7 Calculate a UB matrix

Unless a U or UB matrix has been manually specified, a new UB matrix will be calculated after the second reflection has been found, or whenever one of the first two reflections is changed.

Use the command calcub to force the UB matrix to be calculated from the first two reflections. In case of using lattice orientations instead of reflections, use command orientub to force the UB matrix to be calculated from the first two orientations.

UB matrix can be calculated from any combination of two reflections and/or orientations by providing corresponding reflection/orientation tags or numbers as an argument to calcub. In case of using one reflection and one orientation it is recommended to use tags to avoid ambiguity.

If you have misidentified a reflection used for the orientation the resulting UB matrix will be incorrect. Always use the checkub command to check that the computed reflection indices agree with the estimated values:

>>>	checkub							
	ENERGY	Н	K	L	H_COMP	K_COMP	L_COMP	TAG
1	12.3984	0.00	1.00	1.00	0.0000	1.0000	1.0000	
2	12.3984	0.00	0.00	1.00	0.0000	0.0000	1.0000	

5.8 Calculate a U matrix from crystal mismount

U matrix can be defined from crystal mismount by using a rotation matrix calculated from a provided mismount angle and axis. setmiscut command defines new U matrix by setting it to a rotation matrix calculated from the specified angle and axis parameters. addmiscut command applies the calculated rotation matrix to the existing U matrix, i.e. adds extra mismount to the already existing one:

```
>>> setmiscut 5 [1 0 0]
n_phi: -0.00000 -0.08716 0.99619
n_hkl: 0.00000 0.00000 1.00000 <- set
normal:
    angle: 5.00000
    axis: 1.00000 -0.00000 0.00000</pre>
```

5.9 Manually specify U matrix

Set U matrix manually (pretending sample is squarely mounted):

5.10 Refining UB matrix with one reflection

UB matrix elements can be refined to match diffractometer settings and crystal orientation experimentally found for a given reflection with the corresponding reflection indices. refineub command rescales crystal unit cell dimensions to match with the found scattering angle value and recalculates mismount parameters to update U matrix:

```
>>> refineub [1 0 0]
current pos[y]: y
Unit cell scaling factor: 0.99699
Refined crystal lattice:
  a, b, c: 0.99699 0.99699
                                0.99699
            90.00000 90.00000 90.00000
Update crystal settings?[y]: y
Warning: the old UB calculation has been cleared.
        Use 'calcub' to recalculate with old reflections or
        'orientub' to recalculate with old orientations.
Miscut parameters:
      angle: 2.90000
      axis: -0.00000
                      1.00000 -0.00000
Apply miscut parameters?[y]: y
  n_phi: 0.67043 -0.00000 0.74198
  n_hkl: 0.00000 0.00000 1.00000 <- set
  normal:
     angle: 42.10000
     axis: 0.00000 1.00000
                                0.00000
```

5.11 Calculate UB matrix from multiple reflections

Using fitub command UB matrix can be optimised to find best fit for the selected list of reference reflections. For triclinic crystal system optimal solution is found by solving multivariate linear regression model, while for the higher symmetry systems it is found by running numerical optimiser:

```
>>> fitub 1 2 3 4
Fitting crystal lattice parameters...
Fitting orientation matrix...
Refined crystal lattice:
  a, b, c: 10.56348 10.56348 10.81364
            90.00000 90.00000 90.00000
Update crystal settings?[y]: y
Refined U matrix: 0.94559 -0.32489
                                   0.01762
                 0.32487 0.94575
                                   0.00437
                 -0.01809 0.00160 0.99984
Update U matrix?[y]: y
. . .
. . .
. . .
REFLECTIONS
 ENERGY
           Η
                K
                       L
                               PHI
                                      CHI
                                                EΤΑ
                                                         MU
                                                               DELTA
                                                                          GAM TAG
                                                                      -0.0042 🔒
1 8.000 0.00 0.00 8.00 -26.3000 89.0000 17.3034
                                                      0.0005 33.3569
→None
2 8.000 4.00 4.00 8.00
                           62.4273 53.4451 45.2680
                                                      0.0000 90.0825
                                                                       0.0000
→None
3 8.000 0.00 0.00 8.00
                           13.3485
                                    89.0097 35.0408
                                                      0.0000 69.9326
                                                                       0.0000
→None
4 8.000 4.00 4.00 8.00
                           63.2008 53.4096 44.9007
                                                      0.0000 90.1107
                                                                       0.0000 ...
→None
```

5.12 Set the reference vector

The reference vector can be used to define azimuthal direction within the crystal with which we want to orient the incident or diffracted beam. Orientation of the reference vector w.r.t the incident and diffracted beam is indicated using alpha and beta angles.

By default the reference vector is set parallel to the theta axis. That is, along the x-axis of the laboratory coordinate frame.

The ub command shows the current reference vector at the top its report (or it can be shown by calling setnphi or setnhkl with no args):

```
>>> ub
...
REFERNCE
n_phi: 1.00000 0.00000 0.00000
n_hkl: 1.00000 0.00000 <- set
...</pre>
```

The <- set label here indicates that the reference vector is set in the reciprocal lattice space. In this case, therefore, its direction in the laboratory coordinate frame is inferred from the UB matrix.

To set the reference vector in the phi coordinate frame use:

```
>>> setnphi [1 0 0]
```

To set the reference vector in the crystal's reciprocal lattice space use:

```
>>> setnhkl [1 0 0]
```

5.13 Set the surface normal vector

The orientation of the sample surface can be set using the surface normal vector defined either in laboratory coordinate system or reciprocal space. Orientation of the surface normal vector w.r.t the incident and diffracted beam is indicated using betain and betaout angles.

By default the surface normal vector is set parallel to the phi axis. That is, along the z-axis of the laboratory coordinate frame.

The ub command shows the current surface normal vector at the top its report (or it can be shown by calling surfnphi or surfnhkl with no args):

```
>>> ub
...
SURFACE NORMAL
n_phi: 0.00000 0.00000 1.00000 <- set
n_hkl: 0.00000 0.00000 1.00000
...</pre>
```

The <- set label here indicates that the surface normal vector is set in the laboratory coordinate frame. In this case, therefore, its direction in the crystal's reciprocal lattice space is inferred from the UB matrix.

To set the surface normal vector in the phi coordinate frame use:

```
>>> surfnphi [0 0 1]
```

To set the surface normal vector in the crystal's reciprocal lattice space use:

```
>>> surfnhkl [0 0 1] ...
```

Motion

Once a UB matrix has been calculated, the diffractometer may be driven in hkl coordinates. A given diffractometer setting maps easily into a single hkl value. However for a diffractometer with more than three circles there are excess degrees of freedom when calculating a diffractometer setting from an hkl value. Diffcalc provides many for using up the excess degrees of freedom.

By default Diffcalc selects no mode.

6.1 Constraining solutions for moving in hkl space

To get help and see current constraints:

```
>>> help con
. . .
>>> con
 DET
                REF
                               SAMP
  _____
                _____
                               _____
  delta
                a_eq_b
                              mu
                alpha
                              eta
  qam
  qaz
                beta
                              chi
  naz
                psı
bin_eq_bout
                psi
                              phi
                              mu_is_gam
                betain
                              bisect
                betaout
                              omega
      3 more constraints required
  1
      Type 'help con' for instructions
```

Three constraints can be given: zero or one from the DET and REF columns and the remainder from the SAMP column. Not all combinations are currently available. Use help con to see a summary if you run into troubles.

To configure four-circle vertical scattering:

```
>>> con gam 0 mu 0 a_eq_b
    gam : 0.0000
    a_eq_b
    mu : 0.0000
```

In the following the *scattering plane* is defined as the plane including the scattering vector, or momentum transfer vector, and the incident beam.

DETECTOR COLUMN:

- delta physical delta setting (vertical detector motion) del=0 is equivalent to qaz=0
- gam physical gamma setting (horizontal detector motion) gam=0 is equivalent to qaz=90
- qaz azimuthal rotation of scattering vector (about the beam, from horizontal)
- naz azimuthal rotation of reference vector (about the beam, from horizontal)

REFERENCE COLUMN:

- alpha incident angle to reference vector
- beta exit angle from reference vector
- **psi** azimuthal rotation about scattering vector of reference vector (from scattering plane)
- **a_eq_b** bisecting mode with alpha=beta. Equivalent to psi=90
- betain incident angle to sample surface
- betaout exit angle from sample surface
- **bin_eq_bout** bisecting mode with betain=betaout

SAMPLE COLUMN:

- mu, eta, chi & phi physical settings
- mu_is_gam force mu to follow gamma (results in a 5-circle geometry)
- bisect bisecting mode with scattering vector in chi-circle plane
- omega bisecting mode with omega angle between scattering vector and chi-circle plane

Diffcalc will report two other (un-constrainable) virtual angles:

- theta half of 2theta, the angle through the diffracted beam bends
- tau longitude of reference vector from scattering vector (in scattering plane)

6.2 Example constraint modes

There is sometimes more than one way to get the same effect.

Vertical four-circle mode:

```
>>> con gam 0 mu 0 a_eq_b # or equivalently:
>>> con qaz 90 mu 0 a_eq_b
>>> con alpha 1 # replaces a_eq_b
```

Horizontal four-circle mode:

```
>>> con del 0 eta 0 alpha 1  # or equivalently:
>>> con qaz 0 mu 0 alpha 1
```

Surface vertical mode:

>>> con naz 90 mu 0 betain 1

Surface horizontal mode:

>>> con naz 0 eta 0 betain 1

Z-axis mode (surface horizontal):

```
>>> con chi (-sigma) phi (-tau) betain 1
```

where sigma and tau are the offsets required in chi and phi to bring the surface normal parallel to eta. betain will determine mu directly leaving eta to orient the planes. Or:

>>> con naz 0 phi 0 betain 1 # or any another sample angle

Z-axis mode (surface vertical):

>>> con naz 0 phi 0 betain 1 # or any another sample angle

6.3 Changing constrained values

Once constraints are chosen constrained values may be changed directly:

```
>>> con mu 10
    gam : 0.0000
    a_eq_b
    mu : 10.0000
```

or via the associated scannable:

```
>>> pos mu_con 10
mu_con: 10.00000
```

6.4 Configuring limits and cuts

Diffcalc uses motor limits set in GDA when used from GDA client running on a beamline. The standalone console version maintains its own limits on axes. These limits will be used when choosing solutions. If more than one detector solution exists Diffcalc will ask you to reduce the the limits until there is only one. However if more than one solution for the sample settings is available it will choose one that is closest to the current diffractometer orientation.

Use the hardware command to see the current limits and cuts:

(continues on next page)

(continued from previous page)

```
chi (cut: -180.0)
phi (cut: 0.0)
Note: When auto sector/transforms are used,
cuts are applied before checking limits.
```

To set the limits in standalone Diffcalc session:

```
>>> setmin delta -1
>>> setmax delta 145
```

To set a cut:

```
>>> setcut phi -180
```

This causes requests to move phi to be between the configured -180 and +360 degress above this. i.e. it might dive to -10 degrees rather than 350.

6.5 Moving in hkl space

Configure a mode, e.g. four-circle vertical:

```
>>> con gam 0 mu 0 a_eq_b
gam : 0.0000
a_eq_b
mu : 0.0000
```

Simulate moving to a reflection:

```
>>> sim hkl [0 1 1]
sixc would move to:
    mu : 0.0000
 delta : 90.0000
  gam : 0.0000
   eta: 45.0000
  chi: 45.0000
  phi: 90.0000
 alpha : 45.0000
  beta : 45.0000
betain :
          30.0000
betaout :
          30.0000
   naz :
          35.2644
   psi :
          90.0000
   qaz : 90.0000
   tau : 45.0000
 theta : 45.0000
ttheta : 90.0000
```

Move to reflection:

```
>>> pos hkl [0 1 1]
hkl: h: 0.00000 k: 1.00000 l: 1.00000
```

(continues on next page)

(continued from previous page)

```
>>> pos sixc
sixc: mu: 0.0000 delta: 90.0000 gam: 0.0000 eta: 45.0000 chi: 45.0000 phi: _
$90.0000
```

Simulate moving to a location:

Scanning in hkl space

All scans described below use the same generic scanning mechanism provided by the GDA system or by minigda. Here are some examples.

7.1 Fixed hkl scans

In a 'fixed hkl scan' something (such as energy or Bin) is scanned, and at each step hkl is 'moved' to keep the sample and detector aligned. Also plonk the diffractometer scannable (sixc) on there with no destination to monitor what is actually happening and then throw on a detector (ct) with an exposure time if appropriate:

```
>>> #scan scannable_name start stop step [scannable_name [pos or time]]..
>>> scan en 9 11 .5 hkl [1 0 0] sixc ct 1
>>> scan en 9 11 .5 hklverbose [1 0 0] sixc ct 1
>>> scan betain 4 5 .2 hkl [1 0 0] sixc ct 1
>>> scan alpha_par 0 10 2 hkl [1 0 0] sixc ct 1
```

7.2 Scanning hkl

Hkl, or one component, may also be scanned directly:

```
>>> scan h .8 1.2 .1 hklverbose sixc ct 1
```

At each step, this will read the current hkl position, modify the h component and then move to the resulting vector. There is a danger that with this method k and l may drift. To get around this the start, stop and step values may also be specified as vectors. So for example:

>>> scan hkl [1 0 0] [1 .3 0] [1 0.1 0] ct1

is equivilant to:

>>> pos hkl [1 0 0] >>> scan k 0 .3 .1 ct1

but will not suffer from drifting. This method also allows scans along any direction in hkl space to be performed.

7.3 Multidimension scans

Two and three dimensional scans:

```
>>> scan en 9 11 .5 h .9 1.1 .2 hklverbose sixc ct 1
>>> scan h 1 3 1 k 1 3 1 l 1 3 1 hkl ct 1
```

CHAPTER $\mathbf{8}$

Commands

8.1 Orientation Commands

STATE	
<pre>- newub {'name'}</pre>	start a new ub calculation name
– loadub 'name' num	load an existing ub calculation
– lastub	load the last used ub calculation
– listub	list the ub calculations available to load
– rmub 'name' num	remove existing ub calculation
– saveubas 'name'	save the ub calculation with a new name
LATTICE	
– setlat	interactively enter lattice parameters (Angstroms and Deg)
– setlat name a	assumes cubic
– setlat name a b	assumes tetragonal
– setlat name a b c	assumes ortho
– setlat name a b c gamma	assumes mon/hex with gam not equal to 90
– setlat name a b c alpha beta gamma	arbitrary
– c2th [h k l]	calculate two-theta angle for reflection
– hklangle [h1 k1 l1] [h2 k2 l2]	calculate angle between [h1 k1 l1] and [h2 k2 l2] crystal planes
REFERENCE	
– setnphi {[x y z]}	sets or displays n_phi reference
– setnhkl {[h k l]}	sets or displays n_hkl reference
SURFACE NORMAL	
– surfnphi {[x y z]}	sets or displays surface normal vector in lab space
– surfnhkl {[h k l]}	sets or displays surface normal vector in reciprocal space
REFLECTIONS	
– showref	shows full reflection list
– addref	add reflection interactively
– addref [h k l] { 'tag' }	add reflection with current position and energy
– addref [h k l] (p1,, pN) energy { 'tag' }	add arbitrary reflection

	Table 1 – continued from previous page
– editref num	interactively edit a reflection
– delref num	deletes a reflection (numbered from 1)
– clearref	deletes all the reflections
– swapref	swaps first two reflections used for calculating U matrix
- swapref num1 num2	swaps two reflections (numbered from 1)
CRYSTAL ORIENTATIONS	
– showorient	shows full list of crystal orientations
– addorient	add crystal orientation interactively
- addorient [h k l] [x y z] {'tag'}	add crystal orientation in laboratory frame
– editorient num	interactively edit a crystal orientation
– delorient num	deletes a crystal orientation (numbered from 1)
– clearorient	deletes all the crystal orientations
– swaporient	swaps first two crystal orientations used for calculating U matrix
- swaporient num1 num2	swaps two crystal orientations (numbered from 1)
UB MATRIX	
– fitub ref1 ref2 ref3	fit UB matrix to match list of provided reference reflections
– checkub	show calculated and entered hkl values for reflections
- setu {[[][]]}	manually set u matrix
- setub {[[][]]}	manually set ub matrix
– calcub	(re)calculate u matrix from ref1 and ref2
– calcub idx1 idx2	(re)calculate U matrix from reflections and/or orientations referred by indi
– orientub	(re)calculate U matrix from reflections and/or orientations
– orientub idx1 idx2	(re)calculate U matrix from the first two orientations referred by indices an
– trialub	(re)calculate U matrix from the first reflection (check carefully)
– trialub idx1	(re)calculate U matrix from reflection with index or tag idx only (check ca
– refineub {[h k 1]} {pos}	refine unit cell dimensions and U matrix to match diffractometer angles fo
- addmiscut angle {[x y z]}	apply miscut to U matrix using a specified miscut angle in degrees and a re
- setmiscut angle {[x y z]}	manually set U matrix using a specified miscut angle in degrees and a rota

8.2 Motion commands

CONSTRAINTS	
– con	list available constraints and values
- con <name> {val}</name>	constrains and optionally sets one constraint
- con <name> {val} <name> {val} <name> {val}</name></name></name>	clears and then fully constrains
- uncon <name></name>	remove constraint
HKL	
– allhkl [h k l]	print all hkl solutions ignoring limits
HARDWARE	
– hardware	show diffcalc limits and cuts
<pre>- setcut {name {val}}</pre>	sets cut angle
- setmin {axis {val}}	set lower limits used by auto sector code (None to clear)
- setmax {name {val}}	sets upper limits used by auto sector code (None to clear)
MOTION	
– sim hkl scn	simulates moving scannable (not all)
– sixc	show Eularian position
– pos sixc [mu, delta, gam, eta, chi, phi]	move to Eularian position(None holds an axis still)
- sim sixc [mu, delta, gam, eta, chi, phi]	simulate move to Eulerian positionsixc
– hkl	show hkl position
– pos hkl [h k l]	move to hkl position
– pos {h k 1} val	move h, k or l to val
– sim hkl [h k l]	simulate move to hkl position

Good luck — RobW

References

Bibliography

- [You1999] H. You. Angle calculations for a '4S+2D' six-circle diffractometer. J. Appl. Cryst. (1999). **32**, 614-623. (pdf link).
- [Busing1967] W. R. Busing and H. A. Levy. *Angle calculations for 3- and 4-circle X-ray and neutron diffractometers*. Acta Cryst. (1967). **22**, 457-464. (pdf link).