
EDIT2013 Documentation

Release 0.1.0

Shota Takahashi

June 06, 2016

1	How to obtain samples	3
2	What's inside	5
3	How to take data	7
3.1	Save data	7
3.2	Overwrite data	7
3.3	Change ADC channels	8
3.4	Install CAMAC driver	8
4	How to make plots	9
4.1	Change branch descriptor	9
5	How to use USB	11
6	Indices and tables	13

The Document for DAQ and ROOT macros used in EDIT2013 NEU course

Author Shota TAKAHASHI

Created 2013/03/12

Modified 2013/03/12

How to obtain samples

Just download “sample.zip” and expand it in any place you want.

What's inside

daq/	DAQ programs
data/	data
macros/	ROOT macros
plots/	save plots here if neccessary

How to take data

DAQ programs are in `daq/` directory. Move to this directory first. Then execute program named **adc** or **multiadc**. You need to specify `NEVENT` (event number) and `OFN` (output filename) as arguments. Output file is created in space-separated-text format.

Sample command 1

```
$ cd daq
$ ./adc NEVENT OFN (additional info).
```

Sample command 2 If you run without any argument, it will print usage.

```
$ ./adc
## --> usage will be printed
```

ADVANCED You can also set additional arguments if you want, like below

```
$ ./adc NEVENT OFN MPPC_ID BIAS_VOLTAGE
```

3.1 Save data

It is **REALLY** important to name data with easy-to-remember-and-handle convention. (If you named `test1.txt`, `test2.txt`, ... and so on,

you will get depressed when you came back to check data next week.)

My favorite way is to separate data by day. (it depends by purpose, though) I prefer `OFN` to be named like `../data/YMD/ANYTHING%03d.txt`, where

YMD	date (YYYYMMDD format)
ANYTHING	any words you want
%03d	run number, this makes easy to use ROOT macros below.

3.2 Overwrite data

DAQ program has non-overwriting feature to prevent **ACCIDENTAL** overwriting. In case of overwriting, remove that file first.

```
$ ./adc 100 SAME_NAMED_FILE
## --> Error: 'SAME_NAMED_FILE' already exists.

$ rm SAME_NAMED_FILE
$ ./adc 100 SAME_NAMED_FILE
```

3.3 Change ADC channels

Number of ADCs and its channels are defined as variable(array) named `ModId` and `ChId` in `adc.cc` (`multiadc.cc`). Modify these numbers and re-compile.

Sample command

```
$ emacs adc.cc
$ ... (edit adc.cc)
$ make
```

3.4 Install CAMAC driver

If you restarted PC for some reason, you need to install CAMAC driver before executing DAQ. We use `camdrv` for CAMAC driver. Its source is in `/opt/hep/kinoko/drv/camdrv/`.

Sample command

```
$ cd /opt/hep/kinoko/drv/camdrv
$ su (ask passwd to lecturer)
# make install
# dmesg
# exit
```

How to make plots

ROOT macro demos are in `macros` directory. Move to this directory first. Start ROOT CINT and load MACRO. Then execute as samples below.

*** IMPORTANT NOTICE *** These are just simple demos. I encourage you to look into the codes and improve it, or make your own. (If you know ROOT, you can do it =D)

Sample command 1 `mppcTest.C` to check one histogram

The function `histText()` is defined to return `TH1D*`. You need to create `TH1D*` and draw histogram.

```
$ root
root[] > .L mppcTest.C
root[] > TH1D *h1 = histText("h1", "../data/20130305/mpptest001.txt")
root[] > h1->Draw()
```

Sample command 2 `adcCalibration.C` to check three adc data at one time

The function “`histText`” is defined to return `TCanvas*`. You need to create `TCanvas*`. (Canvas will be drawn automatically)

```
$ root
root[] > .L adcCalibration.C
root[] > TCanvas *c1 = histText(32)
```

Sample command 3 `tracker.C` for eventdisplay

You don’t need to load `tracker.C`, just type

```
$ root tracker.C
```

or, in case you want to start from certain run number,

```
$ 'root tracker.C(3)'
```

4.1 Change branch descriptor

All macros use method `TTree::ReadFile(FILENAME, BRANCH_DESCRIPTOR)()` to create `TTree` from text file. If you changed orders or number of columns of output text format, modify `BRANCH_DESCRIPTOR` as needed.

How to use USB

Unfortunately, PCs used in EDIT2013 does not mount USB automatically. So please mount USB manually.

Sample command

```
$ dmesg
$ su (ask passwd to lecturer)
# mount /dev/sdcl /mnt/usb (or /media/usb)
# ls /mnt/usb
# ... (rsync or cp files)
# umount /mnt/usb
# exit
```

Indices and tables

- `genindex`
- `modindex`
- `search`