
leicascanningtemplate Documentation

Release 0.0.3

Arve Seljebu

April 29, 2015

Contents

```
class leicascanningtemplate.ScanningTemplate(filename)
```

Bases: object

Python object of Leica LAS Matrix Screener Scanning Template XML. Provides easy access to elements via attributes:

```
>>> tmpl = ScanningTemplate('{ScanningTemplate}tmpl.xml')
>>> # attributes of MatrixScreenerTemplate/ScanningTemplate/Properties
>>> print(tmpl.properties.attrib)
```

Parameters `filename` (*str*) – XML to load.

filename

str

Path XML-filename.

root

lxml.objectify.ObjectifiedElement

Objectified root of loaded XML. See <http://lxml.de/objectify.html#the-lxml-objectify-api>

add_well (*well_x*, *well_y*, *start_x*, *start_y*)

Add well with associated scan fields. `self.wells[0]` and `self.fields[0]` will be used as base. ScanWellData will be added to ScanWellArray and ScanFieldData to ScanFieldArray. The amount of fields added is decided by Properties/CountOfScanFields.

Parameters

- **well_x** (*int*) –
- **well_y** (*int*) –
- **start_x** (*int*) – In meters. FieldXCoordinate of first field in well.
- **start_y** (*int*) – In meters. FieldYCoordinate of first field in well.

Raises `ValueError` – If well or fields already exists.

count_of_assigned_jobs

Number of fields that have `attrib['JobAssigned']` set to true.

count_of_wells

Number of wells in x/y-direction of template.

Returns (*xs*, *ys*) number of wells in x and y direction.

Return type tuple

field (*well_x=1*, *well_y=1*, *field_x=1*, *field_y=1*)

ScanFieldData of specified field.

Parameters

- **well_x** (*int*) –
- **well_y** (*int*) –
- **field_x** (*int*) –
- **field_y** (*int*) –

Returns ScanFieldArray/ScanFieldData element.

Return type `lxml.objectify.ObjectifiedElement`

field_array

Short hand for self.root.ScanFieldArray

field_exists (well_x, well_y, field_x, field_y)

Check if field exists ScanFieldArray.

fields

All ScanFieldData elements.

Returns

Return type list of objectify.ObjectifiedElement

move_well (well_x, well_y, start_x, start_y)

Move well and associated scan fields. Spacing between fields will be what Properties/ScanFieldStageDistance is set to.

Parameters

- **well_x (int)** –
- **well_y (int)** –
- **start_x (int)** – In meters. FieldXCoordinate of first field in well.
- **start_y (int)** – In meters. FieldYCoordinate of first field in well.

Raises ValueError – If specified well or associated fields does not exist.

properties

Short hand for self.root.ScanningTemplate.Properties

remove_well (well_x, well_y)

Remove well and associated scan fields.

Parameters

- **well_x (int)** –
- **well_y (int)** –

Raises AttributeError – If well not found.

update_counts ()

Update counts of fields and wells.

update_start_position ()

Set start position of experiment to position of first field.

update_well_positions ()

Set well_attrib['FieldXStartCoordinate'] and well_attrib['FieldYStartCoordinate'] to FieldXCoordinate and FieldYCoordinate of first field in well.

well (well_x=1, well_y=1)

ScanWellData of specific well.

Parameters

- **well_x (int)** –
- **well_y (int)** –

Returns

Return type lxml.objectify.ObjectifiedElement

well_array

Short hand for self.root.ScanWellArray

well_attrib (well_x=1, well_y=1)

Attributes of specific well.

Parameters

- **well_x** (*int*) –
- **well_y** (*int*) –

Returns Attributes of ScanWellArray/ScanWellData.

Return type dict

well_exists (well_x, well_y)

Check if well exists in ScanWellArray.

well_fields (well_x=1, well_y=1)

All ScanFieldData elements of given well.

Parameters

- **well_x** (*int*) –
- **well_y** (*int*) –

Returns All ScanFieldData elements of given well.

Return type list of lxml.objectify.ObjectifiedElement

wells

All ScanWellData elements.

Returns

Return type list of objectify.ObjectifiedElement

write (filename=None)

Save template to xml. Before saving template will update date, start position, well positions, and counts.

Parameters **filename** (*str*) – If not set, XML will be written to self.filename.

A

add_well() (leicascanningtemplate.ScanningTemplate method), [1](#)

C

count_of_assigned_jobs (leicascanningtemplate.ScanningTemplate attribute), [1](#)

count_of_wells (leicascanningtemplate.ScanningTemplate attribute), [1](#)

F

field() (leicascanningtemplate.ScanningTemplate method), [1](#)

field_array (leicascanningtemplate.ScanningTemplate attribute), [1](#)

field_exists() (leicascanningtemplate.ScanningTemplate method), [2](#)

fields (leicascanningtemplate.ScanningTemplate attribute), [2](#)

filename (leicascanningtemplate.ScanningTemplate attribute), [1](#)

L

leicascanningtemplate (module), [1](#)

M

move_well() (leicascanningtemplate.ScanningTemplate method), [2](#)

P

properties (leicascanningtemplate.ScanningTemplate attribute), [2](#)

R

remove_well() (leicascanningtemplate.ScanningTemplate method), [2](#)

root (leicascanningtemplate.ScanningTemplate attribute), [1](#)

S

ScanningTemplate (class in leicascanningtemplate), [1](#)

U

update_counts() (leicascanningtemplate.ScanningTemplate method), [2](#)

update_start_position() (leicascanningtemplate.ScanningTemplate method), [2](#)

update_well_positions() (leicascanningtemplate.ScanningTemplate method), [2](#)

W

well() (leicascanningtemplate.ScanningTemplate method), [2](#)

well_array (leicascanningtemplate.ScanningTemplate attribute), [2](#)

well_attrib() (leicascanningtemplate.ScanningTemplate method), [3](#)

well_exists() (leicascanningtemplate.ScanningTemplate method), [3](#)

well_fields() (leicascanningtemplate.ScanningTemplate method), [3](#)

wells (leicascanningtemplate.ScanningTemplate attribute), [3](#)

write() (leicascanningtemplate.ScanningTemplate method), [3](#)