
PyJibe Documentation

Release 0.0.post50+gb5699d6c3

Paul Müller

Jul 08, 2026

CONTENTS

1	Getting started	3
1.1	Installation	3
1.2	Update	3
1.3	Supported file formats	3
1.4	How to cite	3
2	User Interface	5
2.1	Overview	5
2.2	Force-Distance Analysis	5
2.3	Tools	11
2.4	Advanced options	11
3	Tutorials	13
3.1	T1: Stiffness of two PAAm hydrogels	13
4	Quick Guides	19
4.1	Importing a pre-computed training set	19
4.2	Loading Extensions (nanite fit models)	19
5	Code Signing Policy	21
5.1	Project Integrity	21
5.2	Privacy Policy	21
6	Changelog	23
6.1	version 0.16.4	23
6.2	version 0.16.3	23
6.3	version 0.16.2	23
6.4	version 0.16.1	23
6.5	version 0.16.0	23
6.6	version 0.15.10	24
6.7	version 0.15.9	24
6.8	version 0.15.8	24
6.9	version 0.15.7	24
6.10	version 0.15.6	24
6.11	version 0.15.5	24
6.12	version 0.15.4	24
6.13	version 0.15.3	24
6.14	version 0.15.2	24
6.15	version 0.15.1	24
6.16	version 0.15.0	25
6.17	version 0.14.3	25

6.18	version 0.14.2	25
6.19	version 0.14.1	25
6.20	version 0.14.0	25
6.21	version 0.13.4	25
6.22	version 0.13.3	25
6.23	version 0.13.2	26
6.24	version 0.13.1	26
6.25	version 0.13.0	26
6.26	version 0.12.1	26
6.27	version 0.12.0	26
6.28	version 0.11.2	26
6.29	version 0.11.1	26
6.30	version 0.11.0	27
6.31	version 0.10.1	27
6.32	version 0.10.0	27
6.33	version 0.9.4	27
6.34	version 0.9.3	27
6.35	version 0.9.2	28
6.36	version 0.9.1	28
6.37	version 0.9.0	28
6.38	version 0.8.8	28
6.39	version 0.8.7	28
6.40	version 0.8.6	29
6.41	version 0.8.5	29
6.42	version 0.8.4	29
6.43	version 0.8.3	29
6.44	version 0.8.2	29
6.45	version 0.8.1	29
6.46	version 0.8.0	29
6.47	version 0.7.6	29
6.48	version 0.7.5	30
6.49	version 0.7.4	30
6.50	version 0.7.3	30
6.51	version 0.7.2	30
6.52	version 0.7.1	30
6.53	version 0.7.0	30
6.54	version 0.6.8	30
6.55	version 0.6.7	31
6.56	version 0.6.6	31
6.57	version 0.6.5	31
6.58	version 0.6.4	31
6.59	version 0.6.3	31
6.60	version 0.6.2	31
6.61	version 0.6.1	31
6.62	version 0.6.0	31
6.63	version 0.5.7	32
6.64	version 0.5.6	32
6.65	version 0.5.5	32
6.66	version 0.5.4	32
6.67	version 0.5.3	32
6.68	version 0.5.2	32
6.69	version 0.5.1	32
6.70	version 0.5.0	33
6.71	version 0.4.4	33

6.72	version 0.4.3	33
6.73	version 0.4.2	33
6.74	version 0.4.1	33
6.75	version 0.4.0	33
7	Bilbliography	35
8	Indices and tables	37
	Bibliography	39

PyJibe

PyJibe is a user interface for data analysis in atomic force microscopy (AFM) with an emphasis on biological specimens, such as tissue sections or single cells. This is the documentation of PyJibe version 0.0.post50+gb5699d6c3.

This documentation is also available as a [PDF](#).

GETTING STARTED

1.1 Installation

PyJibe can be installed via multiple channels:

1. **Windows installer:** Download the latest version for your architecture (e.g. `PyJibe_X.Y.Z_win_64bit_setup.exe` for 64bit Windows) from the official [release page](#).
2. **macOS:** Download the latest version (`PyJibeApp_X.Y.Z.dmg`) from the official [release page](#).
3. **Python 3 with pip:** PyJibe can easily be installed with [pip](#):

```
pip install pyjibe
```

To start PyJibe, simply run `python -m pyjibe` or just `pyjibe` in a command shell.

1.2 Update

If you install a newer (or older) version of PyJibe, the previously installed version will be automatically uninstalled.

If you installed `pyjibe` with `pip`, you may upgrade it with:

```
pip install --upgrade pyjibe
```

1.3 Supported file formats

PyJibe relies on the [afmformats](#) package. A list of supported file formats can be found [here](#).

1.4 How to cite

If you use PyJibe in a scientific publication, please cite it with:

Paul Müller and others (2019), PyJibe version X.X.X: Atomic force microscopy data analysis [Software]. Available at <https://github.com/AFM-analysis/PyJibe>.

Please also consider citing Müller et al., *BMC Bioinformatics* (2019) [[MAM+19](#)] (for fitting and rating force-indentation data).


USER INTERFACE

2.1 Overview

When you first start PyJibe, you will be greeted with an empty user interface. To open experimental data, choose one of the corresponding options in the *File* menu.

2.2 Force-Distance Analysis

Force-distance (FD) measurements consist of an approach and a retract curve. The point of maximum indentation is usually defined by the acquisition software as a *setpoint* using an active feedback loop.




The screenshot shows a typical data analysis in PyJibe. A *.jpk-force-map* file has been loaded. The plot at the top shows the experimental FD curve (approach light blue, retract light red), the model fit (dark blue), and the fitting region (white region, bounded by gray background in the upper plot). Below are the fit residuals. On the right are the main controls (discussed in the sections below). At the bottom is the curve list, highlighting which curve is currently shown at the top. At the bottom right are the *curve list controls*.

Note

Notice how the fit residuals become small around the contact point? PyJibe comes with the feature *Suppress residuals near contact point* (in the fit tab on the right). This option suppresses the contributions of the FD curve near the contact point. This feature was designed for biological samples, where the physical tip-sample interactions near the “contact point” are not well-understood [MAM+19]. The contributions near the contact point are suppressed to give the remainder of the indentation curve more weight.

The dataset shown resembles a 2D scan of AFM force-distance curves from a live zebrafish spinal cord microtome section and is available on figshare [MMG20] (same data as in [MKH+20] fig. 1e and [MAM+19] fig. 3a,d,g).

2.2.1 Tab: Preprocess



The screenshot shows the 'Preprocess' tab selected in a software interface. The tab bar at the top includes 'Preprocess', 'Fit', 'E(δ)', 'Plot', 'Info', and 'QMap'. Below the tab bar, there is a 'Load preset:' dropdown menu set to 'Recommended'. The main area contains several options with checkboxes and dropdown menus:


- ☒ tip-sample separation
- ☒ contact point estimation
 - Deviation from baseline
- ☐ baseline slope correction
 - Correct up until contact point
 - Temporal drift correction
- ☒ baseline offset correction
- ☒ segment discovery
- ☐ monotonic height data

At the bottom right, there is a link: [Preprocessing in the nanite docs](#).

Before starting to fit FD data, it is important to perform the correct preprocessing steps. Usually, this includes computing the tip position and correcting for a force and tip position offset.

The preprocessing pipeline is defined via the checkable options in this tab. If you hover over the options, a tool tip is displayed with an explanation. If your analysis pipeline strongly depends on a good estimate of the contact point, then you have multiple options which are described in detail in the [nanite docs](#).

2.2.2 Tab: Fit



Preprocess Fit E(δ) Plot Info QMap

Segment: Approach

Model: parabolic indenter (Hertz)

GCF: 1.000000

Range: absolute

Set indentation depth globally

Left (δ) [μm]: -inf

Right [μm]: 2.3

☒ Suppress residuals near contact point

☒ % tip radius 10.0 ☐ manual [μm] 1.86

Initial Parameters:

	fix	value	min	max
Young's Modulus [Pa]	☐	3000	0.0	inf
Tip Radius [μm]	☒	18.65	0.0	inf
Poisson's Ratio	☒	0.5	0.0	0.5
Contact Point [nm]	☐	0	-inf	inf
Force Baseline [pN]	☐	0	-inf	inf

Results:

Young's Modulus [Pa]	68.881
Contact Point [nm]	-875.32
Force Baseline [nN]	21.88

PyJibe offers a variety of options for fitting FD data. In the upper part, you may select the segment (approach or retract), the axes, the fit model, and the fitting range.

Range:

The fitting range can be set to absolute (relative to the X axis after preprocessing) or relative to the contact point. In the latter case, the fit is repeated four times, updating the new contact point at each iteration. You have three options for setting the indentation depth (left part of the fit interval):

1. globally: For each FD curve, the same indentation depth is used (after preprocessing has been applied).
2. individually: You select the indentation depth for each FD curve separately.
3. guess optimal indentation depth: The optimal indentation depth is guessed using the *E()*-curve.

Suppress residuals near contact point:

This option multiplies the residuals with a linear ramp, left and right of the contact point. The ramp has a value of 0 at the contact point and a value of 1 at the distance defined by the controls. All other residuals are not changed. This feature was designed for biological samples, where the physical tip-sample interactions near the “contact point” are not well-understood [MAM+19]. The contributions near the contact point are suppressed to give the remainder of the indentation curve more weight.

Ancillary Parameters:

Ancillary parameters (not shown in the screenshot here), are parameters that are defined in the fitting model. They are computed prior to the fit and can be set as fitting parameters. Standard models usually do not have ancillary parameters.

Initial Parameters:


These are the parameters set initially during fitting. In the table, you may choose which parameters should remain

fixed during fitting, what the initial values should be, and in which interval this value may be varied.

Results:

The fit results of the varied parameters are shown here.

2.2.3 Tab: E()



The $E()$ curve is used to test whether the fitted Young's modulus E is dependent on the fitting interval. For a reliable fit, the $E()$ curve should exhibit a plateau. This is the case for the present example, where the fitted value of E does not vary much in the range from $-2.5\mu\text{m}$ to $-4.8\mu\text{m}$ (maximum indentation).

The control at the top is identical to the control in the [fitting tab](#). You may choose the number of samples of the $E()$. The controls for setting the indentation depth manually are enabled when you select *set indentation depth individually* in the dropdown menu. Below the plot you have the options to export the $E()$ curve as an image or as a data file (.tsv). If you would like to export all $E()$ curves of the dataset, you can do so via the *Export* menu of the main window.

2.2.4 Tab: Plot

Preprocess Fit **E(δ)** Plot Info QMap

☒ Autoscale X-axis to fitting range

X-axis minimum -5.00

X-axis maximum 2.00

☒ Autoscale Y-axis to fitting range

Y-axis minimum -0.50

Y-axis maximum 7.00

The plotting tab allows you to modify the plot settings. By default, the plot is adapted to the fitting interval. You may, however, modify the axes ranges manually by unchecking the check boxes.

2.2.5 Tab: Info

Preprocess Fit E(δ) Plot **Info** QMap

Acquisition:
 Cantilever spring constant: 10.951 mN/m
 Feedback mode: contact
 Imaging modality: force-distance
 Sensitivity: 84.671 nm/V

Ancillaries:
 Maximum indentation: 3.8619 μm

Dataset:
 Axial piezo range: 39.932 μm
 Dataset index within the experiment: 210
 Duration of approach segment: 5 s
 Duration of retract segment: 2 s
 Duration: 6.8215 s
 Number of segments: 2
 Piezo speed of approach segment: 10 $\mu\text{m/s}$
 Piezo speed of retract segment: 25 $\mu\text{m/s}$
 Sampling rate of approach segment: 1928.6 Hz
 Sampling rate of retract segment: 2000 Hz
 Size of the dataset in points: 13643
 Target indentation force: 4 nN

QMap:
 Horizontal center of grid: -1159.8 μm
 Horizontal grid image size: 500 μm
 Horizontal grid position index: 10
 Horizontal grid shape: 25 px
 Horizontal position: -1199.8 μm
 Vertical center of grid: -1614.1 μm
 Vertical grid image size: 500 μm
 Vertical grid position index: 8
 Vertical grid shape: 25 px
 Vertical position: -1694.1 μm

The info tab shows metadata related to the currently shown curve. Unknown values are indicated as *nan*.

Acquisition:

Settings of the acquisition software

Ancillaries:

Ancillary parameters are computed for each model. Here, the *maximum indentation* is listed, which is the difference between the fitted contact point and the value of the tip position where the fitted curve has its maximum.

Fit models may have their own specific ancillary parameters.

Dataset:

Measurement parameters of the dataset

QMap:

If the curve is part of a quantitative map (2D FD scan), then the scan grid properties and the curve position are listed.


Setup:

Information about the AFM setup used

Storage:

Information about the data file, which are not particularly important for the experiment, but allow to identify a dataset or curve

2.2.6 Tab: QMap



If the dataset consists of a quantitative map, then this map is shown here. It serves as an interactive overview of the dataset. You can choose which data to plot (e.g. piezo height, fitted parameters, curve rating, scan order), which colormap to use, and how the colormap should be scaled. The coordinate axes are identical to those shown in the info tab.

The red square in the plot indicates which curve is currently shown. You may click on the plot to select curves manually (as opposed to using the curve list at the bottom of the window).

Below the plot are controls for exporting the qmap as an image or as a text file.

2.2.7 Curve list controls

The curve list controls perform operations on the entire FD curve list. The *Apply Model and Fit All* button does exactly what it says.

Curve rating allows you to rate each indentation curve. This is useful e.g. when you would like to sort out bad curves automatically. In this example, all curves with a rating below a threshold of 4.5 are excluded from further analysis using the *zef18 + Extra Trees* rating scheme [MAM+19].

Below the curve rating options, you may enable or disable autosaving of the fit results. Only the results for curves that are selected as “use” in the curve list are exported.

For more information on rating, please have a look at the [nanite documentation](#). The button at the bottom starts the PyJibe curve rater which is compatible to the nanite rating workflow. If you would like to import your own training set, please read the quick guide [Importing a pre-computed training set](#).

2.3 Tools

2.3.1 Data conversion

Under *Tools | Convert AFM data...*, you can select raw AFM data files and convert them to the CSV or HDF5 file formats. PyJibe uses the [afmformats](#) library to load and save AFM data.

2.4 Advanced options

2.4.1 Extensions

Under *Edit | Preferences | Extensions*, you may import your own custom fit models. See [Loading Extensions \(nanite fit models\)](#) for more information.

2.4.2 Expert and Developer mode

Under *Edit | Preferences | Advanced*, you can enable expert and developer mode.

- Expert mode:
 - Enables you to manually set the minimizer method and the corresponding keyword arguments applied by `lmfit`.
- Developer mode:
 - Everything that is available in expert mode.
 - Allows you to open experimental data that were not recorded via force-distance modality (e.g. creep-compliance). Related issues are [nanite #11](#) and [afmformats #15](#).

- Slows down loading of large datasets (because the modality has to be determined first).
- Displays and exports hidden model parameters (parameters whose name starts with an underscore _).

3.1 T1: Stiffness of two PAAm hydrogels

3.1.1 Introduction

In this tutorial, you will quantify the stiffness of two polyacrylamide (PAAm) hydrogels. The gels were produced according to the protocol in [RLY+17] with an expected Young’s modulus of 0.8 to 1.2 kPa (compliant gel) and 25 to 28 kPa (stiff gel).


3.1.2 Prerequisites

- PyJibe version 0.6.7 or above
- Complete dataset available on figshare [RMG20]

3.1.3 Data analysis

To make data analysis easier, we first divide the dataset by copying the files of the compliant gel (*PAAm_Compliant_*.jpk-force*) and the stiff gel (*PAAm_Stiff_*.jpk-force*) into two separate folders “compliant” and “stiff”.


In PyJibe, we load each of the folders using *File | Open bulk data*. We now have two subwindows called *Force-Distance #1* and *Force-Distance #2* (screenshot below).




We change the analysis pipeline in each of the subwindows as follows:

- In the *Fit* tab, set model to *spherical indenter (Sneddon, approximative)*. The measurements were performed using a spherical indenter. The *approximative* description indicates, that this model is an analytical representation of the Sneddon model (as used by the JPK analysis software). See [nanite.model.model_sneddon_spherical_approximation](#) for more information.
- In the *Fit* tab, set the right range to 2 μm. This primarily helps with the visualization of the data.
- In the *Fit* tab, uncheck the *Suppress residuals near contact point* check box. This feature is not necessary for PAAm gels.
- In the *Fit* tab, set the initial parameter *Tip Radius* to 5 μm.

Now, click the button *Apply Model and Fit All*. The results should look like this, for the compliant gel...



...and for the stiff gel:



There are several things to note here:

- The fit of the compliant gel looks much better than the fit of the stiff gel: There seems to be a systematic (sine-like) deviation from the fit. One could only speculate about the reason for this deviation. However, compared to the maximal indentation force (setpoint at 32 nN, see *Info* tab), the deviation is comparatively small.
- Furthermore, the stiff gel measurements exhibit a “ringing”, which is clearly visible in the approach part. This is a result of the *force-modulation* feedback mode. The gel is not probed in *contact* mode, but in a mixture between contact mode and intermittent mode (See the NanoWizard User Manual v. 4.2, sec. 5.7.). Maybe this ringing is not visible for the compliant gel, because of the setpoint at 6n N (The setpoint likely defines the amplitude of the force-modulation feedback mode and for the compliant gels, this amplitude might be below the noise level).
- The stiff gel gets a better rating than the compliant gel with the *zef18 + Extra Trees* rating scheme (please see the [nanite rating workflow](#) and [MAM+19] for how rating works). Of course, this observation is misleading - it nicely illustrates a limit of machine learning. The *zef18* training set was created using zebrafish spinal cord section measurements. In the context of hydrogels, it does not make much sense, although it appears to be consistent and gives both, compliant and stiff gels, a “good” rating.

3.1.4 Results

You might have realized that PyJibe creates the file *pyjibe_fit_results_leaf.tsv* in each of the measurement folders (if the *Autosave fit results as .tsv* check box is checked). These files contain (amongst other things) the fit results of each curve. With a simple Python script, we can visualize the Young’s modulus of the two gels:

```
import matplotlib.pyplot as plt
import seaborn as sns

data_compl = pandas.read_table("./compliant/pyjibe_fit_results_leaf.tsv")
data_stiff = pandas.read_table("./stiff/pyjibe_fit_results_leaf.tsv")

sns.set_style("darkgrid")

plt.subplot(121, title="compliant hydrogel")
sns.boxplot("Young's Modulus [Pa]", data=data_compl, fliersize=0,
            color="#d6ff7d")
sns.swarmplot("Young's Modulus [Pa]", data=data_compl,
              size=5, color=".3", linewidth=0)

plt.subplot(122, title="stiff hydrogel")
sns.boxplot("Young's Modulus [Pa]", data=data_stiff, fliersize=0,
            color="#98ff80")
sns.swarmplot("Young's Modulus [Pa]", data=data_stiff,
              size=5, color=".3", linewidth=0)

plt.show()
```




Fig. 3.1: Comparison of the hydrogels. Note that the X axes are scaled differently.


The compliant hydrogel has a Young's modulus of 1090 ± 12 Pa (mean $\pm$ SD) and the stiff hydrogel has a Young's modulus of 27676 ± 270 Pa. These values agree well with the values we expected initially.

QUICK GUIDES

4.1 Importing a pre-computed training set

A training set is a set of text files containing rating features and manual user ratings. Training sets are the basis of the rating scheme employed in PyJibe (and nanite) [MAM+19]. By default, PyJibe comes with the *zef18* training set [MMG18]. For general information about how manual rating and generation of the training set works, please refer to the [nanite rating workflow](#).

In order to import a training set into PyJibe, the text files (`train_feat_*.txt` and `train_response.txt`) must be zipped and named according to the scheme `ts_NAME.zip`, where `NAME` is a descriptive name.



The import in PyJibe is carried out by selecting *Add...* in the *Scheme* dropdown menu, located in the *Curve rating* box (lower right of a Force-Distance analysis). A dialog will ask for the zipped training set (e.g. `ts_organoids19.zip`).

The import is persistent, i.e. the training set is extracted to the user's configuration directory and thus does not have to be imported again the next time PyJibe is run.

4.2 Loading Extensions (nanite fit models)

Since version 0.12.0, PyJibe allows loading [nanite fit models](#).

You may need this feature if you are developing fit models yourself or if you received a fit model from your peer.


Note

If you installed PyJibe via installer (not via *pip*), then many extensions might not work due to software dependencies that those extensions might have. If this happens, please create an issue in the PyJibe repository so we can find a solution.

Warning

Extensions can be harmful. Please only load extensions that you received first-hand from people you trust.

You can load and manage extensions via the *Edit | Preferences* dialog in the *Extensions* tab.



CODE SIGNING POLICY

PyJibe uses free code signing provided by [SignPath.io](#), certificate by [SignPath Foundation](#). These binaries can be downloaded from the [release page](#).

5.1 Project Integrity

To verify the integrity of PyJibe, we manage access to the PyJibe repository via GitHub Team roles:

- **Triage:** The [Members team of AFM-analysis](#) has permission to triage (e.g. modify issues).
- **Core:** The [Core team of AFM-analysis](#) consists of people who are trusted to modify the source code in the project's version control system without additional reviews.
- There is no special role for *Committers*. External contributors or members of the *Triage* team must create a pull request which is reviewed by a *Core* team member.
- Furthermore, we enforce 2FA for every member of the AFM-Analysis GitHub organization.

5.2 Privacy Policy

This program will not transfer any information to other networked systems unless specifically requested by the user or the person installing or operating it.

Specifically, PyJibe does not require an internet connection to function.

CHANGELOG

List of changes in-between PyJibe releases.

6.1 version 0.16.4

- fix: correctly handle params states when swapping to new model ([#51](#))
- fix: ensure identical results when clicking per file and compute button ([#47](#), [#51](#))
- fix: remember user's choice upon curve open ([#44](#), [#45](#))
- enh: show χ^2 annotation and toggle fit line visibility ([#48](#), [#52](#))
- setup: support latest python versions ([#32](#), [#46](#))
- setup: support latest matplotlib versions ([#32](#), [#46](#))
- setup: bump afmformats to 0.18.7

6.2 version 0.16.3

- setup: bump afmformats to 0.18.6
- setup: bump nanite to 4.2.3

6.3 version 0.16.2

- setup: bump afmformats to 0.18.5 (chiaro file format)

6.4 version 0.16.1

- fix: plot extent does not match user settings ([#42](#))

6.5 version 0.16.0

- enh: update curve list even when ancillary computation fails
- ref: migrate from pyqt5 to pyqt6 ([#35](#))
- setup: bump nanite to 4.2.1

6.6 version 0.15.10

- ref: migrate from pkg_resources to importlib.resources (#33)

6.7 version 0.15.9

- build: fix installation directory
- build: disable unsigned installer

6.8 version 0.15.8

- docs: fix broken link to code signing

6.9 version 0.15.7

- build: maintenance

6.10 version 0.15.6

- build: maintenance

6.11 version 0.15.5

- build: maintenance

6.12 version 0.15.4

- build: fix signpath integration with test certificate

6.13 version 0.15.3

- build: finish signpath integration with test certificate

6.14 version 0.15.2

- setup: pin matplotlib>=3.3.7.5 until #32 is resolved
- build: migrate from appveyor to GHActions
- build: introduce code signing via signpath.org

6.15 version 0.15.1

- fix: properly deregister a model
- ref: remove iterative “sneddon_spher” model in developer mode
- setup: bump nanite from 4.0.0 to 4.1.1

6.16 version 0.15.0

- **BREAKING CHANGE:** Due to the upgrade of nanite, the rating values, e.g. for the “zef18” dataset might change slightly (on average by about 0.1). If you are excluding curves from your analysis based on rating, please rerun your analysis.
- setup: bump nanite from 3.7.0 to 4.0.0 (support training sets that contain nan-values for zero responses and infinity values; minor other fixes)

6.17 version 0.14.3

- setup: bump afmformats from 0.16.7 to 0.18.0 (support for specifying custom metadata when opening JPK files)
- ci: maintenance

6.18 version 0.14.2

- enh: support mirroring the input directory structure in the data conversion tool ([#29](#))
- docs: add info about data conversion tool

6.19 version 0.14.1

- maintenance release

6.20 version 0.14.0

- feat: support baseline corrections with linear slope (drift or shift)
- ref: generalize the view of preprocessing insights
- ref: support multiple options for nanite preprocessors
- setup: bump nanite from 3.5.2 to 3.7.0
- setup: bump Python version to 3.10

6.21 version 0.13.4

- ref: don't use deprecated matplotlib.cm.get_cmap
- ref: use packaging.version.parse to compare versions
- ref: fix matplotlib deprecation warning for set_data
- tests: don't be so strict

6.22 version 0.13.3

- setup: bump afmformats from 0.16.4 to 0.16.7
- setup: bump nanite from 3.5.0 to 3.5.2
- setup: require Python 3.9
- ci: update build pipelines

- tests: loosen allclose call

6.23 version 0.13.2

- fix: model parameters defined via expressions were not exported
- fix: only export hidden parameters in developer mode
- enh: export ancillary parameters in autosave

6.24 version 0.13.1

- fix: make sure fit parameters are identified by their label instead of by an index (which can be wrong, now that we have hidden parameters)
- fix: do not show “sneddon_spher” model in expert mode

6.25 version 0.13.0

- feat: allow to specify the geometrical correction factor for non-single-contact experiments (#16)
- fix: export fitting method and method kwargs with initial parameters
- enh: add expert mode which is just a parital developer mode (#22)
- enh: add links to docs in two places (#19)
- enh: show “Axes” in fitting tab only in developer mode (#21)
- setup: bump nanite from 3.4.0 to 3.5.0 (#16)

6.26 version 0.12.1

- build: fixed testing pipeline

6.27 version 0.12.0

- feat: support loading nanite model files as extensions
- feat: added automatic update check during startup
- ref: add preferences dialog
- setup: bump nanite from 3.3.1 to 3.4.0

6.28 version 0.11.2

- maintenance release

6.29 version 0.11.1

- fix: check for identical label when updating parameters (previously, a parameter that starts with the label of another parameter got assigned the wrong boundaries)
- enh: allow to hide model parameters via key names (#20)

- enh: show hidden parameters in developer mode
- ref: code-cleanup and simplification using latest nanite version
- setup: bump nanite from 3.2.1 to 3.3.1

6.30 version 0.11.0

- feat: allow custom fit methods with custom keyword arguments in developer mode (#18)
- fix: correctly handle dependencies in preprocessing tab
- fix: set POC method when loading recommended preprocessing (#17)
- setup: bump nanite from 3.1.3 to 3.2.1 (new POC estimation methods: Frechet distance and linear+polynomial fit)

6.31 version 0.10.1

- enh: use human-readable units in missing-metadata entry dialogs
- enh: change default unit for spring constant to mN/m
- ref: IndentationPreprocessor is deprecated in nanite
- setup: bump nanite to version 3.1.3 (ask for missing spring constant for tip position computation)

6.32 version 0.10.0

- **BREAKING CHANGE:** The default contact point estimation method changed. This means that your fitted contact points and any E(delta)-related results will change!
- feat: allow to select point of contact method (#15)
- enh: restructure preprocessing tab and avoid situations where preprocessing dependencies are not met (#15)
- enh: visualize contact point estimation algorithms in preprocessing tab
- setup: bump nanite from 2.0.0 to 3.1.1 (improved preprocessing)

6.33 version 0.9.4

- enh: only show exact sneddon model “sneddon_spher” in developer mode to avoid confusion
- setup: bump afmformats from 0.15.0 to 0.16.0 (improve support for AFM workshop data)
- setup: bump nanite from 1.7.6 to 2.0.0

6.34 version 0.9.3

- fix: possible fix for TypeError in InfDoubleSpinBox
- setup: bump afmformats from 0.14.3 to 0.15.0 as well as nanite from 1.7.6 to 1.7.8 (requirement for #14)

6.35 version 0.9.2

- ref: improve speed when displaying QMap data
- enh: add developer mode for loading creep-compliance data (#14)
- setup: bump afmformats from 0.13.3 to 0.14.3 (mostly speed, support AFM workshop maps)
- setup: bump nanite from 1.7.4 to 1.7.6 (mostly speed)
- ref: speed-up QMap visualization by caching nanite.QMap instances

6.36 version 0.9.1

- docs: minor update
- fix: regression where left fitting range could not be set
- fix: improved checks for InfDoubleSpinBox which acted up when entering absolutely sane floating point values

6.37 version 0.9.0

- feat: allow users to enter missing metadata during loading of data files (#5)
- fix: automatically append .h5 suffix when rating datasets
- enh: add busy cursor where applicable
- ref: directory for importing training datasets for rating changed (but I think this currently affects nobody)
- setup: dump appdirs dependency (moved to QSettings)
- setup: bump afmformats from 0.12.5 to 0.13.3 (mostly #5 and speed)
- setup: bump nanite from 1.7.3 to 1.7.4 (mostly #5)
- tests: increase coverage

6.38 version 0.8.8

- fix: handle bad curves properly when running bulk analysis (#13)
- ref: make sure there are no duplicate files when loading data (possibly fixes #12)
- ref: cleanup - call overloaded processEvents
- ref: properly decorate PyQtSlots
- tests: add helper function for setting up datasets

6.39 version 0.8.7

- fix: add workaround for macOS where the non-native PyQt5 dialogs did not allow browsing to “/Volumes” (it is now added to the sidebar)
- build: bump afmformats from 0.12.4 to 0.12.5 (OSError: “Too many files open” with the JPK file format)
- build: bump nanite from 1.6.2 to 1.7.3
- build: pin PyQt5 for GitHub releases
- tests: don’t used deprecated weight_cp function from nanite

6.40 version 0.8.6

- build: bump afmformats from 0.12.2 to 0.12.4 (AttributeError when reading HDF5 files)
- tests: deprecate setup.py test

6.41 version 0.8.5

- build: migrate from travisCI to GitHub Actions
- build: overhauled build pipeline (#9)

6.42 version 0.8.4

- build: add pyinstaller hooks for afmformats

6.43 version 0.8.3

- fix: main window not focused after startup
- fix: Windows installation location was confusing and did not coincide configuration file location
- fix: files are sorted before batch loading
- fix: change visualization of fitting region by gray wraparound instead of yellow highlight (not visible on all screens)
- setup: bump afmformats from 0.11.0 to 0.12.2 (improved JPK file format speed)

6.44 version 0.8.2

- setup: bump matplotlib to >=3 (NavigationToolbar2QT modifications)

6.45 version 0.8.1

- fix: correctly handle lmfit models with expressions
- setup: bump nanite from 1.6.0 to 1.6.2 (improve contact point retrieval for bad data, check models during registration)

6.46 version 0.8.0

- fix: exclude misplaced 'available' method from preprocessors
- setup: bump nanite from 1.5.5 to 1.6.0 (improved contact point estimation - this may introduce small changes in the fit result which should not be significant, but which justify a new minor version)

6.47 version 0.7.6

- maintenance release (fix macOS build)
- setup: bump afmformats from 0.10.2 to 0.11.0
- setup: bump nanite from 1.5.4 to 1.5.5

6.48 version 0.7.5

- maintenance release (fix build)

6.49 version 0.7.4

- setup: bump afmfomats from 0.10.0 to 0.10.2
- setup: bump nanite from 1.5.2 to 1.5.4
- build: add command-line option `--version` to just print the version and exit

6.50 version 0.7.3

- fix: update initial parameters when the user edits the left fitting interval boundary
- fix: conical indenter model did not work, because of a wrong (and unnecessary) entry in the *human_units* scheme.

6.51 version 0.7.2

- fix: setting bad fit range lead to `ValueError` during plotting

6.52 version 0.7.1

- fix: wrong SI units assignend to fit results parameters
- setup: bump afmformats from 0.9.0 to 0.10.0 (new formats: JPK (.jpk-qi-data), AFM workshop (.csv), NT-MDT (.txt exported by Nova))
- setup: bump nanite from 1.4.1 to 1.5.2

6.53 version 0.7.0

- feat: add data conversion tool
- enh: allow to select multiple files when loading data
- ref: changed order and clarified meaning of options in file menu
- setup: bump afmformats from 0.7.0 to 0.9.0 (fix HDF5 export, support new file formats)

6.54 version 0.6.8

- setup: bump afmformats from 0.5.1 to 0.7.0 (improved metadata view in info tab)
- setup: bump nanite from 1.4.1 to 1.5.1 (compatibility to afmformats)
- enh: fit results and parameter export dialog now supports new groups in afmformats 0.7.0
- enh: autosave now only stores fit results and rating parameters
- enh: remember scroll position in info tab

6.55 version 0.6.7

- setup: bump nanite from 1.4.0 to 1.4.1 (baseline is now a free parameter by default)
- docs: add first tutorial (PAAm gels)

6.56 version 0.6.6

- fix: typo lead to error when using rating threshold

6.57 version 0.6.5

- build: macOS build failed due to PyInstaller issue 4626

6.58 version 0.6.4

- fix: info tab did not display QMap metadata
- fix: do not fit the entire dataset when applying a rating threshold but take the values from previous ratings
- fix: deselecting curves of a qmap resulted in ValueError if QMap tab was selected
- fix: set correct display range for residuals
- docs: UI section for basics FD analysis
- build: Windows build broken since 0.5.6

6.59 version 0.6.3

- fix: support JPK data recorded in the “force-modulation” feedback mode (bump afmformats from 0.5.0 to 0.5.1)
- fix: make sure ancillary parameters are computed from the initial parameters set in the user interface (not from the default model parameters)

6.60 version 0.6.2

- enh: allow to select which metadata is exported
- fix: do not apply and fit to all before exporting metadata (user may have performed individual fits)
- setup: bump nanite from 1.3.0 to 1.4.0

6.61 version 0.6.1

- build: workaround for Pyinstaller issue 4626

6.62 version 0.6.0

- feat: initial support for ancillary parameters
- fix: always display parameter units
- fix: other minor UI bugs

- setup: bump afmformats from 0.4.1 to 0.5.0
- setup: bump nanite from 1.2.4 to 1.3.0

6.63 version 0.5.7

- enh: allow to also set right part of fitting range individually
- fix: improve layout of FD fitting panel
- setup: bump nanite from 1.2.3 to 1.2.4 (improved default params)

6.64 version 0.5.6

- setup: bump nanite from 1.2.2 to 1.2.3 (fixes issue with fits not being redone when the user changes an initial parameter)
- fix: IndexError when editing the fitting range

6.65 version 0.5.5

- maintenance release

6.66 version 0.5.4

- enh: write maximum indentation to statistics file (#3)
- fix: unicode characters were not read correctly from statistics file by libre office (added UTF-8 BOM)

6.67 version 0.5.3

- enh: improved visualization of meta data parameters in the Info tab
- setup: bump afmformats from 0.3.0 to 0.4.1
- setup: bump nanite from 1.2.0 to 1.2.2

6.68 version 0.5.2

- ref: bump nanite to 1.2.0
- ref: migrate to afmformats (0.3.0)

6.69 version 0.5.1

- fix: missing title in FD window and missing FD number in export menu
- enh: add help menu with link to docs, about, and software used (#1)
- docs: add quick guide for importing a nanite training set
- docs: add black/white logo

6.70 version 0.5.0

- feat: allow to import training sets produced with nanite
- ref: major code refactoring of force-distance interface and cleanup
- setup: bump nanite dependency to $\geq 1.1.1$

6.71 version 0.4.4

- ci: pin 'joblib==0.11.0' on travis-CI (workaround for infinite loop in macOS build, <https://github.com/pyinstaller/pyinstaller/issues/4067>)
- docs: add installation instructions

6.72 version 0.4.3

- maintenance release (migrate to GitHub infrastructure)
- experimental macOS builds

6.73 version 0.4.2

- enh: merge training set and regressor selection into scheme selection (prevents confusion for users that are unsure what to choose)
- enh: user rating is saved during rating, not afterwards
- enh: rating containers for user rating can be reused explicitly
- enh: qmap plot now better distinguishes between unavailable and invalid data
- ref: migrate user rating output to nanite

6.74 version 0.4.1

- fix: autoscaling according to fitting range did not work
- fix: not possible to load .jpk-force-map files
- ref: migrate qmap generation to nanite

6.75 version 0.4.0

- migration to GitHub

BILBLIOGRAPHY

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

BIBLIOGRAPHY

- [MKH+20] Stephanie Möllmert, Maria A. Kharlamova, Tobias Hoche, Anna V. Taubenberger, Shada Abuhattum, Veronika Kusch, Thomas Kurth, Michael Brand, and Jochen Guck. Zebrafish spinal cord repair is accompanied by transient tissue stiffening. *Biophysical Journal*, 118(2):448–463, jan 2020. doi:10.1016/j.bpj.2019.10.044.
- [MMG20] Stephanie Möllmert, Paul Müller, and Jochen Guck. Two-dimensional force-distance afm map of a live zebrafish spinal cord microtome section. *Figshare*, Jan 2020. doi:10.6084/m9.figshare.11709720.v1.
- [MAM+19] Paul Müller, Shada Abuhattum, Stephanie Möllmert, Elke Ulbricht, Anna V. Taubenberger, and Jochen Guck. Nanite: using machine learning to assess the quality of atomic force microscopy-enabled nano-indentation data. *BMC Bioinformatics*, 20(1):1–9, sep 2019. doi:10.1186/s12859-019-3010-3.
- [MMG18] Paul Müller, Stephanie Möllmert, and Jochen Guck. Atomic force microscopy indentation data of zebrafish spinal cord sections. *Figshare*, 11 2018. doi:10.6084/m9.figshare.7297202.v1.
- [RLY+17] Gonzalo Rosso, Ivan Liashkovich, Peter Young, Dominik Röhr, and Victor Shahin. Schwann cells and neurite outgrowth from embryonic dorsal root ganglions are highly mechanosensitive. *Nanomedicine: Nanotechnology, Biology and Medicine*, 13(2):493–501, feb 2017. doi:10.1016/j.nano.2016.06.011.
- [RMG20] Gonzalo Rosso, Paul Müller, and Jochen Guck. Atomic force microscopy indentation data of stiff and compliant polyacrylamide hydrogels. *Figshare*, Jan 2020. doi:10.6084/m9.figshare.11637675.v3.