

---

# **pylinkwrapper Documentation**

*Release 1.0*

**Nick DiQuattro**

March 16, 2016



<b>1</b>	<b>Installation</b>	<b>3</b>
<b>2</b>	<b>Connector Class</b>	<b>5</b>
<b>3</b>	<b>Sample Experiment</b>	<b>9</b>
<b>4</b>	<b>Indices and tables</b>	<b>11</b>
	<b>Python Module Index</b>	<b>13</b>



Contents:



---

## Installation

---

1. **Download** and place the pylinkwrapper folder somewhere you can easily reach in the file system.
2. Add the folder to Psychopy's python path as detailed here: <http://www.psychopy.org/recipes/addCustomModules.html>
3. Update pylink to the latest version provided by the SR Research development kit. It should be installed on the display computer.

That's it! Here's a quick demo of getting connected to the Eyelink:

```
import pylinkwrapper

win = visual.window(monitor='nickMon', fullScr=True, allowGUI=False, color=-1)
tracker = pylinkwrapper.connect(win, '1_nd')
```

Check out the documentation for the [Connector Class](#) class to see what functions are available.



---

## Connector Class

---

This class provides methods for interacting with the EyeLink from psychopy.

**class** `connector.Connect` (*window, edfname*)

Provides functions for interacting with the EyeLink via Pylink.

### Parameters

- **window** – Psychopy window object.
- **edfname** (*str*) – Desired name of the EDF file.

**calibrate** (*cnum=13, paval=1000*)

Calibrates eye-tracker using psychopy stimuli.

### Parameters

- **cnum** (*int*) – Number of points to use for calibration. Options are 3, 5, 9, 13.
- **paval** (*int*) – Pacing of calibration, i.e. how long you have to fixate each target in milliseconds.

**convert\_coords** (*x, y, to='eyelink'*)

Converts from degrees visual angle units to EyeLink Pixel units.

### Parameters

- **x** (*float or int*) – X coordinate in visual angle.
- **y** (*float or int*) – Y coordinate in visual angle.
- **to** – Direction of conversion. Options: 'eyelink' or 'psychopy'.

**Returns** Two values in order x, y

**draw\_ia** (*x, y, size, index, color, name*)

Draws square interest area in EDF and a corresponding filled box on eye-tracker display. Must be called after `set_trialid()` for interest areas to appear in the EDF.

### Parameters

- **x** (*float or int*) – X coordinate in degrees visual angle for center of check area.
- **y** (*float or int*) – Y coordinate in degrees visual angle for center of check area.
- **size** (*float or int*) – length of one edge of square in degrees visual angle.
- **index** (*int*) – number to assign interest area in EDF
- **color** (*int*) – color of box drawn on eye-tracker display (0 - 15)

- **name** (*str*) – Name of interest area in EDF

**draw\_text** (*msg*)

Draws text on eye-tracker screen.

**Parameters** **msg** (*str*) – Text to draw.

**end\_experiment** (*spath*)

Closes and transfers the EDF file.

**Parameters** **spath** (*str*) – Absolute file path of where to save the EDF file.

**fix\_check** (*size, ftime, button*)

Checks that fixation is maintained for certain time.

**Parameters**

- **size** (*float or int*) – Length of one side of box in degrees visual angle.
- **ftime** (*float*) – Length of time to check for fixation in seconds.
- **button** (*char*) – Key to press to recalibrate eye-tracker.

**get\_gaze** ()

Gets current gaze position of eye. Must be called between *record\_on()* and *record\_off()*. Sendlink must be set to True as well.

**Returns** list of coordinates in the form of [x, y].

**record\_off** ()

Stops recording.

**record\_on** (*sendlink=False*)

Starts recording. Waits 50ms to allow eyelink to prepare.

**Parameters** **sendlink** (*bool*) – Toggle for sending eye data over the link to the display computer during recording.

**sac\_detect** (*x, y, radius*)

Checks if current gaze position is outside a circular interest area.

**Parameters**

- **x** (*float or int*) – X coordinate in degrees visual angle for center of circle IA.
- **y** (*float or int*) – Y coordinate in degrees visual angle for center of circle IA.
- **radius** (*float or int*) – Radius of detection circle in degrees visual angle.

**Returns** List of whether saccade was detected and the gaze coordinates

at time of detection. In the form of [bool, [x, y]]. :rtype: bool, list

**send\_command** (*cmd*)

Sends a command to the Eyelink.

**Parameters** **cmd** (*str*) – Command to send.

**send\_message** (*txt*)

Sends a message to the tracker that is recorded in the EDF.

**Parameters** **txt** (*str*) – Message to send.

**send\_var** (*name, value*)

Sends a trial variable to the EDF file.

**Parameters**

- **name** (*str*) – Name of variable.
- **value** (*float, str, or int*) – Value of variable.

**set\_status** (*message*)

Sets status message to appear while recording.

**Parameters** **message** (*str*) – Text object to send, must be < 80 char

**set\_trialid** (*idval=1*)

Sends message that indicates start of trial in EDF.

**Parameters** **idval** – Values to set TRIALID.

**set\_trialresult** (*rval=0, scrcol=0*)

Sends trial result to indicate trial end in EDF and clears screen on EyeLink Display.

**Parameters**

- **rval** (*float, str, or int*) – Value to set for TRIAL\_RESULT.
- **scrcol** (*int*) – Color to clear screen to. Defaults to black.



---

## Sample Experiment

---

This is a very simple experiment that demonstrates use of `pylinkwrapper`. It can be found in the `sample` folder.

```
"""
Pylink Wrapper test experiment
N. DiQuattro - January 2015

This is a simple experiment where a circle appears randomly on the screen. It's
purpose is to provide examples of how to use the pylink wrapper with a psychopy
experiment.

There's help documentaiton available for each of the functions that show the
available parameters.
"""

# Import modules
from psychopy import visual
from psychopy import core, event
import numpy as np

import pylinkwrapper # Here's the special one

# Window set-up
win = visual.Window(monitor='nickMon', units='deg', fullscr=True,
                    allowGUI=False, color=0)

# Initiate eye-tracker link and open EDF
tracker = pylinkwrapper.Connect(win, '1_test')

# Calibrate eye-tracker
tracker.calibrate()

# Stimulus
fix = visual.Circle(win, radius=1, pos=(0, 0), fillColor=[1, 0, 0],
                    lineColor=[1, 0, 0])

cfix = visual.Circle(win, radius=.15, fillColor=-1, lineColor=-1)

# Display stimulus 5 times
for t in range(5):
    # Find random coordinates and set them
    fx = np.random.randint(-10, 10)
    fy = np.random.randint(-10, 10)
    fix.setPos((fx, fy))
```

```
# Eye tracker trial set-up
stxt = 'Trial %d' % t
tracker.set_status(stxt) # Define status message that appears on eye-link
# display

tracker.set_trialid() # Sends trial start message for EDF
tracker.send_message('Circle Trial')

# Draw IA
tracker.draw_ia(fx, fy, 1, 1, 5, 'circle') # Draw interest area and box

# Start recording
tracker.record_on()

# Draw and display circle
cfix.draw()
fix.draw()
win.flip()

# Wait for response
keyp = event.waitKeys()

# Stop Recording
tracker.record_off()

# Send response key to EDF file
tracker.send_var('response', keyp[0][0])

# End trial for EDF
tracker.set_trialresult()

# ISI with fixation check
cfix.draw()
win.flip()
tracker.fixcheck(2, 1, 'z')

# Retrieve EDF
tracker.end_experiment('C:\\\\edfs\\\\') # Closes and retrieves EDF file to
# specified path
```

---

## Indices and tables

---

- `genindex`
- `modindex`
- `search`



**C**

connector, 5



## C

calibrate() (connector.Connect method), 5  
Connect (class in connector), 5  
connector (module), 5  
convert\_coords() (connector.Connect method), 5

## D

draw\_ia() (connector.Connect method), 5  
draw\_text() (connector.Connect method), 6

## E

end\_experiment() (connector.Connect method), 6

## F

fix\_check() (connector.Connect method), 6

## G

get\_gaze() (connector.Connect method), 6

## R

record\_off() (connector.Connect method), 6  
record\_on() (connector.Connect method), 6

## S

sac\_detect() (connector.Connect method), 6  
send\_command() (connector.Connect method), 6  
send\_message() (connector.Connect method), 6  
send\_var() (connector.Connect method), 6  
set\_status() (connector.Connect method), 7  
set\_trialid() (connector.Connect method), 7  
set\_trialresult() (connector.Connect method), 7