
pyMedoc Documentation

Release 0.1.0

Cosan Lab

Feb 21, 2018

Contents

1	Installation	3
1.1	Usage	3
1.2	API Reference	4
2	API	7
3	Usage	9

Python package for wireless network triggering of Medoc Pathway Thermal Stimulation System


```
pip install git+https://github.com/cosanlab/pymedoc
```

1.1 Usage

1.1.1 Basic usage

To use this package, simply initialize the Pathway device and call its methods to interface with the Medoc system. pyMedoc automatically tries to establish a connection on initialization. Below are some example. See the API for a full list of commands.

```
from pymed.devices import Pathway
# Obviously this is a fake IP
my_pathway = Pathway(ip='10.10.10.10',port_number=9991)

# Check status
response = my_pathway.check_status()

# Start protocol 18
response = my_pathway.test_program(18)

# Stop a running program
response = my_path.stop()
```

1.1.2 Complete example of timed stimulation

This is a fully working script that ensures the execution of a program designed to deliver 10s of stimulation. The script assumes that program is numbered '100' on the Medoc system internally and configured for external control.

```
from pymedoc.devices import Pathway
import time

# Get settings from Medoc system
ip = 'X.X.X.X'
port = 20121

# We're going to trigger a protocol that stimulates for 10s
# Stimulation occurs at 46.5 deg but the system baseline is 32 deg
# The fastest possible ramp up is 10 deg/s so we need to account for
# ramp up and ramp down time before the trial ends for real
stimulation_time = 10
ramp_time = (46.5 - 32) / 10.
total_time = stimulation_time + ramp_time

# Establish initial connection test
medoc = Pathway(ip,port,verbose=False)

# Pause for a second
time.sleep(1)

# Issue the command to load the program
print("Starting program")
resp = medoc.program(100)

# The system has a variable length "pre-test" phase before a trigger can be sent.
# Any triggers sent during this time will be completely missed by the system.
# More annoying is that the length of the "pre-test" phase is unknowable ahead of time
# Here we use a convenient method built specifically for this situation
# It repeatedly checks the system 'test_state' until its value is 'RUNNING'
# where it can reliably get a trigger
# Print each poll attempt with the verbose flag
ready = medoc.poll_for_change('test_state', 'RUNNING', verbose=True)

# Trigger the start of the stimulation
print("Triggering")
resp = medoc.trigger()

# Wait the duration of stimulation
time.sleep(total_time)

# Send the stop signal
print("Stopping")
resp = medoc.stop()

# Again check until the system has ACTUALLY stopped to do anything else
ready = medoc.poll_for_change('test_state', 'IDLE', verbose=True)

print("END")
```

1.2 API Reference

This reference provides detailed documentation for all the features in pyMedoc

1.2.1 pymedoc.devices: pyMedoc Devices

class Pathway (*ip, port_number, timeout=5.0, verbose=True, buffer_size=1024*)

Pathway is a class to communicate with the Medoc Pathway thermal stimulation machine.

Parameters

- **ip** (*str*) – device ip address
- **port_number** (*int*) – port the device is listening on
- **timeout** (*float*) – seconds until connection timeouts; default 5s
- **verbose** (*bool*) – flag whether to print responses; default True
- **buffer_size** (*int*) – size of connection buffer; default 1024

call (*command, protocol=None, reuse_socket=False, verbose=False*)

Send command to device.

Parameters

- **command** (*str/int*) – command name or command_id number to send to device
- **protocol** (*str/int*) – protocol number on device to issue command to (only needed for command TEST_PROGRAM)
- **reuse_socket** (*bool*) – try to reuse the last created socket connection; *NOT CURRENTLY FUNCTIONAL*
- **verbose** (*bool*) – whether to print out the device callback

Returns response from Medoc system

Return type response (dict)

poll_for_change (*to_watch, desired_value, poll_interval=0.5, poll_max=-1, verbose=False, server_lag=1.0, reuse_socket=False*)

Poll system for a value change. Useful for waiting until the Medoc system has transitioned to a specific state in order to issue another command, but the transition length is unknowable.

Parameters

- **to_watch** (*str*) – the response field we should be monitoring; most often ‘test_state’ or ‘pathway_state’
- **desired_value** (*str*) – the desired value of the field to wait on, i.e. keep checking until response_field has this value
- **poll_interval** (*float*) – how often to poll; default .5s
- **poll_max** (*int*) – upper limit on polling attempts; default -1 (unlimited)
- **verbose** (*bool*) – print poll attempt number and current state
- **server_lag** (*float*) – sometimes if the socket connection is pinged too quickly after a value change the subsequent command after this method is called can get missed. This adds an additional layer of timing delay before returning from this method to prevent this; default 1s
- **reuse_socket** (*bool*) – try to reuse the last created socket connection; *NOT CURRENTLY FUNCTIONAL*

Returns whether desired_value was achieved

Return type *status* (bool)

status ()
Convenience method.

program (*protocol*)
Convenience method.

start ()
Convenience method.

pause ()
Convenience method.

trigger ()
Convenience method.

stop ()
Convenience method.

abort ()
Convenience method.

yes ()
Convenience method.

no ()
Convenience method.

CHAPTER 2

API

CHAPTER 3

Usage

A

[abort\(\)](#) (Pathway method), 6

C

[call\(\)](#) (Pathway method), 5

N

[no\(\)](#) (Pathway method), 6

P

[Pathway](#) (class in `pymedoc.devices`), 5

[pause\(\)](#) (Pathway method), 6

[poll_for_change\(\)](#) (Pathway method), 5

[program\(\)](#) (Pathway method), 6

S

[start\(\)](#) (Pathway method), 6

[status\(\)](#) (Pathway method), 5

[stop\(\)](#) (Pathway method), 6

T

[trigger\(\)](#) (Pathway method), 6

Y

[yes\(\)](#) (Pathway method), 6