

---

# **SILF Experiment API Documentation**

***Release 0.0***

**SILF Collaboration**

**Apr 18, 2018**

---

## Contents

---

<b>1</b>	<b>device Package</b>	<b>1</b>
<b>2</b>	<b>_const Module</b>	<b>2</b>
<b>3</b>	<b>_device Module</b>	<b>3</b>
<b>4</b>	<b>API fine print</b>	<b>6</b>
<b>5</b>	<b>Power management</b>	<b>8</b>
<b>6</b>	<b>Threading considerations</b>	<b>9</b>
<b>7</b>	<b>Change device state</b>	<b>10</b>
<b>8</b>	<b>How to test the devices according to the API</b>	<b>11</b>
8.1	Use DeviceWorkerWrapper from interpreter . . . . .	11
8.2	Auto result pooling . . . . .	12
<b>9</b>	<b>Device Api examples</b>	<b>13</b>
9.1	Engine driver . . . . .	13
9.2	Engine driver . . . . .	14
9.3	Engine and voltmeter connected . . . . .	14
	<b>Python Module Index</b>	<b>16</b>

# CHAPTER 1

---

## device Package

---

It packages api for single device.

## CHAPTER 2

---

### `_const` Module

---

```
silf.backend.common.device._const.DEVICE_STATES = ('off', 'stand-by', 'ready', 'running',  
    Tuple containing all allowable device states.
```

## CHAPTER 3

---

### \_device Module

---

This is a api for a device.

```
class silf.backend.common.device._device.Device (device_id='default',          con-  
                                              fig_file=None)
```

Bases: object

Defines plugin for a particular device in the experiment.

All methods are blocking, that is should block the current thread until finished.

---

**Note:** Instances of this object don't need to use any synchronization, they will always be called from single thread, This instance will be constructed used and destroyed on single process.

---

**Warning:** All methods should exit relatively fast.

**Warning:** Both method parameters and responses should be pickleable, these will be travelling between process boundaries.

**MAIN\_LOOP\_INTERVAL = 0.1**

Interval between invocations of main loop. Represents number of seconds as float.

**apply\_settings** (*settings*)

Applies some set of settings to this device.

**Parameters settings** (*dict*) – Settings to be applied, it is already validated by the IDeviceManager.

**Raises**

- **InvalidStateException** – If device is in invalid state (that is not STAND\_BY or RUNNING)

- **DeviceRuntimeException** – If any exception occurs

**Returns** None

**Return type** `None`

### **logger**

**Returns** Logger instance attached to this device. Utility method, you may use whatsoever logger you want

### **loop\_iteration()**

Perform an iteration of main experiment loop. Should terminate quickly,

**Raises** **DeviceRuntimeException** – If any exception occurs

**Returns** If returned value is *False* or *None* next iteration of this method will be scheduled after `MAIN_LOOP_INTERVAL` seconds, if result is true it will be sheduled earlier (after at most one command from controller was performed);

### **perform\_diagnostics** (*diagnostics\_level='short'*)

Performs diagnostics on the device. Can be ran if this device is OFF. or STAND\_BY.

**Parameters** **diagnostics\_level** (*str*) – Whether diagnostisc should be thororough or not, must be in `DEVICE_STATES`

**Raises**

- **DiagnosticsException** – If there is error in diagnostics.
- **InvalidStateException** – If device is in invalid state (that is not OFF)
- **DeviceRuntimeException** – If any exception occurs.

### **pop\_results()**

This method returns list of recently acquired points, it should clear this list so next calls won't return the same result points.

**Raises**

- **InvalidStateException** – If device is in invalid state (that is not RUNNING
- **InvalidStateException** – If device is in invalid state (that is not READY
- **DeviceRuntimeException** – If any exception occurs

**Returns** Returns results for (possibly) many points.

**Return type** `list` (or any other iterable) of *dict*.

### **post\_power\_up\_diagnostics** (*diagnostics\_level='short'*)

### **power\_down()**

Call to this method moves this class to OFF state.

**Raises**

- **InvalidStateException** – If device is in invalid state (that is not STAND\_BY
- **DeviceRuntimeException** – If any exception occurs

### **power\_up()**

Call to this method enables consecutive `apply_settings()`.

It also should power up the device (if this action makes any sense for this particular device see also: *Power management*).

**Raises**

- **InvalidStateException** – If device is in invalid state (that is not OFF)
- **DeviceRuntimeException** – If any exception occurs

**pre\_power\_up\_diagnostics** (*diagnostics\_level='short'*)

**start** ()

Starts the acquisition on the device (that is starts the measurements).

Blocks until this device is started.

**Raises** **InvalidStateException** – If device is in invalid state (that is not READY)

**Returns** None

**Return type** None

**state = None**

State of this device should be in DEVICE\_STATES, full state chart is available in: [Device state chart](#).

**stop** ()

Stops the acquisition on the device (that is stops the measurements).

Blocks until this device is started.

**Raises**

- **InvalidStateException** – If device is in invalid state (that is not RUNNING)
- **DeviceRuntimeException** – If any exception occurs

**Returns** None

**Return type** None

**tearDown** ()

**tear\_down** ()

Called when current process is being disabled.

This method can be called multiple times.

---

**Note:** do not override this method, override `_tear_down()`.

---

**Raises** **DeviceRuntimeException** – If any exception occurs

**Returns** None

**Return type** None

**exception** `silf.backend.common.device._device.InvalidCallToAssertState`

Bases: `Warning`

## CHAPTER 4

---

API fine print

---






Fig. 4.1: Device state chart

---

## Power management

---

---

**Note:** If your device does not need to power itself up or down, please just ignore `power_up()` and `power_down()` methods.

---

Devices should be powered up when we start call `power_up()`, but needn't do so, they must be powered up when after we exit from `start()`. So there are three methods in which devices should power up:

- `power_up()`, this method is called relatively early in during the experiment, and should allow plenty of time to initialize everything
- `apply_settings()`, use this method if your device powers up quickly.
- `start()`, if your device is volatile and you want to minimize the time it is powered up use this.

You can power down the device when following methods are called:

- `power_down()`
- `stop()`

## CHAPTER 6

---

### Threading considerations

---

Devices are accessed from single thread. All methods should exit relatively fast, **you should not use loops that are infinite** (or **can be infinite** — for example if hardware will not respond).

## CHAPTER 7

---

### Change device state

---

It is quite important to change state of your device after appropriate method calls.

---

## How to test the devices according to the API

---

There are two ways in which you can test it: start ipython interpreter create device and manage it by hand:

### 8.1 Use DeviceWorkerWrapper from interpreter

Import classes:

```
>>> from silf.backend.common_test.device.test_device import *
>>> from silf.backend.common.device_manager import start_worker_interactive
```

Start the device:

```
>>> work = start_worker_interactive('foo', MockDevice,
... configure_logging=False, auto_pull_results=False)

>>> work.state
'off'

>>> work.power_up()
UUID(...)
```

Let's setup the device:

```
>>> work.apply_settings({"foo": 3, "bar": 2})
UUID(...)

>>> work.start()
UUID(...)
```

This device will perform own acquisition in separate process, well wait for results to be acquired:

```
>>> time.sleep(1.2)
```

First `pop_results()` will return stale data, and schedule acquisition of new data:

```
>>> work.state
'running'
>>> work.pop_results() == []
True
```

Wait for results to get processed (will be faster on server!)

```
>>> time.sleep(0.5)
>>> results = work.pop_results()
>>> results == [{'foo_result': 3, 'bar_result': 2}]
True
```

Kill it without waiting;

```
>>> work.kill(wait_time=None)
```

## 8.2 Auto result pooling

You can configure this to auto poll for results:

```
>>> work = start_worker_interactive('foo', MockDevice,
... configure_logging=False, auto_pull_results=True)

>>> work.power_up()
UUID(...)
```

As in last test:

```
>>> work.apply_settings({"foo": 3, "bar": 2})
UUID(...)
>>> work.start()
UUID(...)
```

Wait for results to be gathered

```
>>> time.sleep(2)
```

Notice that results are available at once (no need to query)

```
>>> results = work.pop_results()
>>> results == [{'foo_result': 3, 'bar_result': 2}]
True

>>> work.kill(wait_time=None)
>>> results == [{'foo_result': 3, 'bar_result': 2}]
True

>>> work.kill(wait_time=None)
```

This is pseudocode

### 9.1 Engine driver

This imaginary device implements an engine. This is not actual experiment code, experiment will not be doing any waiting!

```
engine = ImaginaryDriver()

assert engine.state == 'off'

engine.power_up() # Powers up the device

assert engine.state == 'stand-by'

engine.apply_settings({"position" : 512})

assert engine.state == 'ready'

engine.start() # Start the engine

assert engine.state == 'acquiring'

# Silnik ruszył i teraz jest w stanie `acquiring`

# .. wait

while engine.state != 'ready':
    time.sleep(0.1)

# Silnik doszedł do końca i jest w stanie `ready`

# Następny punkt
```

```
engine.apply_settings({"position" : 1024})

engine.start() # Start the engine
```

## 9.2 Engine driver

Imaginary voltmeter

```
volt = ImaginaryVoltmeter()

assert volt.state == 'off'

volt.power_up() # Powers up the device

assert volt.state == 'stand-by'

volt.apply_settings({'range' : 15})

assert volt.state == 'ready'

volt.start() # Start the volt

assert volt.state == 'acquiring'

while volt.state != 'ready':
    time.sleep(0.1)

assert volt.pop_results() == [{'voltage' : 243.11}]
```

## 9.3 Engine and voltmeter connected

It works that so voltmeter measures single point after position is set by the engine.

```
engine = ImaginaryDriver()
volt = ImaginaryVoltmeter()

engine.power_up() # Powers up the device
volt.power_up() # Powers up the device

engine.apply_settings({"position" : 512})

engine.start();

while engine.state != 'ready':
    time.sleep(0.1)

volt.apply_settings({'range' : 15})

volt.start() # Start the volt

while volt.state != 'ready':
    time.sleep(0.1)
```



```
assert volt.pop_results() == [{'voltage' : 243.11}]

engine.apply_settings({"position" : 1024})

engine.start();

while engine.state != 'ready':
    time.sleep(0.1)

while volt.state != 'ready':
    time.sleep(0.1)

assert volt.pop_results() == [{'voltage' : 123.123}]
```

### S

`silf.backend.common.device`, [1](#)  
`silf.backend.common.device._const`, [2](#)  
`silf.backend.common.device._device`, [3](#)

## A

`apply_settings()` (`silf.backend.common.device._device.Device` method), 3

## D

`Device` (class in `silf.backend.common.device._device`), 3

`DEVICE_STATES` (in module `silf.backend.common.device._const`), 2

## I

`InvalidCallToAssertState`, 5

## L

`logger` (`silf.backend.common.device._device.Device` attribute), 4

`loop_iteration()` (`silf.backend.common.device._device.Device` method), 4

## M

`MAIN_LOOP_INTERVAL` (`silf.backend.common.device._device.Device` attribute), 3

## P

`perform_diagnostics()` (`silf.backend.common.device._device.Device` method), 4

`pop_results()` (`silf.backend.common.device._device.Device` method), 4

`post_power_up_diagnostics()` (`silf.backend.common.device._device.Device` method), 4

`power_down()` (`silf.backend.common.device._device.Device` method), 4

`power_up()` (`silf.backend.common.device._device.Device` method), 4

`pre_power_up_diagnostics()` (`silf.backend.common.device._device.Device` method), 5

## S

`silf.backend.common.device` (module), 1

`silf.backend.common.device._const` (module), 2

`silf.backend.common.device._device` (module), 3

`start()` (`silf.backend.common.device._device.Device` method), 5

`state` (`silf.backend.common.device._device.Device` attribute), 5

`stop()` (`silf.backend.common.device._device.Device` method), 5

## T

`tear_down()` (`silf.backend.common.device._device.Device` method), 5

`tearDown()` (`silf.backend.common.device._device.Device` method), 5