

---

# **pyCDB Documentation**

***Release 0.5.0-alpha***

**J. Urban et al.**

**Jul 28, 2017**



---

## Contents

---

<b>1</b>	<b>Description</b>	<b>3</b>
1.1	The data model . . . . .	3
1.2	Data acquisition management . . . . .	4
1.3	Database structure . . . . .	5
<b>2</b>	<b>Installation</b>	<b>7</b>
2.1	Dependencies . . . . .	7
2.2	Clone CDB repository . . . . .	7
2.3	pyCDB . . . . .	7
2.4	cyCDB . . . . .	8
2.5	Database . . . . .	8
<b>3</b>	<b>Configuration</b>	<b>9</b>
<b>4</b>	<b>Usage</b>	<b>11</b>
4.1	Reading data . . . . .	11
4.1.1	Basic example in Python . . . . .	11
4.1.2	Matlab . . . . .	11
4.1.3	Signal id's . . . . .	12
4.1.4	Linear signals with get_signal_data . . . . .	12
4.2	Writing data . . . . .	13
<b>5</b>	<b>CDB primer for (future) admins</b>	<b>15</b>
5.1	Populate SQL database . . . . .	15
5.2	Work with CDB . . . . .	16
<b>6</b>	<b>pyCDB reference</b>	<b>21</b>
6.1	pyCDB.client . . . . .	21
6.2	pyCDB.DAQClient . . . . .	31
6.3	pyCDB.pyCDBBase . . . . .	33
6.4	pyCDB.logbook . . . . .	36
6.5	pyCDB.CodeGeneration . . . . .	37
<b>7</b>	<b>Matlab CDB reference</b>	<b>39</b>
7.1	Matlab errors . . . . .	39
<b>8</b>	<b>Installation of FireSignal with pyCDB</b>	<b>41</b>

8.1	Installation steps . . . . .	41
8.1.1	INSTALLATION OF POSTGRES DB: . . . . .	41
8.1.2	HOW TO CREATE NEW CDB DATABASE? . . . . .	41
8.1.3	INSTALLATION OF FIRESIGNAL: . . . . .	42
8.1.4	TO RUN CENTRAL SERVER: . . . . .	42
8.1.5	INSTALLATION OF FIRESIGNAL GUI CLIENT: . . . . .	42
8.1.6	HDF5 plugin: . . . . .	42
8.2	Nota Bene: How many instances of FSs are running? . . . . .	43
<b>9</b>	<b>Compilation of FireSignal:</b>	<b>45</b>
9.1	Example: LAST RECORD NUMBER IN GUI CLIENT . . . . .	45
<b>10</b>	<b>JyCDB</b>	<b>47</b>
10.1	Javadoc . . . . .	47
10.1.1	cz.cas.ipp.compass.jycdb . . . . .	47
10.1.2	cz.cas.ipp.compass.jycdb.plotting . . . . .	67
10.1.3	cz.cas.ipp.compass.jycdb.test . . . . .	68
10.1.4	cz.cas.ipp.compass.jycdb.util . . . . .	68
<b>11</b>	<b>Indices and tables</b>	<b>85</b>
	<b>Python Module Index</b>	<b>87</b>

Contents:



The Compass DataBase (CDB) is a lightweight system designed for storing the COMPASS (IPP Prague) tokamak experimental data. It can equally well be used for any tokamak or generally any device that repeatedly produces experimental data.

CDB uses HDF5 (or NetCDF 4) files to store numerical data and a relational database, actually MySQL, to store metadata. The core application is implemented in Python (pyCDB). Cython is used to wrap the Python code in a C API. Matlab, IDL etc. clients can then be built using the C API.

There are several major advantages of this scheme:

- Vast of the required functionality is readily implemented and available for numerous operating systems and applications (high/low level data input/output, database functionality).
- Data can be stored on any file system (local or remote), no need for specific protocol.
- Rapid, platform independent development in Python.

CDB has a possibility to store the information about the data acquisition sources (channels) of the data. The database contains information about DAQ channels associations to physical quantities.

An important point in CDB is its *never overwrite* design. Anything stored in CDB cannot be overwritten (at least using the standard API); instead, revisions are possible as corrective actions.

## The data model

A relational database is used to store metadata of the numerical data. Metadata include physical quantities names, units, information about axes (which themselves are physical quantities and are the same entities as any other quantities).

**generic\_signals** Describe physical quantities stored in the database. In particular, contain names, units, axes id's (axis are treated as any other signals), description, signal type (FILE or LINEAR), record numbers validity range, and data source.

**data\_sources** Used as primary grouping criterion. Contain name, description and the directory name of the data files.

**data\_files** List of all data files in the system. Each file can contain one or more signals. Data files are stored in subdirectories specified in the data source under the main CDB directory. **files\_status** states whether a file is ready for reading.

**data\_signals** In fact, data signals are instances of generic signals. A data signal either points to data in a data file or contains only coefficients for linear function (LINEAR signal). Data signal contains record number to which the data belong. Revisions can be created when a correction to a signal is needed.

Each data signal contains *offset* and *coefficient* used either for a linear signal construction or for linear transformation of data stored in a file. See [Linear signals with get\\_signal\\_data](#) for details. *time0* specifies the time of the first data point for time-dependent signals.

**shot\_database** Contains record numbers—unique numbers characterizing a data set. This is mostly a tokamak shot, can however be a simulation, a DAQ system test etc. Tokamak shots also have **shot\_numbers**. Data files for a particular shot are stored in **record\_directories**.

**FireSignal\_Event\_IDs** CDB can be used as storage system for FireSignal. In this case, this table relates FireSignal id's and CDB record numbers.

## Data acquisition management

CDB has a possibility to track information about DAQ A/D channels. Each data acquisition system can be associated (“attached”) with a physical quantity (generic\_signal) it outputs.


**DAQ\_channels** List of all DAQ channels available. Unique identification consist of *computer\_id*, *board\_id*, *channel\_id*. In case CDB is used with FireSignal, *nodeuniqueid*, *hardwareuniqueid* and *parameteruniqueid* relates the two databases.

Each DAQ channel has a *default\_generic\_signal\_id*, which is the id of a (unique to the channel) generic signal associated with the channel if no other generic signal is associated. The reason for this is that, in CDB, every **data\_signal** must have a *generic\_signal\_id*.

**channel\_setup** This table tells to which *generic\_signal\_id* is a DAQ channel associated (attached). It's an event-style table containing association events (date, time, user id and a note of a DAQ channel association to a generic signal.)



## Database structure





### Dependencies

For pyCDB

- Python 2.7+
- MySQLdb or pymysql
- numpy
- h5py
- pint
- matplotlib (for the example only)

For JyCDB - Java, Matlab, IDL clients

- Jython

For cyCDB (the C-interface)

- Cython (tested with versions 0.14 and 0.15)

### Clone CDB repository

CDB's Mercurial repository is hosted at [bitbucket.org/compass-tokamak/cdb](https://bitbucket.org/compass-tokamak/cdb). Simply clone it via

```
hg clone https://bitbucket.org/compass-tokamak/cdb
```

### pyCDB

Just place pyCDB on Python path and `import pyCDB`

If pyCDB is not on Python path, it is possible to modify the path:

```
export CDB_PATH=/path/to/cdb/src && export PYTHONPATH=$PYTHONPATH:/path/to/cdb/src
```

## cyCDB

Run `make` in the `src` directory. This step is optional.

## Database

CDB needs MySQL database. Import *database/mysql\_setup.sql* to create the *CDB* database structure, run

```
mysql -u root -p < database/CDB_schema.sql
```

To create CDB users and the testing database (modify passwords for production use), run

```
mysql -u root -p < CDB_users.sql
```

## CHAPTER 3

---

### Configuration

---

CDB configuration parameters are sought in:

1. environment variables
2. `./CDBrc`
3. `$HOME/.CDBrc`

CDBrc file explanation, environment variables have the same names:

```
# CDB configuration file
#
# Edit and save to .CDBrc
#
# Possible directories, ordered by ascending priority
# - $CDB_PATH/..
# - user home directory
# - current directory
# - alternatively, CDB can be configured with environment variables

# Database configuration
[database]
# Root data directory
CDB_DATA_ROOT = /var/local/CDB

# MySQL host name or IP (:port)
CDB_HOST = localhost
# MySQL user name
CDB_USER = CDB_test
# MySQL password
CDB_PASSWORD = cdb_data
# MySQL database name
CDB_DB = CDB_test

# Logging configuration
[logging]
```

```
# Screen logging level: INFO, DEBUG, WARNING, ERROR, NOTSET
CDB_LOG_LEVEL = WARNING
# File log file name
CDB_LOG_FILE =
# File log level
CDB_FILE_LOG_LEVEL = INFO

# logbook configuration (optional)
[logbook]
CDB_LOGBOOK_HOST = localhost
CDB_LOGBOOK_USER = CDB
CDB_LOGBOOK_PASSWORD = cdb_data
CDB_LOGBOOK_DB = logbook
```

## Reading data

### Basic example in Python

First import pyCDB and connect to the database:

```
from pyCDB.client import CDBClient
cdb = CDBClient()
```

Retrieve references for some generic signal:

```
generic_signal_name = 'electron density'
# get the full generic signal reference (all columns)
generic_signal_refs = cdb.get_generic_signal_references(generic_signal_name = generic_
↳ signal_name)
# get signal references
signal_refs = cdb.get_signal_references(record_number=-1, generic_signal_id=generic_
↳ signal_refs[0]['generic_signal_id'])
```

Now get the signal data, including description and axes, using `pyCDB.client.CDBClient.get_signal()`:

```
sig = cdb.get_signal(signal_ref=signal_refs[-1])
```

## Matlab

Use `cdb_client` to instantiate a client:

```
cdb = cdb_client();
```

To get the signal data including axes use `cdb_client.cdb_get_signal()`

```
signal = cdb.get_signal(str_id)
```

## Signal id's

*Generic signal* is uniquely identified either by

- numeric id (`generic_signal_id`)
- alias + record number
- name (`generic_signal_name`) + data source id + record number

CDB interface supports generic signal string id's in the following forms:

- *alias\_or\_name* - search by alias first, if no match is found search by name (>1 results possible)
- *name/data\_source\_id* - generic signal name followed by '/' and `data_source_id` (data source name or numeric id)
- *generic\_signal\_id* - numeric id of the generic signal

See also `pyCDB.client.decode_generic_signal_strid()`.

*Data signal* (single data set) is uniquely identified by

- generic signal id (numeric) + record number + revision

Data signal string id is in the form of:

```
id_type:generic_signal_string_id:record_number:revision[units]
```

For example, this refers to EFIT psi 2D, shot 4073, last revision:

```
CDB:psi_2D:4073:-1[default]
```

which is equivalent to:

```
psi_2D:4073
```

See also `pyCDB.client.decode_signal_strid()`.

## Linear signals with get\_signal\_data

Described below is the logic for linear signals, implemented in `pyCDB.client.CDBClient.get_signal_data()`.

- number of axes > 1 → **unresolved**
- number of axes = 1
  - axis is linear
    - \* axis is time\_axis
      - time\_limit provided: `n_samples = int(ceil((time_limit - offset) / coefficient))`
      - n\_samples not provided → **unresolved**
      - **result** = `offset + coefficient*(axis_data[0..n_samples] - time0)`
    - \* axis is not time\_axis



- `n_samples` not provided → **unresolved**
- **result** = `offset + coefficient*axis_data[0..n_samples]`
- axis contains data
  - \* axis is `time_axis`
    - **result** = `offset + coefficient*(axis_data - time0)`
  - \* axis is not `time_axis`
    - **result** = `offset + coefficient*axis_data`
- number of axes = 0
  - `n_samples` and `x0` provided (`x0` is optional, defaults to 0)
    - \* **result** = `x0 + [0,1, ..., n_samples-1]*coefficient`
  - `x0` and `x1` provided
    - \* **result** = `[x0, x0+coefficient, x0+2*coefficient, ..., x1]`
  - otherwise
    - \* **result = function:** `f(i) = offset + coefficient * i`

## Writing data

First create a new data file record (go to the last step for liner signals), providing data source id, record number and collection (a base name for the data file chosen by the user):

```
file_ref = cdb.new_data_file(collection_name, data_source_id =
    data_source_id, record_number = record_number, \
    file_format = "HDF5")
```

Now you can create the file and fill with data, e.g.:

```
import h5py
fh5 = h5py.File(file_ref['full_path'], 'w')
grp_name = 'raw data'
data_file_key = grp_name + '/' + generic_signal_name
f_grp = fh5.create_group(grp_name)
f_grp.create_dataset(generic_signal_name, data=numpy_data)
fh5.close()
```

One has to say when the file is ready (for reading):

```
cdb.set_file_ready(file_ref['data_file_id'])
```

Finally, the CDB database must know what signals are stored in the file. Linear signals are created solely by this step:

```
cdb.store_signal(generic_signal_id, \
    record_number=record_number, \
    data_file_id=file_ref['data_file_id'], \
    data_file_key = data_file_key, \
    offset=0.5, coefficient=1.3, time0=-0.2, \
    computer_id=1, board_id=1, channel_id=1, \
    note='my first stored signal')
```



---

## CDB primer for (future) admins

---

This tutorial demonstrates some basics of CDB when starting from an empty database. This is useful for testing a new installation, particularly by (future) admins and developers.

```
In [1]: %matplotlib inline
In [2]: import numpy as np
        import matplotlib.pyplot as plt
```

Import CDB and connect

```
In [3]: from pyCDB.client import CDBClient
        cdb = CDBClient()
```

## Populate SQL database

For simplicity, we will populate the database for this tutorial by SQL directly. Login to MySQL / MariaDB with enough permissions (e.g. as root) and execute:

```
INSERT INTO `da_computers` (`computer_id`, `computer_name`, `location`,
→ `description`) VALUES
(1, 'daq_comp_1', NULL, NULL);

INSERT INTO `data_sources` (`data_source_id`, `name`, `description`, `subdirectory`)
→ VALUES
(1, 'MAGNETICS', NULL, 'magnetics'),
(2, 'EFIT', NULL, 'efit'),
(3, 'DAQ', NULL, 'daq');

INSERT INTO `generic_signals` (`generic_signal_id`, `generic_signal_name`, `alias`,
→ `first_record_number`, `last_record_number`, `data_source_id`, `time_axis_id`,
→ `axis1_id`, `axis2_id`, `axis3_id`, `axis4_id`, `axis5_id`, `axis6_id`, `units`,
→ `description`, `signal_type`) VALUES
(1, 'daq_time', NULL, 1, -1, 3, NULL, NULL, NULL, NULL, NULL, NULL, NULL, 's', NULL,
→ 'LINEAR'),
```

```

(2, 'daq_comp_1_channel_1_1', NULL, 1, -1, 3, 1, NULL, NULL, NULL, NULL, NULL, NULL,
→ NULL, NULL, 'FILE'),
(3, 'I_plasma', 'Ip', 1, -1, 1, 1, NULL, NULL, NULL, NULL, NULL, NULL, 'A', NULL,
→ 'FILE'),
(4, 't_efit', NULL, 1, -1, 2, NULL, NULL, NULL, NULL, NULL, NULL, NULL, 's', NULL,
→ 'FILE'),
(5, 'I_plasma', NULL, 1, -1, 2, 4, NULL, NULL, NULL, NULL, NULL, NULL, 'A', NULL,
→ 'FILE');

INSERT INTO `DAQ_channels` (`computer_id`, `board_id`, `channel_id`, `default_generic_
→ signal_id`, `note`,
`nodeuniqueid`, `hardwareuniqueid`, `parameteruniqueid`) VALUES
(1, 1, 1, 2, NULL, NULL, NULL, NULL);

INSERT INTO `shot_database` (`record_number`, `record_time`, `record_type`,
→ `description`) VALUES
('1', 'NOW()', 'EXP', NULL);

INSERT INTO `record_directories` (`record_number`, `data_directory`) VALUES
('1', '1');

```

## Work with CDB

## List all generic signals

```
In [4]: gs_list = cdb.get_generic_signal_references()
        for gs in gs_list:
            ds = cdb.get_data_source_references(data_source_id=gs.data_source_id)[0]
            print("id: {generic_signal_id}, name/source: {generic_signal_name}/{data_source_name}".format(
                **gs, data_source_name=ds.name))

id: 1, name/source: daq_time/DAQ
id: 2, name/source: daq_comp_1_channel_1_1/DAQ
id: 3, name/source: I_plasma/MAGNETICS
id: 4, name/source: t_efit/EFIT
id: 5, name/source: I_plasma/EFIT

In [5]: cdb.get_data_source_references(data_source_id=gs.data_source_id)

Out[5]: [OrderedDict([('data_source_id', 2),
                      ('name', 'EFIT'),
                      ('description', None),
                      ('subdirectory', 'efit')])]
```

## Populate our first shot with data

```
In [6]: shot = cdb.last_record_number()
        print(shot)

1

In [7]: # example data + noise
        n = 100000
        dt = 1e-3
        t_data = np.r_[0:n*dt:dt]
        Ip_data = 1e6 * (1 - np.linspace(-1, 1, n)**8)**0.5
        Ip_data_noisy = Ip_data * (1 + (0.05*np.random.rand(n) - 0.025))
```

Use `put_signal` to store data. `put_signal` <[http://cdb.readthedocs.io/en/latest/reference.html?highlight=put\\_signal#pyCDB.client.CDBClient.put\\_signal](http://cdb.readthedocs.io/en/latest/reference.html?highlight=put_signal#pyCDB.client.CDBClient.put_signal)> is a convenience wrapper around the full (and more flexible) store signal procedure, described in <http://cdb.readthedocs.io/en/latest/usage.html#writing-data>.

```
In [8]: gs_id = 3
        cdb.put_signal(gs_id, shot, Ip_data_noisy, time0=0, time_axis_coef=dt)

Out[8]: OrderedDict([('record_number', 1),
                    ('generic_signal_id', 3),
                    ('revision', 1),
                    ('variant', ''),
                    ('timestamp', datetime.datetime(2017, 7, 28, 11, 50, 24)),
                    ('data_file_id', 1),
                    ('data_file_key', 'I_plasma'),
                    ('time0', 0.0),
                    ('coefficient', 1.0),
                    ('offset', 0.0),
                    ('coefficient_V2unit', 1.0),
                    ('time_axis_id', 1),
                    ('time_axis_revision', 1),
                    ('axis1_revision', 1),
                    ('axis2_revision', 1),
                    ('axis3_revision', 1),
                    ('axis4_revision', 1),
                    ('axis5_revision', 1),
                    ('axis6_revision', 1),
                    ('time_axis_variant', ''),
                    ('axis1_variant', ''),
                    ('axis2_variant', ''),
                    ('axis3_variant', ''),
                    ('axis4_variant', ''),
                    ('axis5_variant', ''),
                    ('axis6_variant', ''),
                    ('note', ''),
                    ('computer_id', None),
                    ('board_id', None),
                    ('channel_id', None),
                    ('data_quality', 'UNKNOWN'),
                    ('deleted', 0)])
```

```
In [9]: cdb.put_signal(5, shot, Ip_data, time_axis_data=t_data)

Out[9]: OrderedDict([('record_number', 1),
                    ('generic_signal_id', 5),
                    ('revision', 1),
                    ('variant', ''),
                    ('timestamp', datetime.datetime(2017, 7, 28, 11, 50, 32)),
                    ('data_file_id', 2),
                    ('data_file_key', 'I_plasma'),
                    ('time0', 0.0),
                    ('coefficient', 1.0),
                    ('offset', 0.0),
                    ('coefficient_V2unit', 1.0),
                    ('time_axis_id', 4),
                    ('time_axis_revision', 1),
                    ('axis1_revision', 1),
                    ('axis2_revision', 1),
                    ('axis3_revision', 1),
                    ('axis4_revision', 1),
                    ('axis5_revision', 1),
```

```
('axis6_revision', 1),
('time_axis_variant', ''),
('axis1_variant', ''),
('axis2_variant', ''),
('axis3_variant', ''),
('axis4_variant', ''),
('axis5_variant', ''),
('axis6_variant', ''),
('note', ''),
('computer_id', None),
('board_id', None),
('channel_id', None),
('data_quality', 'UNKNOWN'),
('deleted', 0)])
```

Get the data back from CDB using the *string id*.

```
In [10]: str_id = 'Ip:{}'.format(shot)
        str_id
```

```
Out[10]: 'Ip:1'
```

```
In [11]: ip_sig = cdb.get_signal(str_id)
        ip_sig
```

```
Out[11]: { Generic signal name: I_plasma
          Generic signal alias: Ip
          Generic signal id: 3
          Record number: 1
          Revision: 1
          Units: A
          Data: shape = (100000,), dtype = float64
          More details in: ref, gs_ref, daq_attachment, file_ref}
```


Now get the data using keyword arguments.

```
In [12]: ip_efit = cdb.get_signal(generic_signal_id=5, shot=shot)
        ip_efit
```

```
Out[12]: { Generic signal name: I_plasma
          Generic signal alias: None
          Generic signal id: 5
          Record number: 1
          Revision: 1
          Units: A
          Data: shape = (100000,), dtype = float64
          More details in: ref, gs_ref, daq_attachment, file_ref}
```

```
In [13]: fig, ax = plt.subplots(figsize=(14, 8))
        ax.plot(ip_sig.time_axis.data, ip_sig.data, label="magnetics")
        ax.plot(ip_efit.time_axis.data, ip_efit.data, c='r', lw=2, label="EFIT")
        ax.legend(loc='best')
```

```
Out[13]: <matplotlib.legend.Legend at 0x7fed7ace5128>
```



In [ ]:





## pyCDB.client

COMPASS database main module

Contains functions to connect, retrieve and store data using the COMPASS database.

```
class pyCDB.client.CDBClient (host=None, user=None, passwd=None, db=None, port=3306,  
                             log_level=None, data_root=None)
```

Fundamental connector class for CDB

```
CDB2FS_ref (computer_id, board_id, channel_id)
```

Return the CDB DAQ\_channel reference of a CDB channel

```
FS2CDB_ref (nodeuniqueid, hardwareuniqueid, parameteruniqueid)
```

Return the CDB DAQ\_channel reference of a FireSignal channel

```
add_unit (unit_system_id, unit_name, auto_dim=False, **dimensions)
```

Add unit to unit system.

**Parameters** **auto\_dim** – if true, dimensions are inferred automatically.

**Returns** unit\_id of the unit

**Return type** int

```
apply_unit_factor_tree (signal, unit_tree)
```

Apply unit conversion of signal data tree.

**Parameters**

- **signal** – a CDBSignal object with dependent axes and fully read data.
- **unit\_tree** – a tree of units to convert to as returned by get\_unit\_factor\_tree.

```
create_generic_signal (generic_signal_name=None, alias=None, first_record_number=1,
                        last_record_number=-1, data_source_id=None, time_axis_id=None,
                        axis1_id=None, axis2_id=None, axis3_id=None, axis4_id=None,
                        axis5_id=None, axis6_id=None, units=None, description=None,
                        signal_type='FILE', **kwargs)
```

Creates a generic signal in database.

Returns the id of the new generic signal.

```
create_record (record_type='EXP', record_time=None, description='', data_root=None,
                FS_event_number=0, FS_event_id='', record_number=None, autonumber=False)
```

Creates a new record in the database

#### Parameters

- **record\_type** – 'EXP', 'VOID' or 'MODEL'
- **record\_time** – default current, can be integer time stamp or a datetime.datetime object
- **description** – record description
- **data\_root** – CDB data root path, see [get\\_data\\_root\\_path\(\)](#)
- **FS\_event\_number** – FireSignal event number (optional) # TODO: Remove?
- **FS\_event\_id** – FireSignal event id (optional) # TODO: Remove?
- **record\_number** – new record number (None for first available)
- **autonumber** – automatic numbering (effective only if record\_number is not None)

**Return type** new record number

```
delete_signal (str_id='', **kwargs)
```

Delete a signal

#### Parameters

- **str\_id** – string id of the signal
- **kwargs** – alternative signal id (signal\_reference)

kwargs - fields from data signal reference:

#### Parameters

- **timestamp** –
- **note** –

```
find_generic_signals (alias_or_name='', regexp=False)
```

Find generic signals by a name or alias

#### Parameters

- **alias\_or\_name** – search string
- **regexp** – use regular expressions (MySQL REGEXP)

**Return type** tuple of signal references (like [get\\_generic\\_signal\\_references](#))

The search criterion in MySQL is “SELECT ... WHERE generic\_signal\_name LIKE %%alias\_or\_name%% OR alias LIKE %%alias\_or\_name%%”

```
get_FS_signal_reference (nodeuniqueid, hardwareuniqueid, parameteruniqueid)
```

Get the currently attached generic signal reference for a FireSignal channel.

Returns None if not found.

**get\_attachment\_table** (*timestamp=None, computer\_id=None, board\_id=None, channel\_id=None, generic\_signal\_id=None, record\_number=None*)  
Returns tuple of rows from a “view” corresponding to channel attachments based on specified criteria.

**Parameters**

- **timestamp** – reference time (default current)
- **computer\_id** – DAQ computer id
- **board\_id** – DAQ board id
- **channel\_id** – DAQ channel id
- **generic\_signal\_id** – attached generic signal id
- **record\_number** – timestamp of this record number will be used

**get\_channel\_references** (*computer\_name=None, computer\_id=None, board\_id=None, channel\_id=None*)  
Get DAQ channel reference(s)

**get\_coefficient\_V2unit** (*generic\_signal\_id*)  
Get coefficient\_V2unit (needed for FireSignal operation)

**get\_coefficient\_lev2V** (*generic\_signal\_id*)  
Get coefficient\_lev2V (needed for FireSignal operation)

**get\_computer\_id** (*computer\_name=None, description=None*)  
Returns computer ID by computer name or description

**get\_data\_file\_reference** (*\*\*kwargs*)  
Returns a tuple of dictionaries containing complete file references containing the signal. Returns one reference for each file revision.

Typical parameters: data\_file\_id: unique id generic\_signal\_name, data\_source\_id (name), data\_source\_id, generic\_signal\_id record\_number

If data\_source\_id is not specified the signal is first sought in generic signal names.

**get\_data\_root\_path** (*data\_root=None*)  
Return the CDB data root path

**Parameters** **data\_root** – user specified root directory, if None given by CDB\_DATA\_ROOT

**Return type** the full path with expanded variables

**get\_data\_source\_id** (*data\_source\_name*)  
Returns data source id of a given data source name

**get\_data\_source\_references** (*\*\*kwargs*)  
Returns a tuple with references to data sources

Typical parameters: data\_source\_name, data\_source\_id, description, record\_number

**get\_data\_source\_subdir** (*data\_source\_id*)  
Returns the data storage subdirectory for a given data source id

**get\_generic\_signal\_id** (*generic\_signal\_name, data\_source\_id=None*)  
Returns generic signal id given the signal name

**get\_generic\_signal\_references** (*str\_id='', \*\*kwargs*)  
Returns a tuple of references to generic signal, based on specified criteria

**Parameters** **str\_id** – generic signal string id – see `decode_generic_signal_strid()`

Keyword parameters: `generic_signal_id` (or `signal_id`): id of the generic signal `generic_signal_ref` or `generic_signal_reference` `generic_signal_name` + (optionally) `data_source_id` or `data_source_name` + `record_number` (for validity) `alias_or_name`: first find by alias; if no match find by name `generic_signal_strid`: deprecated - use `str_id`

**classmethod** `get_max_revisions` (*signals*, *deleted=False*)

From a set of signals, take only those that have maximum revision number for a combination of (`record_number`, `generic_signal_id`).

**Parameters** `deleted` – If true, keep deleted, otherwise filter them out.

**get\_record\_datetime** (*\*\*kwargs*)

Return record datetime (timestamp)

`kwargs`: :param `record_number`: record number

**get\_record\_path** (*\*\*kwargs*)

Returns data storage path (from the database) for a given record number

**get\_signal** (*str\_id=None*, *squeeze=True*, *deleted=False*, *dtype=None*, *\*\*kwargs*)

Get signal information and data.

**Parameters**

- **str\_id** – signal string id, see `decode_signal_strid`
- **squeeze** – remove singleton dimensions in the received data
- **deleted** – return or not deleted signals (this will not work anyway ...)
- **dtype** – force a certain numpy data type

Keyword parameters: :param `units`: raw / dav / default (case insensitive, conflict with `str_id`) :param `signal_ref`: signal reference :param `generic_signal_ref`: generic signal reference :param `n_samples`: number of samples :param `time_limit`: time limit (instead of `n_samples`) :param **slice**: standard Python slice :rtype: CDBSignal instance

If `str_id` argument is not present, keyword arguments are used in the search.

**get\_signal\_base\_tree** (*signal\_ref*)

Create a reference tree of dependent signals.

**Parameters** **signal\_ref** – A valid signal reference

It reads parameters from the database and sets some signal properties as well.

**get\_signal\_calibration** (*str\_id=None*, *sig\_ref=None*, *units=None*, *gs\_ref=None*)

Calculate signal offset and coefficient based on signal reference

**Parameters**

- **str\_id** – signal string id
- **sig\_ref** – signal reference
- **units** – signal units, None uses default

The calibration transformations is: (data + offset) \* coefficient

If `str_id` is input, all the references are fetched from the database. Note: Applies only to FILE signals, LINEAR cannot be transformed this way

**get\_signal\_data** (*signal\_ref=None*, *generic\_signal\_ref=None*, *signal=None*, *squeeze=True*, *raw=False*, *\*\*kwargs*)

Get signal data.

**Parameters** **squeeze** – Reduce  $N \times 1 \times \dots \times M$  data\_sets into  $N \times \dots \times M$  (numpy squeeze function).

For FILE data signals, it returns the data without any modifications (apart from squeeze or slicing). For LINEAR data signals, it returns: - if `n_samples` is set: an array of linear values - else: a linear function

kwargs:

**Parameters** **x0** – if set, it is used instead of `signal_ref[‘offset’]`

Parameters: `time_limit`, `n_samples` `x0`, `x1` - axis limits raw - if True, data won’t be scaled

**get\_signal\_data\_tree** (*signal*, *units*=‘default’, *n\_samples*=None, *squeeze*=True, *x0*=None, *slice*=Ellipsis, *is\_this\_axis*=False, *dtype*=None)

Recursively read data for signal and all its axes.

**Parameters**

- **signal** – the signal object to set data to.
- **x0** – additive constant

**get\_signal\_parameters** (*str\_id*=‘’, *parse\_json*=True, *\*\*kwargs*)

Get data signal parameters from data\_signal\_parameters table

**Parameters**

- **str\_id** – string id
- **kwargs** – keyword arguments for signal identification
- **parse\_json** – if True, the method tries to parse fields as JSON.

**get\_signal\_references** (*str\_id*=‘’, *deleted*=False, *limit*=10000, *dereference*=False, *\*\*kwargs*)

Get signal references by specified criteria

**Parameters**

- **str\_id** – signal string id (see `decode_signal_strid()`)
- **deleted** – return or not deleted references for revision=-1
- **limit** – limit the number of results
- **dereference** – return explicit `gs_ref`, `data_file_ref` and `data_source_ref`

keyword arguments offer an alternative to `str_id`: :param `record_number`: record number :param `revision`: revision :param `signal_ref`: signal reference :param `signal_reference`: the same as `signal_ref`

**get\_signal\_setup** (*gs\_str\_id*, *record\_number*=None, *timestamp*=None)

Get signal setup for a generic signal

**Parameters** **gs\_str\_id** – generic signal alias, name or numeric id

**get\_unit** (*\*\*kwargs*)

Get a unit.

kwargs:

**Parameters**

- **unit\_id** – Direct ID of the unit.
- **(or unit\_system\_id or unit\_system\_name)** (*unit\_system*) – The unit system to use.
- **construct** – whether to construct unit from basic units (default: True)

For unit system, specify any of the following keyword arguments: m, kg, s, A, K, mol, cd (as dimensionality)

You can specify either `unit_id` OR `unit_system` (`unit_system_id/unit_system_name`) and dimensionality.

**get\_unit\_factor\_tree** (*signal*, *units=None*, *unit\_system=None*)

Get tree of factors and unit names for a signal object.

#### Parameters

- **signal** – CDBSignal + its axes signals in a tree
- **units** – name of the unit or unit system to convert to
- **unit\_system** – unit system to convert to (takes precedence)

Works recursively but only unit system is passed along the recursion. Example return value: { 'units': 'mm', 'factor' : 0.001 }

**get\_unit\_system** (*unit\_system\_name=None*, *unit\_system\_id=None*, *load\_parent=True*, *\*\*ignored\_args*)

Get a unit system.

#### Parameters

- **unit\_system\_name** – Name of the system.
- **unit\_system\_id** – The numerical ID of the system.
- **load\_parent** – if true, automatically load parent unit systems.

You have to specify either the system name (takes precedence) or `unit_system_id`.

**is\_file\_ready** (*data\_file\_id*)

Set file\_ready to True

**last\_record\_number** (*record\_type='EXP'*)

Returns the “record number” of the last available record

**Parameters** **record\_type** – record type: EXP, VOID, MODEL or empty string (any type)

**Return type** record number (integer)

**last\_shot\_number** ()

The same as `last_record_number` for `record_type = 'EXP'`

**Return type** record number (integer)

**new\_data\_file** (*collection\_name*, *file\_format='HDF5'*, *file\_extension=None*, *create\_subdir=None*, *data\_root=None*, *\*\*kwargs*)

Creates a record in the COMPASS database for a new data file. Returns the created data file reference.

#### Parameters

- **collection\_name** – signal collection name (used for file name)
- **create\_subdir** – automatically create the file subdirectory (default if `data_root` is not set)
- **data\_root** – use user-supplied `data_root` instead of the default one
- **kwargs** – keyword parameters \* `record_number`: `resolve_record_number` \* `data_source_id`, `data_source_ref`, or `data_source` (name)

First looks for a record with the input record number, data source and collection. If found, a new revision is created. Otherwise a new id with revision 1 is created.

```
put_signal(gs_alias_or_id, record_number, data, variant='', time0=None, data_coef=None,
            data_offset=None, data_coef_V2unit=None, note=None, no_debug=None,
            data_quality='UNKNOWN', time_axis_data=None, time_axis_coef=None,
            time_axis_offset=None, time_axis_note=None, time_axis_id=None, axis1_data=None,
            axis1_coef=None, axis1_offset=None, axis1_note=None, axis2_data=None,
            axis2_coef=None, axis2_offset=None, axis2_note=None, axis3_data=None,
            axis3_coef=None, axis3_offset=None, axis3_note=None, axis4_data=None,
            axis4_coef=None, axis4_offset=None, axis4_note=None, axis5_data=None,
            axis5_coef=None, axis5_offset=None, axis5_note=None, axis6_data=None,
            axis6_coef=None, axis6_offset=None, axis6_note=None, gs_parameters=None,
            daq_parameters=None)
```

Write CDB signal to HDF5 and database

`cdb_signal = cdb_put_signal(gs_alias_or_id, record_number, data, ...)` writes signal data into CDB, including HDF5 file.

#### Parameters

- **gs\_parameters** (*JSON string or dict*) – generic signal parameters
- **daq\_parameters** (*JSON string or dict*) – data acquisition channel parameters

```
record_exists(record_number=None, **kwargs)
```

Returns True if a record exists, False otherwise.

**Parameters** **kwargs** – named parameters interpretable by `resolve_record_number` (see)

**Return type** bool

```
resolve_record_number(argv)
```

Resolves record number from input parameters dictionary `argv`

**Parameters** **argv** – dictionary with parameters

**Return type** record number

Recognized `argv` keys:

- **record\_number** (-1 returns last record number)
- **record\_type**

```
set_file_ready(data_file_id, chmod_ro=True)
```

Set `file_ready` to True

#### Parameters

- **data\_file\_id** – data file id
- **chmod\_ro** – make the file read-only (strongly recommended)

```
classmethod signal_dim(signal_ref)
```

Signal dimensionality from signal dict reference

If the signal has no axes, it is assumed to be an axis and expected to have dimension 1. Otherwise the expected dimension is equal to the number of signal axes.

```
store_channel_data_as_signal(computer_id, board_id, channel_id, generic_signal_id,
                              record_numbers, time_axis_id=None, time0=None, off-
                              set=None, coefficient=None, coefficient_V2unit=None,
                              delete=True, note=None, enable_deleted=False, re-
                              store_default=True)
```

Store data signal(s) from a particular DAQ channel as a given generic signal

#### Parameters

- **computer\_id** – numeric id or computer name
- **board\_id** –
- **channel\_id** –
- **generic\_signal\_id** –
- **record\_numbers** – list (or any iterable) of record numbers (or a single record number)
- **time\_axis\_id** – time axis (generic signal) id, latest revision will be used
- **time0** – time0, None keeps the last value
- **offset** – set offset (in raw levels), None keeps the last value
- **coefficient** – DAQ levels to Volts coefficient, None keeps the last value
- **coefficient\_V2unit** – DAQ Volts to units coefficient, None keeps the last value
- **delete** – delete the previous data signal
- **enable\_deleted** – do it even if there is no generic signal that is not deleted.
- **restore\_default** – if another data signal belongs to the g.s., it is restored to default g.s.

**store\_signal** (*generic\_signal\_id*, *\*\*kwargs*)

Insert new signal into SQL database. Does not create/write to the data file!

**Parameters** **generic\_signal\_id** – generic signal id

**kwargs** - fields from data signal reference:

**Parameters**

- **record\_number** – must be specified
- **data\_file\_id** –
- **data\_file\_key** –
- **time0** –
- **coefficient** –
- **offset** –
- **coefficient\_V2unit** –
- **variant** –
- **time\_axis\_id** –
- **time\_axis\_revision** –
- **axis1\_revision** –
- **axis2\_revision** –
- **axis3\_revision** –
- **axis4\_revision** –
- **axis5\_revision** –
- **axis6\_revision** –
- **time\_axis\_variant** –
- **axis1\_variant** –



- **axis2\_variant** –
- **axis3\_variant** –
- **axis4\_variant** –
- **axis5\_variant** –
- **axis6\_variant** –
- **note** – text note (new revision reason etc.)
- **computer\_id** – specify originating A/D
- **computer\_name** – alternative to computer\_id
- **board\_id** –
- **channel\_id** –
- **data\_quality** – enum('UNKNOWN', 'POOR', 'GOOD', 'VALIDATED')
- **deleted** –
- **daq\_parameters** – DAQ parameters (setup)
- **gs\_parameters** – generic signal parameters (setup)
- **ignore\_if\_exists** – Will not write if such signal exists (useful when multiple sources want to write it at the same time)

Attachment information (computer, board & channel) is automatically filled in if not specified (using *get\_attachment\_table* method).

**update\_signal** (*str\_id*='', *sig\_ref*=None, *\*\*kwargs*)

Create a new revision of a signal with some properties changed.

Use either *str\_id* or *sig\_ref* as signal identifier.

#### Parameters

- **str\_id** – signal string id
- **sig\_ref** – signal reference (dict)
- **kwargs** – new values of the fields

**class** pyCDB.client.CDBSignal (*data*=None, *name*='', *units*='', *axes*=None, *description*='', *ref*=None, *gs\_ref*=None, *file\_ref*=None, *daq\_attachment*=None)

CDB signal class, contains description and data

**get\_axes\_info** (*remove\_empty\_dim*=True)

Return list of dicts describing each viable axis

**get\_log\_level** ()

Get logging level

**info** ()

Information string about the signal

**log\_level**

Get logging level

**plot** (*fig*=None, *subplot*=111, *mplbackend*=None, *down\_sample*=1, *start\_time*=-1, *stop\_time*=-1, *y\_start*=-1, *y\_stop*=-1, *xmin*=None, *xmax*=None, *color\_set*=None, *plot\_kwargs*={}, *contour*=True, *colormap*=True, *\*\*kwargs*)

Plots the signal

**Parameters**

- **fig** – existing figure to add plot to
- **subplot** – subplot to plot to
- **mplbackend** – matplotlib backend to use (None - use default), use 'Agg' for no screen output
- **color\_set** – Color set of the figure, currently supported None (default) and logbook
- **plot\_kwargs** – dictionary with keyword arguments for the matplotlib plotting command
- **(only 2D) (colormap)** – show contour map
- **(only 2D)** – show color map

**set\_log\_level** (*value*)  
Set logging level

**url**  
WebCDB URL.

`pyCDB.client.apply_calibration` (*data*, *calibration*)  
Apply calibration to RAW data.

**Parameters**

- **data** – numpy array (or anything that can be multiplied and added)
- **calibration** – dictionary having keys [offset, coefficient] (signal\_ref can be used)

**Equation:**  $(data + offset) * coefficient$

`pyCDB.client.decode_generic_signal_strid` (*generic\_signal\_strid*)  
Decode generic signal string reference and return a dictionary

**Parameters** **generic\_signal\_strid** – is defined as 'alias or generic\_signal\_name' or 'generic\_signal\_name/data\_source\_name' or 'generic\_signal\_name/data\_source\_id'

**Return type** dictionary containing 'alias\_or\_name' or ('generic\_signal\_name' and ('data\_source\_id' or 'data\_source\_name'))

`pyCDB.client.decode_signal_strid` (*str\_id*)  
Decode signal string id

**Parameters** **str\_id** – string identifier, defined as id\_type:identification[units\_or\_slice][units\_or\_slice]

units\_or\_slice can be either units or slice specification. If two brackets are supplied, the second section in brackets has to differ from the first one. id\_type := CDB | FS | DAQ

**identification := (based on id\_type)**

- CDB -> generic\_signal\_strid:record\_number:variant:revision
- FS -> nodeuniqueid/hardwareuniqueid/parameteruniqueid:record\_number:variant:revision
- DAQ -> computer\_id/board\_id/channel\_id:record\_number:variant:revision \* if computer\_id is a string (contains non-numeric characters), it's understood as computer\_name \* alternatively, computer\_id and channel\_id can be prefixed by "board\_" and "channel\_" \* in DAQ and FS, you can use "." instead of "/"

**special units:**

- [RAW] or [raw] - get raw signal

- [DAV] - signal in Data Acquisition Volts (applied `get_coefficient_lev2V`)
- [] or [DEFAULT] or [default] - default units
- [unit\_name] or [unit\_system\_name] - conversion to unit stored in the database

certain parts can be omitted, the defaults are:

- `id_type = CDB`
- `revision = -1`
- `variant = ''`
- `record_number = -1`
- default units

i.e., `"ne/thomson:100" == "CDB:ne/thomson:100:-1"`

```
pyCDB.client.format_record_path(record_number, record_type)
```

Return data path for a given record number

**Return type** string with the subdirectory name

```
pyCDB.client.get_h5_dataset(file_path='', file_key='', file_format='HDF5', slice_=Ellipsis)
```

Get the HDF5 dataset.

```
pyCDB.client.get_h5_dataset_shape(file_path='', file_key='', file_format='HDF5')
```

Get the HDF5 dataset shape.

```
pyCDB.client.read_file_data(file_path='', file_key='', file_format='HDF5', slice_=Ellipsis)
```

Read data from a file

#### Parameters

- **file\_path** – file path of the file
- **file\_key** – data location in the file (groupA/groupB/.../dataset for HDF5)
- **file\_format** – data format of the file
- **file\_handle** – file handle if file is already open
- **slice** – slice specification

```
pyCDB.client.validate_variant(variant)
```

Check that the supplied variant follows the rules.

## pyCDB.DAQClient

@package DAQClient COMPASS Data Acquisition module

Contains class for DAQ management.

```
class pyCDB.DAQClient.CDBDAQClient(host=None, user=None, passwd=None, db=None,
                                     port=3306, log_level=None, data_root=None)
```

CDB data acquisition class - extends CDBClient

```
attach_channel(pc_id, board_id, channel_id, generic_signal_id, offset=None, coefficient_lev2V=None, coefficient_V2unit=None, uid=None, note='', parameters=None, timestamp=None, force_detach=False, **kwargs)
```

Attach a channel to a generic signal (specified by id)

Parameters that are not specified (default values None) will be inherited from previous channel attachment record.

parameters is a JSON text that describes channel parameters

**attached\_channel\_reference** (*pc\_id, board\_id, channel\_id*)

Checks whether a physical channel is attached

**change\_calibration** (*generic\_signal\_id=None, computer\_id=None, board\_id=None, channel\_id=None, \*\*kwargs*)

Shortcut method for changing the calibration of channel / generic signal)

Specify channel / generic signal by either choice (overspecification will be rejected): 1) *generic\_signal\_id*  
2) *computer\_id, board\_id, channel\_id*

Arguments that are specified, will be changed, other won't.

**create\_DAQ\_channels** (*nodeuniqueid, hardwareuniqueid, parameteruniqueid, data\_source\_id=None, note=', time\_axis\_id=None, axis1\_id=None, axis2\_id=None, axis3\_id=None, axis4\_id=None, axis5\_id=None, axis6\_id=None, \*\*kwargs*)

Creates a new DAQ channel.

It enables automatic creation of new channels and boards (if requested). New computers cannot be created automatically. Use *create\_computer* method in such case.

**Parameters data\_source\_id** – Not required only if you specify *generic\_signal\_id*

Keyword arguments:

#### Parameters

- **default\_generic\_signal\_id** – If set, it is used as default generic signal instead of creating a new one.
- **computer\_id** – Valid computer id that will be used if there is no channel for the computer yet.
- **board\_id** – Board id for a non-existent board (if found, stored *board\_id* is always used). If set to -1, the number will be evaluated automatically. If not set, exception is thrown. Please note that two *hardwareuniqueid*'s can be used as one board (although not recommended).
- **channel\_id** – If set, it will be used. Otherwise (or value of -1), a new one will be created.
- **coefficient\_lev2V** – *coefficient\_lev2V* for the default channel attachment.
- **offset** – offset for the default channel attachment.

**create\_computer** (*computer\_name, description=None, location=None*)

Create new computer in the database.

**detach\_channel** (*pc\_id, board\_id, channel\_id, timestamp=None*)

Detach channel in the channel setup

**find\_daq\_signals** (*record\_number=-1*)

Find all data signals that originate in any DAQ channel.

**find\_duplicate\_signals** (*\*record\_numbers*)

Find duplicate signals from DAQ channels in a set of records.

Example: *cdb.find\_duplicate\_signals(6228, 6229)*

**Parameters \*record\_numbers** – two or more record numbers

**get\_FS\_attachement** (*nodeuniqueid, hardwareuniqueid, parameteruniqueid*)  
Return the currently attached generic signal reference and the corresponding line from the channel\_attachments view for a FireSignal channel

**get\_board\_id** (*description='', description2=''*)  
Returns board ID by computer and board name

**get\_record\_number\_from\_FS\_event** (*event\_number, event\_id='0x0000'*)  
Return record number for a specific FireSignal event.

**Event\_number** Number of the FS event.

**Event\_id** string with the ID of the event. Default is "0x0000", shot event.

**is\_channel\_attached** (*pc\_id, board\_id, channel\_id*)  
Checks whether a physical channel is attached

**is\_generic\_signal\_attached** (*generic\_signal\_id*)  
Checks whether a generic signal is attached to a physical channel

**nodiff\_signal\_data** (*signal1, signal2*)  
Compare data of two signals. Returns True for identical signals.

#### Parameters

- **signal1** – first signal str\_id or reference
- **signal2** – second signal str\_id or reference

**signal\_setup** (*generic\_signal\_id, parameters='', timestamp=None, uid=None*)  
Write generic signal parameters (setup) into signal\_setup

#### Parameters

- **generic\_signal\_id** – numerical generic\_signal\_id
- **parameters** – signal parameters (typically JSON/YAML formatted text)
- **timestamp** – time, default is now
- **uid** – user id

**switch\_channels** (*pc\_id\_1, board\_id\_1, channel\_id\_1, pc\_id\_2, board\_id\_2, channel\_id\_2, uid=None, note='', timestamp=None*)  
Switch 2 channels attached generic signals

Keeps the offset and coefficient from the original channel, i.e., the offset and the coefficient are switched as well.

## pyCDB.pyCDBBase

Created on Jun 13, 2012

@author: jurban

**exception** pyCDB.pyCDBBase.CDBException  
A generic CDB exception class

**class** pyCDB.pyCDBBase.OrderedDict (\*args, \*\*kwargs)  
Customized ordered dictionary

**diff** (*other, mode='norm'*)  
Find differences to another OrderedDict, ignoring the keys order

**Parameters** **other** – OrderedDict object to compare to

**classmethod** **load\_h5** (*filename*)

Deserialize from an HDF5 file (created by save\_h5)

**Parameters** **filename** – input file name

**save\_h5** (*filename, mode='w'*)

Serialize to an HDF5 file

**Parameters**

- **filename** – output file name
- **mode** – file open mode, typically 'w' or 'a'

`pyCDB.pyCDBBase.cdb_base_check` (*func*)

Basic checking (database connection in particular) before calling CDB functions.

To be used as decorator for client methods. Raises an exception if problem occurs.

`pyCDB.pyCDBBase.isintlike` (*obj*)

Is object int-like?

In Python 2, this means int or long. In Python 3, this means just int.

`pyCDB.pyCDBBase.isstringlike` (*obj*)

Is object string-like?

In Python 2, this means string or unicode. In Python 3, this means just string.

**class** `pyCDB.pyCDBBase.pyCDBBase` (*host=None, user=None, passwd=None, db=None, port=3306, log\_level=None, data\_root=None*)

PyCDB base class for other pyCDB classes

**check\_connection** (*reconnect=True*)

Checks if connection is open, optionally (re)connects.

**close** ()

Close all connections (SQL)

**db**

Connection to the database.

- each thread has its local connection
- connection is automatically opened

**delete** (*table\_name, fields*)

Delete a row from the table.

**Parameters**

- **table\_name** – name of table to insert to
- **fields** – dictionary of fields identifying the row

Deletes at most one row: First, it checks for its existence in the table. If more than one row is found, raises exception. Returns True if something was deleted, False if not.

**error** (*exception*)

Log an error and throw the exception

**Parameters** **exception** – can be an exception or string with the message.

**file\_log\_level**

Get logging level

**classmethod filter\_sql\_str** (*sql\_str*, *extra\_chars*='\_', *quiet*=False)

Filter a string intended to be stored to SQL

**Parameters** *sql\_str* – input string

Keep only white listed characters: letters, digits, \_.

**classmethod format\_file\_name** (*collection\_name*, *revision*, *file\_format*='HDF5',  
*file\_extension*=None)

Returns generic file name

**classmethod get\_conf\_value** (*var\_name*, *section*='database', *required*=False)

Get configuration value for a given variable name

First looks for environment variables, then to ./CDBrc, then to ~/.CDBrc

**Parameters**

- **var\_name** – one of CDB\_HOST, CDB\_USER, CDB\_PASSWD, CDB\_DB
- **required** – if True, exception is raised telling the user that he/she has to set the option.

**Return type** variable value

**get\_file\_log\_level** ()

Get logging level

**get\_log\_level** ()

Get logging level

**get\_table\_struct** (*table\_name*)

Get a table structure (as a dict)

**insert** (*table\_name*, *fields*, *return\_inserted\_id*=False, *check\_structure*=True)

Insert a new row to a table.

**Parameters**

- **table\_name** – name of table to insert to
- **fields** – dictionary of ( field name, value ) to insert
- **return\_inserted\_id** – if true, returns the id of the inserted row
- **check\_structure** – if true, checks, whether fields agree with table layout,

Values are processed to be DB-friendly by `mysql_values` (see).

**classmethod is\_json** (*text*)

Checks whether text has a valid JSON syntax

**Parameters** *text* – text (string) to check

**log\_level**

Get logging level

**classmethod makedirs** (*path*[, *mode*=2047 ])

Super-mkdir; create a leaf directory and all intermediate ones. Works like mkdir, except that any intermediate path segment (not just the rightmost) will be created if it does not exist. This is recursive.

**classmethod mysql\_str** (*val*, *float\_format*='%.17g', *datetime\_format*='%Y-%m-%d %H:%M:%S',  
*quote*=True)

Convert Python value into MySQL string and quote if necessary

Use as the first character of MySQL builtin functions to avoid quoting

Lists and tuples are converted into parentheses-wrapped lists (using recursion) to work with IN operator.

**classmethod** `mysql_values` (*val\_list*, *float\_format*='%.17g', *datetime\_format*='%Y-%m-%d %H:%M:%S', *quote*=True)

Construct whole VALUES MySQL INSERT construct for a list of Python values

**query** (*sql*, *error\_if\_empty*=False, *warn\_if\_empty*=False, *error\_if\_multiple*=False)

Execute SQL query and returns a list of rows as dictionaries.

#### Parameters

- **sql** – valid SQL query
- **error\_if\_empty** – raises exception if the result set is empty
- **warn\_if\_empty** – issues warning if the result set is empty but successfully returns it
- **error\_if\_multiple** – the result has to include at maximum 1 row

**classmethod** `row_as_dict` (*row*, *description*)

Convert SQL row tuple to python dict

**set\_file\_log\_level** (*value*)

Set logging level

**set\_log\_level** (*value*)

Set logging level

**classmethod** `sql_to_python_type` (*sql\_type\_str*)

Transforms MySQL type definition to Python type and its length

**update** (*table\_name*, *id\_field*, *id\_value*, *fields*)

Update a row in a table.

#### Parameters

- **table\_name** – name of table to insert to
- **id\_field** – name of the unique key field
- **id\_value** – value of the unique key
- **fields** – dictionary of ( field name, value ) to insert

## pyCDB.logbook

Created on Jun 13, 2012

@author: jurban

**class** `pyCDB.logbook.logbook` (*host*=None, *user*=None, *passwd*=None, *db*=None, *port*=3306, *log\_level*=None)

logbook client

**campaign\_name** (*cid*)

Translates cid to the campaign name

**get\_shot\_comments** (*shot\_number*)

Get shot comments

**get\_shot\_info** (*shot\_number*)

Get shot info

**get\_shot\_linear\_profile\_value** (*shot\_number*, *prof\_name*)

Get a profile as an array



**get\_shot\_linear\_profiles** (*shot\_number*)  
Get all shot parameters (in a dictionary)

**get\_shot\_param\_value** (*shot\_number*, *param\_name*)  
Get the parameter's value

**get\_shot\_params** (*shot\_number*)  
Get all shot parameters (in a dictionary)

**get\_tags** (*shot\_number*)  
Get shot tags

**has\_tag** (*shot\_number*, *tag*)  
Check for a tag

**shot\_param\_name** (*param\_id*)  
Translate param\_id to param name

**uid\_name** (*uid*)  
Translates uid to (LDAP) name

## pyCDB.CodeGeneration

`pyCDB.CodeGeneration.generate_pxi` (*filename*='cyCDB\_api.pxi', *conn*=None)  
Generates .pxi file with C-api definitions

`pyCDB.CodeGeneration.get_table_struct` (*table\_name*, *conn*=None)  
Reads table structure and returns a dict containing row\_name : type



**class** `cdb_client`

CDB client class.

**get\_signal** (*str\_id*)

**Parameters** `str_id` – CDB string id (see *Signal id's*)

Get CDB signal by string id.

Returns a structure with the CDB signal, including the data and the description (references).

## Matlab errors

- `CDB:JyCDBError`
- `CDB:SignalError`
- `CDB:InputError`
- `CDB:StoreError`



---

## Installation of FireSignal with pyCDB

---

### Installation steps

#### INSTALLATION OF POSTGRES DB:

Exactly follow the instruction at IPFN IST web page: [http://metis.ipfn.ist.utl.pt/index.php?title=CODAC/FireSignal/Databases/PostgreSQL\\_%2b\\_FileSystem](http://metis.ipfn.ist.utl.pt/index.php?title=CODAC/FireSignal/Databases/PostgreSQL_%2b_FileSystem)

- Useful commands for work with POSTGRES:
- `sudo -u postgres psql`
- `\d`
- `\l`
- `\c genericdb`
- `SELECT * FROM events;`
- `\q`

#### HOW TO CREATE NEW CDB DATABASE?

```
mysql -u root -p
```

```
GRANT ALL PRIVILEGES ON . TO 'CDB'@'localhost' WITH GRANT OPTION; CREATE USER 'CDB'@'localhost' IDENTIFIED BY 'cmprSQLdata' ;
```

```
mysql -h localhost -u CDB -p < localhost.sql with password in ".CDBrc" file
```

**TO CHECK IT: connect CDB** show tables;

copy file .CDBrc to your home directory

```
mkdir /home/rrr/CDB_data
```

FINAL CHECK of SUCCESSFUL OPERATION OF pyCDB: `python pyCDB.py` → should some nice draw graph

## INSTALLATION OF FIRESIGNAL:

Source: `svn co svn://baco.ipfn.ist.utl.pt/scad/trunk`

`/trunk/java/dist$ sudo ./firesignal-1.1-linux-installer.bin`

IN INSTALLER: installed to: `/home/rrr/FS_PT` locate `javac` `/home/rrr/FS/jdk1.6.0_27/jre/` ..better to use sun's implementation 10.136.245.226 1050 PostgreSQL + filesystem `/home/rrr/CDB_myDATA` Database(Posgres) `genericdbadmin xxx genericdb localhost 5432`

## TO RUN CENTRAL SERVER:

`sudo ./fsignal start`

etc. (e.g `sudo ./fsignal_test_node start`)

## INSTALLATION OF FIRESIGNAL GUI CLIENT:

`sudo ./firesignaljws-1.1-linux-installer.bin`

`/usr/share/mini-httpd/html/XXX`

<http://localhost/XXX> – improve setting of http server 127.0.0.1 or localhost 1050 ok to the rest

TO RUN IT:

`/home/rrr/FS/jdk1.6.0_27/jre/bin/javaws FireSignalClient.jnlp`

`sudo ./fsignal_test_node start` Burn button → OK → you should see the data in the “Data” panel

## HDF5 plugin:

jep-2.4 compilovat s original java

**change scripts :** `/home/rrr/firesignal-1.1/dbcontroller/DBScript` `/home/rrr/firesignal-1.1/server/FSServerScript`

4x change IP 2x export `LD_PRELOAD=/usr/lib/libpython2.6.so.1.0` ( P by which was compiled jep, ldd..) 2x `FS_PT` path 1x `#../libs/libjhdf5.so` Xx parth `firesignal-1.1` ...

**copy libraries:** `sudo cp /home/rrr/firesignal-1.1/libs/jhdf5.jar /home/rrr/FS_PT/libs/` `sudo cp /home/rrr/firesignal-1.1/libs/libjhdf5.so /home/rrr/FS_PT/libs/` `sudo cp /home/rrr/firesignal-1.1/libs/hdf5test.jar /home/rrr/FS_PT/libs/` → `fsPqsqlHDF5.jar` in mail 19.10.2011 `sudo cp /home/rrr/firesignal-1.1/libs/jep.jar /home/rrr/FS_PT/libs/`

`rrr@rrr-laptop:~/FS_PT/init$ sudo ./fsignal start`

TO GET INFOS ABOUT START OF FS: either see the logs

or

`sudo gedit ../server/StartFSServer` & comment lines: `/home/rrr/FS_PT/server/FSServerScript \ #1>>/dev/null \ #2>>/dev/null`

`sudo gedit ../dbcontroller/StartDBController` &

## Nota Bene: How many instances of FSs are running?

If you experience some very strange/unexpected behaviour of FS maybe the error can be caused by that more than one instance of firesignal are running.

```
ps aux | grep java
```





---

### Compilation of FireSignal:

---

in NetBeans IDE -> new project from resources -> choose all java structure

You will need many libraries. Usually they are common and publicly accessible.

SDAS Core Libraries (SDAS.jar) and SDAS Client (SDASClient.jar) you can find here: <http://metis.ipfn.ist.utl.pt/CODAC/SDAS/Download>

### Example: LAST RECORD NUMBER IN GUI CLIENT

```
export LD_PRELOAD=/usr/lib/libpython2.6.so.1.0 ..needs (at least) pymysql export PYTHON-  
PATH=$PYTHONPATH:/home/rrr/workspace/pyCDB/src/
```

```
java -Djava.library.path=/home/rrr/FS_PT/libs/./usr/local/lib/./usr/lib -jar ./dist/FS.jar ..location of c-libs
```

NOTE: Done just in FS.jar X not as FireSignalClient.jpnl



### Javadoc

**cz.cas.ipp.compass.jycdb**

#### CDBClient

public class **CDBClient** extends *PythonAdapter*

Jython-based adapter of python CDBClient. Threading: Each thread should acquire an instance of this class using getInstance(). Multiple calls in one thread result in returning the same object. There is no need of closing the connections, it happens automatically in garbage collector (this is forced when we run out of connections). However, you should not pass reference to one instance among different threads (if you cannot assure the threads will not use it concurrently). The methods map to Client/DAQClient python methods as closely as possible. Several methods have overloaded versions: 1) a totally generic one that accepts ParameterList (a “dictionary” of values) 2,...) more specific versions with frequently used sets of parameters. See pyCDB documentation.

**Author** pipek

#### Methods

##### FS2CDB\_ref

public void **FS2CDB\_ref** (*String nodeUniqueId*, *String hardwareUniqueId*, *String parameterUniqueId*)

Return the CDB DAQ\_channel reference of a FireSignal channel

##### Parameters

- **nodeUniqueId** –
- **hardwareUniqueId** –
- **parameterUniqueId** –

### checkConnection

```
public boolean checkConnection ()
```

### createComputer

```
public Long createComputer (String computerName, String description, String location)
```

### createComputer

```
public Long createComputer (String computerName, String description)
```

### createComputer

```
public Long createComputer (String computerName)
```

### createDAQChannel

```
public void createDAQChannel (String nodeId, String hardwareId, String parameterId, long dataSourceId,  
                               Long[] axisIds)
```

### createDAQChannel

```
public void createDAQChannel (ParameterList parameters)
```

### createGenericSignal

```
public void createGenericSignal (String name, String alias, long dataSourceId, Long[] axisIds, String  
                                  units, String signalType)
```

### createGenericSignal

```
public void createGenericSignal (ParameterList parameters)
```

### createGenericSignalWhichReturnsLong

```
public Long createGenericSignalWhichReturnsLong (String name, String alias, long dataSourceId, Long[] axisIds, String units, String  
                                                    signalType)
```

### createGenericSignalWhichReturnsLong

```
public Long createGenericSignalWhichReturnsLong (ParameterList parameters)
```

### createRecord

public long **createRecord** (*Map*<*String*, *PyObject*> *parameters*)

### createRecord

public long **createRecord** (*String* *recordType*)

### createRecord

public long **createRecord** (long *fsEventNumber*, *String* *fsEventID*)

Creates new FireSignal record with forced record number equivalent to *fsEventNumber*.

#### Parameters

- **fsEventID** – event identification “0x0000” -
- **fsEventNumber** – firesignal shot number

### createRecord

public long **createRecord** (long *fsEventNumber*, *String* *fsEventID*, *String* *recordType*)

### deleteSignal

public void **deleteSignal** (*ParameterList* *parameters*)

### finalize

public void **finalize** ()

Even if finalize is not always called, try to close the connection.

### findGenericSignals

public *List*<*GenericSignal*> **findGenericSignals** (*String* *aliasOrName*, boolean *useRegexp*)

### getAttachment

public *ChannelAttachment* **getAttachment** (long *genericSignalId*)

### getAttachmentTable

public *List*<*ChannelAttachment*> **getAttachmentTable** (*ParameterList* *parameters*)

### getAttachmentTable

public *ChannelAttachment* **getAttachmentTable** (long *computerId*, long *boardId*, long *channelId*)

### getComputerId

public Long **getComputerId** (String *computerName*, String *description*)

### getComputerId

public Long **getComputerId** (String *computerName*)

### getDataFile

public *DataFile* **getDataFile** (*ParameterList* *parameters*)

### getDataFile

public *DataFile* **getDataFile** (long *dataFileId*)

### getDataSourceId

public Long **getDataSourceId** (String *dataSourceName*)

### getFSSignal

public *GenericSignal* **getFSSignal** (String *nodeId*, String *hardwareId*, String *parameterId*)

Get current generic signal attached to node, board & parameter.

**Returns** null if not found. Use Firesignal names.

### getGenericSignal

public *GenericSignal* **getGenericSignal** (long *genericSignalId*)

### getGenericSignal

public *GenericSignal* **getGenericSignal** (String *strid*)

### getGenericSignals

public List<*GenericSignal*> **getGenericSignals** (*ParameterList* *parameters*)

## getGenericSignals

```
public List<GenericSignal> getGenericSignals (String strid)
```

## getInstance

```
public static synchronized CDBClient getInstance ()
```

Get CDB client specific for current thread. It can serve a few hundred connections before problem arises. However, at that moment, there is a slight pause during which we try to close all unused connections (by hinting garbage connection) and obtain the connection for the second time. Only after this it eventually fails. Note: you should not pass CDBClient instances among threads or undefined behaviour results. But it is possible if you know what you are doing (e.g. in FS database controller).

**Returns** null if not connected.

## getInstance

```
public static synchronized CDBClient getInstance (ParameterList parameters)
```

Get CDB Client with specific construction parameters.

### Parameters

- **parameters** – (host, user, passwd, db, port, log\_level, data\_root). If any of the parameters is not specified, default value (from env. variables, etc.) is taken. Warning: This version of constructor is not recommended for general use. It doesn't use any clever method of preserving connections, memory, thread safety etc.

## getRecordNumberFromFSEvent

```
public long getRecordNumberFromFSEvent (long eventNumber, String eventId)
```

## getSignal

```
public Signal getSignal (String strid, String variant)
```

## getSignal

```
public Signal getSignal (String strid)
```

## getSignal

```
public Signal getSignal (long recordNumber, long genericSignalId)
```

## getSignal

```
public Signal getSignal (ParameterList parameters)
```

### getSignalCalibration

public *SignalCalibration* **getSignalCalibration** (*String strId*)

### getSignalCalibration

public *SignalCalibration* **getSignalCalibration** (*ParameterList* parameters)

### getSignalParameters

public *SignalParameters* **getSignalParameters** (*ParameterList* parameters)

### getSignalReference

public *SignalReference* **getSignalReference** (long *recordNumber*, long *genericSignalId*, long *revision*)

### getSignalReference

public *SignalReference* **getSignalReference** (long *recordNumber*, long *genericSignalId*)

### getSignalReference

public *SignalReference* **getSignalReference** (*String strId*)

### getSignalReferences

public *List*<*SignalReference*> **getSignalReferences** (*ParameterList* parameters)

### insert

public *Long* **insert** (*String tableName*, *ParameterList* fields, *Boolean* returnInsertedId, *Boolean* checkStructure)

### insert

public *Long* **insert** (*String tableName*, *ParameterList* fields, *Boolean* returnInsertedId)

### insert

public *Long* **insert** (*String tableName*, *ParameterList* fields)

### lastRecordNumber

public long **lastRecordNumber** (*String recordType*)



### lastShotNumber

public long **lastShotNumber** ()  
Last shot (or experimental record) number.

### newDataFile

public *DataFile* **newDataFile** (*ParameterList* parameters)

### newDataFile

public *DataFile* **newDataFile** (*String* collectionName, long recordNumber, long dataSourceId)

### recordExists

public boolean **recordExists** (long recordNumber)  
Check whether a record exists.

### setFileReady

public void **setFileReady** (long fileId)

### storeLinearSignal

public SignalReference **storeLinearSignal** (long genericSignalId, long recordNumber, double time0,  
double coefficient, double offset)  
Store signal of type LINEAR.

### storeSignal

public SignalReference **storeSignal** (*ParameterList* parameters)

### storeSignal

public SignalReference **storeSignal** (long genericSignalId, long recordNumber, *String* dataFileKey, long  
dataFileId, double time0, double coefficient, double offset, double co-  
efficient\_V2unit)  
Store signal of type FILE. Underlying Python method ensures that computer, board and channel ids are correctly set.

### storeSignal

public SignalReference **storeSignal** (long genericSignalId, long recordNumber, *String* dataFileKey, long  
dataFileId, double time0)  
Store signal of type FILE. Underlying Python method ensures that computer, board and channel ids are correctly set.

## updateSignal

public void **updateSignal** (*ParameterList* parameters)

Update signal.

### Parameters

- **parameters** – As few parameters as needed.

## CDBConnectionException

public class **CDBConnectionException** extends *CDBException*

The connection to database does not work. This means that the pyCDB was correctly loaded but there is a problem inside it or between it and MySQL.

**Author** pipek

## Constructors

### CDBConnectionException

public **CDBConnectionException** (*String* message)

### CDBException

public class **CDBException** extends *Exception*

Some problem with CDB. Common parent for more specific exception classes.

**Author** pipek

## Constructors

### CDBException

public **CDBException** (*String* message)

## ChannelAttachment

public class **ChannelAttachment** extends *DictionaryAdapter*

## Methods

### create

public static *ChannelAttachment* **create** (PyObject *object*)

### getAttachedGenericSignalId

public long **getAttachedGenericSignalId** ()

**getCoefficientLev2V**

public double **getCoefficientLev2V** ()

**getCoefficientV2Unit**

public double **getCoefficientV2Unit** ()

**getParameters**

public [Object](#) **getParameters** ()

**getParametersString**

public [String](#) **getParametersString** ()

**DataFile**

public class **DataFile** extends [DictionaryAdapter](#)  
Wrapper around rows of *data\_files* table.

**Author** honza

**Methods****create**

public static [DataFile](#) **create** (PyObject *object*)

**getCollectionName**

public [String](#) **getCollectionName** ()

**getFileName**

public [String](#) **getFileName** ()

**getFullPath**

public [String](#) **getFullPath** ()  
Full path to the file containing data.

**getId**

public long **getId** ()

## FSNodeWriter

public class **FSNodeWriter**

Methods for FireSignal nodes. It is mostly a facade over CDBClient methods to simplify development of nodes that communicate directly with CDB. How to write a signal? 1) For each record, create an instance of FSNodeWriter (you have to know generic signal id as well). 2) Create a file using createDataFile(). (see variants of this method) 3) Write data (as many times as you want). a) Using writeData(). Please, notice that you have to be aware of the data offset in HDF5 (or write them in one step). b) If you have an existing HDF5 file, use copyHdf5(). 4) Tell DB that the file is ready using setFileReady(). (if you changed it) 5) Write all axes using writeAxis(index). index=0 for time axis, index=1,2,3,... for other axes. (can be called multiple times for the same axis). 6) Write signal using writeSignal().

**Author** pipek

## Fields

### gson

Gson **gson**

## Constructors

### FSNodeWriter

public **FSNodeWriter** (long *genericSignalId*, long *recordNumber*, String *collectionName*, String *fileKey*)

#### Parameters

- **collectionName** – Name of the file without prefixes.
- **fileKey** – Name of the dataset in HDF5 file.

#### Throws

- **cz.cas.ipp.compass.jycdb.CDBException** –

### FSNodeWriter

public **FSNodeWriter** (long *computerId*, long *boardId*, long *channelId*, long *recordNumber*, String *collectionName*, String *fileKey*)

## Methods

### addDaqParameter

public void **addDaqParameter** (String *key*, Object *value*)

### addGenericSignalParameter

public void **addGenericSignalParameter** (String *key*, Object *value*)

### copyHdf5

public void **copyHdf5** (*String path*)

Copy an existing HDF5 file to the correct destination.

#### Parameters

- **path** – Local path of the HDF5 file.

#### Throws

- **IOException** – If the destination exist, copying is prevented and exception thrown.

### createDataFile

public void **createDataFile** ()

Create new data file (row in *data\_files*).

#### Throws

- **CDBException** – In this default version, an already existing file with requested properties is taken as error.

### createDataFile

public boolean **createDataFile** (boolean *okIfExists*, boolean *createIfExists*)

Create new data file (row in *data\_files*).

#### Parameters

- **okIfExists** – If false, an exception is thrown if a file with the same collection name, data source and record number exists
- **createIfExists** – If true and a file already exists, create a new one (otherwise the existing one is returned).

#### Throws

- **CDBException** –

**Returns** true if a file was created, false if nothing happens.

### getDataFile

public *DataFile* **getDataFile** ()

### getFileKey

public *String* **getFileKey** ()

### getGenericSignal

public *GenericSignal* **getGenericSignal** ()

### getRecordNumber

public long **getRecordNumber** ()

### setFileReady

public void **setFileReady** ()

Tell CDB that you have finished writing data. After this, you cannot write more data.

### setRequireChannelAttachment

public void **setRequireChannelAttachment** (boolean *require*)

Set whether writing should succeed (assuming value 1.0 for coefficients) when channel attachment not found. It has to be specified explicitly because it is quite probable that non-existence marks a bug.

### writeAxis

public void **writeAxis** (int *axis*, double *coefficient*, double *offset*)

Write one of the axes.

#### Parameters

- **axis** – Index of the axis (0=time, 1,2,3,...=other axes)

### writeData

public void **writeData** (*Data* *data*, long[] *dataOffset*)

Write data (with possible offset to append).

#### Parameters

- **data** – Data to be written.
- **dataOffset** – Zero-based index of the first data element. Size is calculated automatically. If data Offset is null, starts from beginning (fails with existing file and dataset).

#### Throws

- *cz.cas.ipp.compass.jycdb.util.Hdf5WriteException* –

### writeNotAttachedSignal

public void **writeNotAttachedSignal** (double *time0*, double *coefficient*)

Write signal info to CDB for signals that have no DAQ attachment. Usage is similar to writeSignal.

### writeSignal

public void **writeSignal** (double *time0*)

Write signal info to CDB. Call this after you have finished with writing data and axes. Can be used only for generic signal with valid DAQ attachment.

## FSWriter

public class **FSWriter** extends *PythonAdapter*

Java wrapper for quick storing of signal data to CDB in database controller. See: fs\_writer.py Warning: This class is meant to be used only in FireSignal database controller! Not elsewhere, even the nodes. Unexpected behaviour may result.

**Author** pipek

## Constructors

### FSWriter

public **FSWriter** (*CDBClient* client, *String* nodeId, *String* hardwareId, *String* parameterId, long recordNumber, double time0, double sampleLength)

## Methods

### getFileKey

public *String* **getFileKey** ()  
Get the name of the dataset to write to.

### getFilePath

public *String* **getFilePath** ()  
Get the full path where to write the file.

### readSignalInfo

public void **readSignalInfo** ()  
1st step - read all information we need for file storage.

### signalExists

public boolean **signalExists** ()  
Is the signal already stored in the database?

### storeSignalAndFile

public void **storeSignalAndFile** ()  
2nd step - write all that is to be written to database. This step should be undertaken only if signal does not exist.

## GenericSignal

public class **GenericSignal** extends *DictionaryAdapter*  
Wrapper around rows of *generic\_signals* table.

**Author** honza

## Fields

### FILE

public static final String **FILE**

### LINEAR

public static final String **LINEAR**

### MAX\_AXES

public static int **MAX\_AXES**

### VALID\_TYPES

public static final String[] **VALID\_TYPES**

## Methods

### create

public static *GenericSignal* **create** (PyObject *object*)

### getAlias

public String **getAlias** ()

### getAxis

public *GenericSignal* **getAxis** (int *index*)

### getAxisId

public long **getAxisId** (int *index*)  
Get id of one of the axes (time or spatial).

#### Parameters

- **index** – Number of the axis (starting with 1). If 0, time axis is returned.



**Throws**

- *CDBException* –

**Returns** 0 if there is no axis.

**getDataSourceId**

public long **getSourceId** ()

**getDescription**

String **getDescription** ()

**getId**

public long **getId** ()

**getLastRecordNumber**

public long **getLastRecordNumber** ()

**getName**

public String **getName** ()

**getSignalType**

public String **getSignalType** ()

**getTimeAxis**

public *GenericSignal* **getTimeAxis** ()

**getTimeAxisId**

public long **getTimeAxisId** ()

**getUnits**

public String **getUnits** ()

**isFile**

public boolean **isFile** ()

## isLinear

public boolean **isLinear** ()

## isValidType

public static boolean **isValidType** (*String signalType*)  
Check if a signal type is valid. Only “LINEAR” and “FILE” are accepted now.

## Signal

public class **Signal** extends *PythonAdapter*  
Class mimicking the CDBSignal Python class. In a few aspects, it is more object-oriented and lazy-evaluated to simplify manipulation in MATLAB/IDL.  
**Author** pipek

## Constructors

### Signal

public **Signal** (SignalReference *signalRef*, *String units*)

### Signal

protected **Signal** (*Signal parent*, int *nthChild*, *PythonAdapter signalTree*, *DictionaryAdapter unitFactorTree*)  
Constructor for dependent objects from signal tree with unit factor.

#### Parameters

- **parent** – Signal that uses this as an axis.
- **nthChild** – The order of this signal among axes in the parent signal.
- **signalTree** – Signal subtree.
- **unitFactorTree** – Unit factor subtree.

### Signal

protected **Signal** (*Signal parent*, int *nthChild*, *PythonAdapter signalTree*, *String units*)  
Constructor for dependent objects from signal tree without unit factor.

#### Parameters

- **parent** – Signal that uses this as an axis.
- **nthChild** – The order of this signal among axes in the parent signal.
- **signalTree** – Signal subtree.
- **units** – Most times, this will be “default”, but any other value can used.

## Methods

### getAxes

public [SortedMap](#)<[String](#), [Signal](#)> **getAxes** ()

A map of all existing axes.

**Throws**

- [CDBException](#) – Non-existent axes are omitted.

**Returns** A dictionary of [ axis\_name, Signal object ] values.

### getAxis

public [Signal](#) **getAxis** (int *i*)

An axis as Signal instance.

**Parameters**

- *i* – 0-time axis, N-axisN

**Throws**

- [CDBException](#) –

**Returns** null if not found

### getData

public synchronized [Data](#) **getData** ()

Numerical data of the signal. Note: this method is not used in bindings using JyCDB (IDL/MATLAB) due to efficiency. Try using native implementations of file reading.

### getDataFile

public [DataFile](#) **getDataFile** ()

### getDataFileId

public long **getDataFileId** ()

### getDescription

public [String](#) **getDescription** ()

Generic signal description.

### getGenericSignal

public [GenericSignal](#) **getGenericSignal** ()

### getGenericSignalId

public long **getGenericSignalId** ()

### getGenericSignalReference

public *GenericSignal* **getGenericSignalReference** ()

### getName

public *String* **getName** ()  
Generic signal name.

### getParameters

public *SignalParameters* **getParameters** ()

### getPythonObject

public PyObject **getPythonObject** ()  
Return the Python object of CDBSignal class.

**Returns** null if CDBSignal not found. Note: This method catches exceptions thus hiding problems.

### getRecordNumber

public long **getRecordNumber** ()

### getSignalCalibration

public *SignalCalibration* **getSignalCalibration** (boolean *includeUnitFactor*)  
Signal calibration for default/DAV/RAW.

**Throws**

- *CDBException* –

### getSignalReference

public *SignalReference* **getSignalReference** ()

### getTimeAxis

public *Signal* **getTimeAxis** ()  
Time axis as *Signal* instance. Just shortcut to *getAxis* method.

### getUnit

public **String** **getUnit** ()  
Unit name.

#### Throws

- *CDBException* – RAW/DAQ/default GS units or units from factor tree (if requested).

### getUnitFactor

public double **getUnitFactor** ()  
Unit factor, by which the default signal has to be multiplied.

#### Throws

- *CDBException* – If DAV/RAW/default are selected, this factor is 1.0. If another unit/unit system is selected, this is a conversion from the default unit to the selected one. Internally, this reflects `unit_factor_tree` returned from Python.

### hasUnitFactor

public boolean **hasUnitFactor** ()  
Whether the signal has associated unit factor tree. This is true only if unit or unit system were requested.

### isDAV

public boolean **isDAV** ()

### isFile

public boolean **isFile** ()

### isLinear

public boolean **isLinear** ()

### isRaw

public boolean **isRaw** ()

### SignalCalibration

public class **SignalCalibration** extends *DictionaryAdapter*  
Signal calibration. Meaning:  $\text{value} = (\text{data} + \text{offset}) * \text{coefficient}$

**Author** pipek

## Fields

### DAQ\_VOLTS

public static final *String* **DAQ\_VOLTS**

### RAW

public static final *String* **RAW**

## Methods

### apply

public double[] **apply** (double[] *array*)  
If you have an array of data, apply calibration to it. Returns new array.

### create

public static *SignalCalibration* **create** (PyObject *object*)

### getCoefficient

public double **getCoefficient** ()

### getOffset

public double **getOffset** ()

### getUnits

public *String* **getUnits** ()

## SignalParameters

public class **SignalParameters** extends *DictionaryAdapter*  
Signal parameters (both GS & DAQ-based). These are strings interpreted as JSON. Strings are accessible using `getDAQParametersString()` & `getGenericSignalParametersString()`, Parsed JSON objects (as trees) are accessible using `getDAQParameters()` & `getGenericSignalParameters()`.

**Author** pipek

## Methods

### create

public static *SignalParameters* **create** (PyObject *object*)

### getDAQParameters

public Object **getDAQParameters** ()

### getDAQParametersString

public String **getDAQParametersString** ()

### getGenericSignalParameters

public Object **getGenericSignalParameters** ()

### getGenericSignalParametersString

public String **getGenericSignalParametersString** ()

## WrongSignalTypeException

public class **WrongSignalTypeException** extends *CDBException*

Signal type is not valid for its intended use. It can be either LINEAR or FILE with different operations available.

**Author** pipek

## Constructors

### WrongSignalTypeException

public **WrongSignalTypeException** (String *message*)

## cz.cas.ipp.compass.jycdb.plotting

### DataPlotter

public class **DataPlotter** extends JFrame

## Constructors

### DataPlotter

public **DataPlotter** ()

## Methods

### appendChart

public static JFreeChart **appendChart** (JFreeChart *oldChart*, *Signal* *signal*)

### createChart

public static JFreeChart **createChart** (*Signal* *signal*)

### main

public static void **main** (String[] *args*)

## cz.cas.ipp.compass.jycdb.test

### JyCDBTest

public class **JyCDBTest**  
    **Author** pipek

## Methods

### fsWriterTest

public static void **fsWriterTest** ()

### instanceTest

public static void **instanceTest** ()

### main

public static void **main** (String[] *args*)

## cz.cas.ipp.compass.jycdb.util

### ArrayUtils

public class **ArrayUtils**  
    Utils for quick manipulation with (numeric) arrays. Warning: most methods work only on rectangular arrays.  
    **Author** pipek



## Methods

### add

public static double[] **add** (double[] *array*, double *offset*)

Add a constant to all elements of an array.

**Returns** a new array

### asDoubleArray

public static double[] **asDoubleArray** (Object *array*)

Cast all values of an array to double (runtime version).

#### Parameters

- **array** – Array of allowed type.

#### Throws

- *UnknownDataTypeException* –

### asDoubleArray

public static double[] **asDoubleArray** (long[] *array*)

Cast all values of an array to double.

### asDoubleArray

public static double[] **asDoubleArray** (int[] *array*)

Cast all values of an array to double.

### asDoubleArray

public static double[] **asDoubleArray** (short[] *array*)

Cast all values of an array to double.

### asDoubleArray

public static double[] **asDoubleArray** (float[] *array*)

Cast all values of an array to double.

### asIntArray

public static int[] **asIntArray** (long[] *array*)

Cast all values of an array to int.

## asLongArray

public static long[] **asLongArray** (int[] *array*)  
Cast all values of an array to long.

## dimensions

public static long[] **dimensions** (Object *array*)  
Get the dimensions along all axes of an array. Works only on rectangular arrays.

## get

public static Object **get** (Object *array*, int[] *index*)  
Get an element from a multidimensional array. Works on all arrays.

## linearTransform

public static double[] **linearTransform** (double[] *array*, double *multiplyBy*, double *add*)  
Apply a linear transformation  $y = ax + b$  on an array.

### Parameters

- **multiplyBy** – Multiplication constant
- **add** – Addition constant

**Returns** a new array

## linearTransform

public static Object **linearTransform** (Object *array*, double *multiplyBy*, double *add*)  
A general linear transformation of multidimensional array.

### Parameters

- **array** –
  - N-dimensional rectangular array of double/long/int/short/float
- **multiplyBy** –
  - multiplicative factor
- **add** –
  - additive factor

### Returns

- N-dimensional array of the same shape as input

## max

public static double **max** (double[] *array*)  
Maximum value found in the array.

## min

public static double **min** (double[] *array*)  
Minimum value found in the array.

## multiply

public static double[] **multiply** (double[] *array*, double *coefficient*)  
Multiply all elements of an array by a constant.

**Returns** a new array

## product

public static long **product** (int[] *array*)

## product

public static long **product** (long[] *array*)

## range

public static int[] **range** (int *xmin*, int *xmax*)

## rank

public static int **rank** (Object *array*)  
Get the number of dimensions of an array. Works only on rectangular arrays.

## reshape

public static Object **reshape** (Object *array*, int[] *newDimensions*)  
Reshape array. Create a new array with the same total size but with different dimensions. Last index is moving fastest, while the elements are copied one after another. Works only on rectangular arrays.

### Parameters

- **array** –
- **newDimensions** –

**Returns** new array

## set

public static void **set** (Object *array*, int[] *index*, Object *value*)  
Set an element in a multidimensional array. Works on all arrays.

## size

public static long **size** (*Object array*)

Total number of items in array. Works only on rectangular arrays.

## Data

public class **Data**

Wrapper around data. It enables us to live with just one copy of most methods and pass data around independently of its dimension and numerical type. It can be constructed using: 1) one of the 6x3 constructors with explicit array types 2) a general constructor using Object and dimensionality data specified by you.

**Author** pipek

## Fields

### DOUBLE

public static final int **DOUBLE**

### FLOAT

public static final int **FLOAT**

### INT16

public static final int **INT16**

### INT32

public static final int **INT32**

### INT64

public static final int **INT64**

### INT8

public static final int **INT8**

### MAX\_RANK

public static final int **MAX\_RANK**

## Constructors

### Data

public **Data** (byte[] *data*)

### Data

public **Data** (byte[][] *data*)

### Data

public **Data** (byte[][][] *data*)

### Data

public **Data** (short[] *data*)

### Data

public **Data** (short[][] *data*)

### Data

public **Data** (short[][][] *data*)

### Data

public **Data** (int[] *data*)

### Data

public **Data** (int[][] *data*)

### Data

public **Data** (int[][][] *data*)

### Data

public **Data** (long[] *data*)

### Data

```
public Data (long[][] data)
```

### Data

```
public Data (long[][][] data)
```

### Data

```
public Data (double[] data)
```

### Data

```
public Data (double[][] data)
```

### Data

```
public Data (double[][][] data)
```

### Data

```
public Data (float[] data)
```

### Data

```
public Data (float[][] data)
```

### Data

```
public Data (float[][][] data)
```

### Data

```
public Data (Object data, int type, long[] dimensions)
```

Generic constructor.

#### Parameters

- **data** – The correct array object (or null)
- **type** – Type of data in terms of this class constants.
- **dimensions** – Length along all dimensions. Gets copied. This is useful if we obtain data from external source, we know its properties but we don't want to cast them unnecessarily (like read from HDF5). Data needn't be specified (in such case a new array is created.)

## Methods

### asDoubleArray1D

public double[] **asDoubleArray1D** ()

Get a 1-D double array representation of data. Works for 1-D data. For doubles, it simply returns, for others, it converts them using `ArrayUtils.asDoubleArray()` (see).

#### Throws

- *WrongDimensionsException* – if data are not 1-D.
- *UnknownDataTypeException* – if conversion to double is not supported by underlying procedure.

### createRawObject

public static *Object* **createRawObject** (int *type*, long[] *dimensions*)

Create a multidimensional array object of a selected type and dimensions.

### flatten

public boolean **flatten** ()

Remove dimensions that have length 1. Preserves 1D arrays. Makes changes only if there is a trivial dimension.

**Returns** true if there was a change, false otherwise. Motivation: Our HDF5 file sometimes have Nx1 arrays instead of N arrays.

### from1DArray

public static *Data* **from1DArray** (*Object* *array*, int *type*, long[] *dimensions*)

Reshapes 1-D array to a multidimensional array of correct dimensions. Motivation: HDF5 library returns only 1-D arrays. In C order.

### getDimensions

public long[] **getDimensions** ()

Get length of the data along all dimensions.

### getRank

public int **getRank** ()

Get the number of dimensions. i.e. `double[] => 1`, `short[][] => 2` etc.

### getRawData

public *Object* **getRawData** ()

Get the original array used for the initialization of this object. No type information, you have to cast it yourself.

## getType

public int **getType** ()

Get the type of stored array (in terms of constants defined in this class). Our data types differ from HDF5 so that this class is not dependent on HDF5 library. Hdf5Utils class contains conversion routines.

## linearTransform

public *Data* **linearTransform** (double *mul*, double *add*)

## DebugUtils

public class **DebugUtils**

Utilities helping with debugging.

**Author** pipek

## Methods

### logWithDuration

public static void **logWithDuration** (*String message*)

Print a message alongside with the number of nanoseconds since last call of this method.

#### Parameters

- **message** – Message to log

### logWithDuration

public static void **logWithDuration** (*String message*, *Logger logger*)

Print a message alongside with the number of nanoseconds since last call of this method.

#### Parameters

- **message** – Message to log
- **logger** – Logger to use

## DictionaryAdapter

public class **DictionaryAdapter** extends *PythonAdapter*

Object that is represented by a python dictionary (or OrderedDict) and usually obtained as a row from database. Each of the classes should implement static method create(PyObject) which returns null if underlying Python object is None. There should be no public constructor in child classes.

## Constructors

### DictionaryAdapter

public **DictionaryAdapter** (PyObject *object*)



## Methods

### asMap

public [SortedMap](#) **asMap** ()  
Map python dictionary to java map. Motivation: Binding for MATLAB.

### get

public [PyObject](#) **get** ([String](#) key)  
Get dictionary value indexed by key as python object.

### getDictKeys

public [String](#)[] **getDictKeys** ()  
Get list of dictionary keys. Motivation: Binding for IDL.  
**Returns** Array of keys in the order they are stored in by Python.

### getDictValue

public [Object](#) **getDictValue** ([String](#) key)  
Get dictionary value indexed by key as “Java native” object. Motivation: Binding for IDL.  
**See also:** [PythonUtils.asNative](#) ([org.python.core.PyObject](#))

### getDictValueAsDouble

public double **getDictValueAsDouble** ([String](#) key)  
Get dictionary value indexed by key as double. Motivation: Binding for IDL.

### getDictValueAsLong

public long **getDictValueAsLong** ([String](#) key)  
Get dictionary value indexed by key as long. Motivation: Binding for IDL.

### getDictValueAsString

public [String](#) **getDictValueAsString** ([String](#) key)  
Get dictionary value indexed by key as string. Motivation: Binding for IDL.

## Hdf5ReadException

public class **Hdf5ReadException** extends [CDBException](#)  
**Author** pipek

## Constructors

### Hdf5ReadException

public **Hdf5ReadException** (*String message*)

### Hdf5Utils

public class **Hdf5Utils**

Utilities for reading/writing HDF5 files. Methods readData, writeData and writeData

**Author** pipek

## Methods

### readFileData

public static *Data* **readFileData** (*String filePath*, *String fileKey*)

Read data from file. Up to 3-D arrays are supported. (Higher rank would require reimplementations of Data class)

#### Parameters

- **filePath** – Path to the file.
- **fileKey** – Name of the dataset in the file.

### writeData

public static void **writeData** (*String filePath*, *String fileKey*, *Data data*)

Write (from beginning) new HDF5 dataset.

#### Parameters

- **filePath** – Path to the file (may or may not exist).
- **fileKey** – Name of the dataset (error if it exists).
- **data** – Wrapped data (see class Data).

### writeData

public static void **writeData** (*String filePath*, *String fileKey*, *Data data*, long[] *offset*)

Write (or append) data to an existing HDF5 dataset.

#### Parameters

- **filePath** – Path to the file.
- **fileKey** – Name of the (existing) dataset.
- **data** – Wrapped data (see class Data).
- **offset** – Offset (in dimensions) from which start with writing data. If offset == null, new data set is created. You have to be sure that your offset is correct. The dataset is enlarged if possible but you can overwrite existing data if you are not careful.

## writeData

```
public static void writeData (String filePath, String fileKey, int dataType, long[] dimensions, Object rawData, long[] offset)
```

Write data to a HDF5 dataset (generic version). This generic version is available as a public method, still it is more convenient (and practical) to use `writeData(String, String, Data)` or `writeData`.

### Parameters

- **filePath** – Path to the file (may or may not exist).
- **fileKey** – Name of the dataset.
- **dataType** – HDF5 data type.
- **dimensions** – The dimensions of the data array.
- **rawData** – The data array of any format.
- **offset** – Offset of data (null => create new dataset)

## Hdf5WriteException

```
public class Hdf5WriteException extends CDBException
```

An error when writing to HDF5.

**Author** pipek

## Constructors

### Hdf5WriteException

```
public Hdf5WriteException (String message)
```

## JsonUtils

```
public class JsonUtils
```

Utility to read and write data from/to JSON format.

**Author** pipek

## Methods

### asNative

```
public static Object asNative (JsonElement element)
```

Transform JSON element into native object. Works recursively: - double, boolean, string - native - arrays => Object[] - objects => TreeMap

### parseNative

```
public static Object parseNative (String source)
```

Take JSON string and turn it into native objects.

## ParameterList

public class **ParameterList** extends [HashMap<String, PyObject>](#)

List of parameters for python methods that accepts (via put method) some of the basic types in addition to default PyObject. Various put methods just simplify adding of native objects.

## Methods

### put

public [String](#) **put** ([String](#) key, [String](#) value)

### put

public long **put** ([String](#) key, long value)

### put

public int **put** ([String](#) key, int value)

### put

public [Date](#) **put** ([String](#) key, [java.sql.Date](#) value)

### put

public boolean **put** ([String](#) key, boolean value)

### put

public double **put** ([String](#) key, double value)

### put

public [Timestamp](#) **put** ([String](#) key, [Timestamp](#) value)

### put

public short **put** ([String](#) key, short value)

### put

public [ParameterList](#) **put** ([String](#) key, [ParameterList](#) value)

## PythonAdapter

public class **PythonAdapter**

Adapter of a PyObject that offers some shortcut methods for invoking etc. Motivation: The Jython API is not as elegant as would be desirable. However, with proper understanding of PyObject, this class would not be necessary.

**Author** pipek

## Fields

### PyObject

protected PyObject **pyObject**

## Constructors

### PythonAdapter

public **PythonAdapter** (PyObject *object*)

## Methods

### getAttribute

public PyObject **getAttribute** (String *name*)

Get the attribute of PyObject.

#### Parameters

- **name** – Name of the attribute.

**Returns** The attribute as PyObject or null (if attribute not present).

### getPythonObject

public PyObject **getPythonObject** ()

Get the unwrapped PyObject.

### hasAttribute

public boolean **hasAttribute** (String *name*)

### invoke

public PyObject **invoke** (String *method*, Map<String, PyObject> *parameters*)

Call a method on PyObject with parameters.

## invoke

public PyObject **invoke** (*String method*)  
Call a method on PyObject without parameters.

### Parameters

- **method** – Name of the method (callable attribute) to call on the object.

**Returns** Result of the call as PyObject

## PythonUtils

public class **PythonUtils**  
Utilities for easier manipulation with Python through Jython.

**Author** pipek

## Methods

### asNative

public static Object **asNative** (PyObject *object*)  
Interpret the python object as a native Java one.

**Returns** Object of the correct type. It has to be cast to be useful. Numbers converts to numbers (long or double). Dicts converts to TreeMap. Tuples and lists converts to Object[]. Dates converts to Date. NoneType converts to null. In case it does not know a type, it returns its string representation along with ! and type name. (e.g. "!representation of my weird type{ Weirdtype}")

### getInterpreter

public static synchronized PythonInterpreter **getInterpreter** ()  
Get a single copy of python interpreter. More of them should not be needed.

## invoke

public static PyObject **invoke** (PyObject *object*, *String method*, Map<*String*, PyObject> *parameters*)  
Invoke a method on a Python object with named parameters. It calls *invoke* method present in Jython, it only allows easier manipulation with parameters (enabling this to be called with ParameterList).

## newObject

public static PyObject **newObject** (*String className*, Map<*String*, PyObject> *parameters*)  
Create a new instance of a named python class.

## UnknownDataTypeException

public class **UnknownDataTypeException** extends *CDBException*

**Author** pipek

## Constructors

### UnknownDataTypeException

public **UnknownDataTypeException** (*String message*)

### WrongDimensionsException

public class **WrongDimensionsException** extends *CDBException*

**Author** pipek

## Constructors

### WrongDimensionsException

public **WrongDimensionsException** (*String message*)





# CHAPTER 11

---

## Indices and tables

---

- `genindex`
- `modindex`
- `search`



### p

- `pyCDB.client`, [21](#)
- `pyCDB.CodeGeneration`, [37](#)
- `pyCDB.DAQClient`, [31](#)
- `pyCDB.logbook`, [36](#)
- `pyCDB.pyCDBBase`, [33](#)



## Symbols

`__main__`  
module, 41

## A

`add(double[], double)` (Java method), 69  
`add_unit()` (pyCDB.client.CDBClient method), 21  
`addDaqParameter(String, Object)` (Java method), 56  
`addGenericSignalParameter(String, Object)` (Java method), 56  
`appendChart(JFreeChart, Signal)` (Java method), 68  
`apply(double[])` (Java method), 66  
`apply_calibration()` (in module pyCDB.client), 30  
`apply_unit_factor_tree()` (pyCDB.client.CDBClient method), 21  
`ArrayUtils` (Java class), 68  
`asDoubleArray(float[])` (Java method), 69  
`asDoubleArray(int[])` (Java method), 69  
`asDoubleArray(long[])` (Java method), 69  
`asDoubleArray(Object)` (Java method), 69  
`asDoubleArray(short[])` (Java method), 69  
`asDoubleArray1D()` (Java method), 75  
`asIntArray(long[])` (Java method), 69  
`asLongArray(int[])` (Java method), 70  
`asMap()` (Java method), 77  
`asNative(JsonElement)` (Java method), 79  
`asNative(PyObject)` (Java method), 82  
`attach_channel()` (pyCDB.DAQClient.CDBDAQClient method), 31  
`attached_channel_reference()` (pyCDB.DAQClient.CDBDAQClient method), 32

## C

`campaign_name()` (pyCDB.logbook.logbook method), 36  
`CDB2FS_ref()` (pyCDB.client.CDBClient method), 21  
`cdb_base_check()` (in module pyCDB.pyCDBBase), 34  
`cdb_client` (built-in class), 39  
`CDBClient` (class in pyCDB.client), 21

`CDBClient` (Java class), 47  
`CDBConnectionException` (Java class), 54  
`CDBConnectionException(String)` (Java constructor), 54  
`CDBDAQClient` (class in pyCDB.DAQClient), 31  
`CDBException`, 33  
`CDBException` (Java class), 54  
`CDBException(String)` (Java constructor), 54  
`CDBSignal` (class in pyCDB.client), 29  
`change_calibration()` (pyCDB.DAQClient.CDBDAQClient method), 32  
`ChannelAttachment` (Java class), 54  
`check_connection()` (pyCDB.pyCDBBase.pyCDBBase method), 34  
`checkConnection()` (Java method), 48  
`close()` (pyCDB.pyCDBBase.pyCDBBase method), 34  
`copyHdf5(String)` (Java method), 57  
`create(PyObject)` (Java method), 54, 55, 60, 66, 67  
`create_computer()` (pyCDB.DAQClient.CDBDAQClient method), 32  
`create_DAQ_channels()` (pyCDB.DAQClient.CDBDAQClient method), 32  
`create_generic_signal()` (pyCDB.client.CDBClient method), 21  
`create_record()` (pyCDB.client.CDBClient method), 22  
`createChart(Signal)` (Java method), 68  
`createComputer(String)` (Java method), 48  
`createComputer(String, String)` (Java method), 48  
`createComputer(String, String, String)` (Java method), 48  
`createDAQChannel(ParameterList)` (Java method), 48  
`createDAQChannel(String, String, String, long, Long[])` (Java method), 48  
`createDataFile()` (Java method), 57  
`createDataFile(boolean, boolean)` (Java method), 57  
`createGenericSignal(ParameterList)` (Java method), 48  
`createGenericSignal(String, String, long, Long[], String, String)` (Java method), 48  
`createGenericSignalWhichReturnsLong(ParameterList)` (Java method), 48

`createGenericSignalWhichReturnsLong(String, String, long, Long[], String, String)` (Java method), 48  
`createRawObject(int, long[])` (Java method), 75  
`createRecord(long, String)` (Java method), 49  
`createRecord(long, String, String)` (Java method), 49  
`createRecord(Map)` (Java method), 49  
`createRecord(String)` (Java method), 49  
`cz.cas.ipp.compass.jycdb` (package), 47  
`cz.cas.ipp.compass.jycdb.plotting` (package), 67  
`cz.cas.ipp.compass.jycdb.test` (package), 68  
`cz.cas.ipp.compass.jycdb.util` (package), 68

## D

`DAQ_VOLTS` (Java field), 66  
`Data` (Java class), 72  
`Data(byte[])` (Java constructor), 73  
`Data(byte[][])` (Java constructor), 73  
`Data(byte[][][])` (Java constructor), 73  
`Data(double[])` (Java constructor), 74  
`Data(double[][])` (Java constructor), 74  
`Data(double[][][])` (Java constructor), 74  
`Data(float[])` (Java constructor), 74  
`Data(float[][])` (Java constructor), 74  
`Data(float[][][])` (Java constructor), 74  
`Data(int[])` (Java constructor), 73  
`Data(int[][])` (Java constructor), 73  
`Data(int[][][])` (Java constructor), 73  
`Data(long[])` (Java constructor), 73  
`Data(long[][])` (Java constructor), 74  
`Data(long[][][])` (Java constructor), 74  
`Data(Object, int, long[])` (Java constructor), 74  
`Data(short[])` (Java constructor), 73  
`Data(short[][])` (Java constructor), 73  
`Data(short[][][])` (Java constructor), 73  
`DataFile` (Java class), 55  
`DataPlotter` (Java class), 67  
`DataPlotter()` (Java constructor), 67  
`db` (pyCDB.pyCDBBase.pyCDBBase attribute), 34  
`DebugUtils` (Java class), 76  
`decode_generic_signal_strid()` (in module pyCDB.client), 30  
`decode_signal_strid()` (in module pyCDB.client), 30  
`delete()` (pyCDB.pyCDBBase.pyCDBBase method), 34  
`delete_signal()` (pyCDB.client.CDBClient method), 22  
`deleteSignal(ParameterList)` (Java method), 49  
`detach_channel()` (pyCDB.DAQClient.CDBDAQClient method), 32  
`DictionaryAdapter` (Java class), 76  
`DictionaryAdapter(PyObject)` (Java constructor), 76  
`diff()` (pyCDB.pyCDBBase.OrderedDict method), 33  
`dimensions(Object)` (Java method), 70  
`DOUBLE` (Java field), 72

## E

`error()` (pyCDB.pyCDBBase.pyCDBBase method), 34  
execution  
    context, 41

## F

`FILE` (Java field), 60  
`file_log_level` (pyCDB.pyCDBBase.pyCDBBase attribute), 34  
`filter_sql_str()` (pyCDB.pyCDBBase.pyCDBBase class method), 34  
`finalize()` (Java method), 49  
`find_daq_signals()` (pyCDB.DAQClient.CDBDAQClient method), 32  
`find_duplicate_signals()` (pyCDB.DAQClient.CDBDAQClient method), 32  
`find_generic_signals()` (pyCDB.client.CDBClient method), 22  
`findGenericSignals(String, boolean)` (Java method), 49  
`flatten()` (Java method), 75  
`FLOAT` (Java field), 72  
`format_file_name()` (pyCDB.pyCDBBase.pyCDBBase class method), 35  
`format_record_path()` (in module pyCDB.client), 31  
`from1DArray(Object, int, long[])` (Java method), 75  
`FS2CDB_ref()` (pyCDB.client.CDBClient method), 21  
`FS2CDB_ref(String, String, String)` (Java method), 47  
`FSNodeWriter` (Java class), 56  
`FSNodeWriter(long, long, long, long, String, String)` (Java constructor), 56  
`FSNodeWriter(long, long, String, String)` (Java constructor), 56  
`FSWriter` (Java class), 59  
`FSWriter(CDBClient, String, String, String, long, double, double)` (Java constructor), 59  
`fsWriterTest()` (Java method), 68

## G

`generate_pxi()` (in module pyCDB.CodeGeneration), 37  
`GenericSignal` (Java class), 60  
`get(Object, int[])` (Java method), 70  
`get(String)` (Java method), 77  
`get_attachment_table()` (pyCDB.client.CDBClient method), 22  
`get_axes_info()` (pyCDB.client.CDBSignal method), 29  
`get_board_id()` (pyCDB.DAQClient.CDBDAQClient method), 33  
`get_channel_references()` (pyCDB.client.CDBClient method), 23  
`get_coefficient_lev2V()` (pyCDB.client.CDBClient method), 23  
`get_coefficient_V2unit()` (pyCDB.client.CDBClient method), 23

[get\\_computer\\_id\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_conf\\_value\(\)](#) (pyCDB.pyCDBBase.pyCDBBase class method), [35](#)  
[get\\_data\\_file\\_reference\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_data\\_root\\_path\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_data\\_source\\_id\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_data\\_source\\_references\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_data\\_source\\_subdir\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_file\\_log\\_level\(\)](#) (pyCDB.pyCDBBase.pyCDBBase method), [35](#)  
[get\\_FS\\_attachement\(\)](#) (pyCDB.DAQClient.CDBDAQClient method), [32](#)  
[get\\_FS\\_signal\\_reference\(\)](#) (pyCDB.client.CDBClient method), [22](#)  
[get\\_generic\\_signal\\_id\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_generic\\_signal\\_references\(\)](#) (pyCDB.client.CDBClient method), [23](#)  
[get\\_h5\\_dataset\(\)](#) (in module pyCDB.client), [31](#)  
[get\\_h5\\_dataset\\_shape\(\)](#) (in module pyCDB.client), [31](#)  
[get\\_log\\_level\(\)](#) (pyCDB.client.CDBSignal method), [29](#)  
[get\\_log\\_level\(\)](#) (pyCDB.pyCDBBase.pyCDBBase method), [35](#)  
[get\\_max\\_revisions\(\)](#) (pyCDB.client.CDBClient class method), [24](#)  
[get\\_record\\_datetime\(\)](#) (pyCDB.client.CDBClient method), [24](#)  
[get\\_record\\_number\\_from\\_FS\\_event\(\)](#) (pyCDB.DAQClient.CDBDAQClient method), [33](#)  
[get\\_record\\_path\(\)](#) (pyCDB.client.CDBClient method), [24](#)  
[get\\_shot\\_comments\(\)](#) (pyCDB.logbook.logbook method), [36](#)  
[get\\_shot\\_info\(\)](#) (pyCDB.logbook.logbook method), [36](#)  
[get\\_shot\\_linear\\_profile\\_value\(\)](#) (pyCDB.logbook.logbook method), [36](#)  
[get\\_shot\\_linear\\_profiles\(\)](#) (pyCDB.logbook.logbook method), [36](#)  
[get\\_shot\\_param\\_value\(\)](#) (pyCDB.logbook.logbook method), [37](#)  
[get\\_shot\\_params\(\)](#) (pyCDB.logbook.logbook method), [37](#)  
[get\\_signal\(\)](#) (cdb\_client method), [39](#)  
[get\\_signal\(\)](#) (pyCDB.client.CDBClient method), [24](#)  
[get\\_signal\\_base\\_tree\(\)](#) (pyCDB.client.CDBClient method), [24](#)  
[get\\_signal\\_calibration\(\)](#) (pyCDB.client.CDBClient method), [24](#)  
[get\\_signal\\_data\(\)](#) (pyCDB.client.CDBClient method), [24](#)  
[get\\_signal\\_data\\_tree\(\)](#) (pyCDB.client.CDBClient method), [25](#)  
[get\\_signal\\_parameters\(\)](#) (pyCDB.client.CDBClient method), [25](#)  
[get\\_signal\\_references\(\)](#) (pyCDB.client.CDBClient method), [25](#)  
[get\\_signal\\_setup\(\)](#) (pyCDB.client.CDBClient method), [25](#)  
[get\\_table\\_struct\(\)](#) (in module pyCDB.CodeGeneration), [37](#)  
[get\\_table\\_struct\(\)](#) (pyCDB.pyCDBBase.pyCDBBase method), [35](#)  
[get\\_tags\(\)](#) (pyCDB.logbook.logbook method), [37](#)  
[get\\_unit\(\)](#) (pyCDB.client.CDBClient method), [25](#)  
[get\\_unit\\_factor\\_tree\(\)](#) (pyCDB.client.CDBClient method), [26](#)  
[get\\_unit\\_system\(\)](#) (pyCDB.client.CDBClient method), [26](#)  
[getAlias\(\)](#) (Java method), [60](#)  
[getAttachedGenericSignalId\(\)](#) (Java method), [54](#)  
[getAttachment\(long\)](#) (Java method), [49](#)  
[getAttachmentTable\(long, long, long\)](#) (Java method), [50](#)  
[getAttachmentTable\(ParameterList\)](#) (Java method), [49](#)  
[getAttribute\(String\)](#) (Java method), [81](#)  
[getAxes\(\)](#) (Java method), [63](#)  
[getAxis\(int\)](#) (Java method), [60](#), [63](#)  
[getAxisId\(int\)](#) (Java method), [60](#)  
[getCoefficient\(\)](#) (Java method), [66](#)  
[getCoefficientLev2V\(\)](#) (Java method), [55](#)  
[getCoefficientV2Unit\(\)](#) (Java method), [55](#)  
[getCollectionName\(\)](#) (Java method), [55](#)  
[getComputerId\(String\)](#) (Java method), [50](#)  
[getComputerId\(String, String\)](#) (Java method), [50](#)  
[getDAQParameters\(\)](#) (Java method), [67](#)  
[getDAQParametersString\(\)](#) (Java method), [67](#)  
[getData\(\)](#) (Java method), [63](#)  
[getDataFile\(\)](#) (Java method), [57](#), [63](#)  
[getDataFile\(long\)](#) (Java method), [50](#)  
[getDataFile\(ParameterList\)](#) (Java method), [50](#)  
[getDataFileId\(\)](#) (Java method), [63](#)  
[getDataSourceId\(\)](#) (Java method), [61](#)  
[getDataSourceId\(String\)](#) (Java method), [50](#)  
[getDescription\(\)](#) (Java method), [61](#), [63](#)  
[getDictKeys\(\)](#) (Java method), [77](#)  
[getDictValue\(String\)](#) (Java method), [77](#)  
[getDictValueAsDouble\(String\)](#) (Java method), [77](#)  
[getDictValueAsLong\(String\)](#) (Java method), [77](#)  
[getDictValueAsString\(String\)](#) (Java method), [77](#)  
[getDimensions\(\)](#) (Java method), [75](#)  
[getFileKey\(\)](#) (Java method), [57](#), [59](#)  
[getFileName\(\)](#) (Java method), [55](#)  
[getFilePath\(\)](#) (Java method), [59](#)

getFSSignal(String, String, String) (Java method), 50  
getFullPath() (Java method), 55  
getGenericSignal() (Java method), 57, 63  
getGenericSignal(long) (Java method), 50  
getGenericSignal(String) (Java method), 50  
getGenericSignalId() (Java method), 64  
getGenericSignalParameters() (Java method), 67  
getGenericSignalParametersString() (Java method), 67  
getGenericSignalReference() (Java method), 64  
getGenericSignals(ParameterList) (Java method), 50  
getGenericSignals(String) (Java method), 51  
getId() (Java method), 55, 61  
getInstance() (Java method), 51  
getInstance(ParameterList) (Java method), 51  
getInterpreter() (Java method), 82  
getLastRecordNumber() (Java method), 61  
getName() (Java method), 61, 64  
getOffset() (Java method), 66  
getParameters() (Java method), 55, 64  
getParametersString() (Java method), 55  
getPyObject() (Java method), 64, 81  
getRank() (Java method), 75  
getRawData() (Java method), 75  
getRecordNumber() (Java method), 58, 64  
getRecordNumberFromFSEvent(long, String) (Java method), 51  
getSignal(long, long) (Java method), 51  
getSignal(ParameterList) (Java method), 51  
getSignal(String) (Java method), 51  
getSignal(String, String) (Java method), 51  
getSignalCalibration(boolean) (Java method), 64  
getSignalCalibration(ParameterList) (Java method), 52  
getSignalCalibration(String) (Java method), 52  
getSignalParameters(ParameterList) (Java method), 52  
getSignalReference() (Java method), 64  
getSignalReference(long, long) (Java method), 52  
getSignalReference(long, long, long) (Java method), 52  
getSignalReference(String) (Java method), 52  
getSignalReferences(ParameterList) (Java method), 52  
getSignalType() (Java method), 61  
getTimeAxis() (Java method), 61, 64  
getTimeAxisId() (Java method), 61  
getType() (Java method), 76  
getUnit() (Java method), 65  
getUnitFactor() (Java method), 65  
getUnits() (Java method), 61, 66  
gson (Java field), 56

## H

has\_tag() (pyCDB.logbook.logbook method), 37  
hasAttribute(String) (Java method), 81  
hasUnitFactor() (Java method), 65  
Hdf5ReadException (Java class), 77  
Hdf5ReadException(String) (Java constructor), 78

Hdf5Utils (Java class), 78  
Hdf5WriteException (Java class), 79  
Hdf5WriteException(String) (Java constructor), 79

## I

info() (pyCDB.client.CDBSignal method), 29  
insert() (pyCDB.pyCDBBase.pyCDBBase method), 35  
insert(String, ParameterList) (Java method), 52  
insert(String, ParameterList, Boolean) (Java method), 52  
insert(String, ParameterList, Boolean, Boolean) (Java method), 52  
instanceTest() (Java method), 68  
INT16 (Java field), 72  
INT32 (Java field), 72  
INT64 (Java field), 72  
INT8 (Java field), 72  
invoke(PyObject, String, Map) (Java method), 82  
invoke(String) (Java method), 82  
invoke(String, Map) (Java method), 81  
is\_channel\_attached() (py-CDB.DAQClient.CDBDAQClient method), 33  
is\_file\_ready() (pyCDB.client.CDBClient method), 26  
is\_generic\_signal\_attached() (py-CDB.DAQClient.CDBDAQClient method), 33  
is\_json() (pyCDB.pyCDBBase.pyCDBBase class method), 35  
isDAV() (Java method), 65  
isFile() (Java method), 61, 65  
isintlike() (in module pyCDB.pyCDBBase), 34  
isLinear() (Java method), 62, 65  
isRaw() (Java method), 65  
isstringlike() (in module pyCDB.pyCDBBase), 34  
isValidType(String) (Java method), 62

## J

JsonUtils (Java class), 79  
JyCDBTest (Java class), 68

## L

last\_record\_number() (pyCDB.client.CDBClient method), 26  
last\_shot\_number() (pyCDB.client.CDBClient method), 26  
lastRecordNumber(String) (Java method), 52  
lastShotNumber() (Java method), 53  
LINEAR (Java field), 60  
linearTransform(double, double) (Java method), 76  
linearTransform(double[], double, double) (Java method), 70  
linearTransform(Object, double, double) (Java method), 70



load\_h5() (pyCDB.pyCDBBase.OrderedDict class method), 34  
 log\_level (pyCDB.client.CDBSignal attribute), 29  
 log\_level (pyCDB.pyCDBBase.pyCDBBase attribute), 35  
 logbook (class in pyCDB.logbook), 36  
 logWithDuration(String) (Java method), 76  
 logWithDuration(String, Logger) (Java method), 76

## M

main(String[]) (Java method), 68  
 mkdirs() (pyCDB.pyCDBBase.pyCDBBase class method), 35  
 max(double[]) (Java method), 70  
 MAX\_AXES (Java field), 60  
 MAX\_RANK (Java field), 72  
 min(double[]) (Java method), 71  
 module  
   \_\_main\_\_, 41  
   search path, 41  
   sys, 41

multiply(double[], double) (Java method), 71  
 mysql\_str() (pyCDB.pyCDBBase.pyCDBBase class method), 35  
 mysql\_values() (pyCDB.pyCDBBase.pyCDBBase class method), 35

## N

new\_data\_file() (pyCDB.client.CDBClient method), 26  
 newDataFile(ParameterList) (Java method), 53  
 newDataFile(String, long, long) (Java method), 53  
 newObject(String, Map) (Java method), 82  
 nodiff\_signal\_data() (py-CDB.DAQClient.CDBDAQClient method), 33

## O

OrderedDict (class in pyCDB.pyCDBBase), 33

## P

ParameterList (Java class), 80  
 parseNative(String) (Java method), 79  
 path  
   module search, 41  
 plot() (pyCDB.client.CDBSignal method), 29  
 product(int[]) (Java method), 71  
 product(long[]) (Java method), 71  
 put(String, boolean) (Java method), 80  
 put(String, double) (Java method), 80  
 put(String, int) (Java method), 80  
 put(String, java.sql.Date) (Java method), 80  
 put(String, long) (Java method), 80  
 put(String, ParameterList) (Java method), 80  
 put(String, short) (Java method), 80

put(String, String) (Java method), 80  
 put(String, Timestamp) (Java method), 80  
 put\_signal() (pyCDB.client.CDBClient method), 26  
 pyCDB.client (module), 21  
 pyCDB.CodeGeneration (module), 37  
 pyCDB.DAQClient (module), 31  
 pyCDB.logbook (module), 36  
 pyCDB.pyCDBBase (module), 33  
 pyCDBBase (class in pyCDB.pyCDBBase), 34  
 PyObject (Java field), 81  
 PythonAdapter (Java class), 81  
 PythonAdapter(PyObject) (Java constructor), 81  
 PythonUtils (Java class), 82

## Q

query() (pyCDB.pyCDBBase.pyCDBBase method), 36

## R

range(int, int) (Java method), 71  
 rank(Object) (Java method), 71  
 RAW (Java field), 66  
 read\_file\_data() (in module pyCDB.client), 31  
 readFileData(String, String) (Java method), 78  
 readSignalInfo() (Java method), 59  
 record\_exists() (pyCDB.client.CDBClient method), 27  
 recordExists(long) (Java method), 53  
 reshape(Object, int[]) (Java method), 71  
 resolve\_record\_number() (pyCDB.client.CDBClient method), 27  
 row\_as\_dict() (pyCDB.pyCDBBase.pyCDBBase class method), 36

## S

save\_h5() (pyCDB.pyCDBBase.OrderedDict method), 34  
   search  
     path, module, 41  
 set(Object, int[], Object) (Java method), 71  
 set\_file\_log\_level() (pyCDB.pyCDBBase.pyCDBBase method), 36  
 set\_file\_ready() (pyCDB.client.CDBClient method), 27  
 set\_log\_level() (pyCDB.client.CDBSignal method), 30  
 set\_log\_level() (pyCDB.pyCDBBase.pyCDBBase method), 36  
 setFileReady() (Java method), 58  
 setFileReady(long) (Java method), 53  
 setRequireChannelAttachment(boolean) (Java method), 58  
 shot\_param\_name() (pyCDB.logbook.logbook method), 37  
 Signal (Java class), 62  
 Signal(Signal, int, PythonAdapter, DictionaryAdapter) (Java constructor), 62  
 Signal(Signal, int, PythonAdapter, String) (Java constructor), 62

Signal(SignalReference, String) (Java constructor), [62](#)  
 signal\_dim() (pyCDB.client.CDBClient class method), [27](#)  
 signal\_setup() (pyCDB.DAQClient.CDBDAQClient method), [33](#)  
 SignalCalibration (Java class), [65](#)  
 signalExists() (Java method), [59](#)  
 SignalParameters (Java class), [66](#)  
 size(Object) (Java method), [72](#)  
 sql\_to\_python\_type() (pyCDB.pyCDBBase.pyCDBBase class method), [36](#)  
 store\_channel\_data\_as\_signal() (pyCDB.client.CDBClient method), [27](#)  
 store\_signal() (pyCDB.client.CDBClient method), [28](#)  
 storeLinearSignal(long, long, double, double, double) (Java method), [53](#)  
 storeSignal(long, long, String, long, double) (Java method), [53](#)  
 storeSignal(long, long, String, long, double, double, double, double) (Java method), [53](#)  
 storeSignal(ParameterList) (Java method), [53](#)  
 storeSignalAndFile() (Java method), [59](#)  
 switch\_channels() (pyCDB.DAQClient.CDBDAQClient method), [33](#)  
 sys module, [41](#)

## U

uid\_name() (pyCDB.logbook.logbook method), [37](#)  
 UnknownDataTypeException (Java class), [82](#)  
 UnknownDataTypeException(String) (Java constructor), [83](#)  
 update() (pyCDB.pyCDBBase.pyCDBBase method), [36](#)  
 update\_signal() (pyCDB.client.CDBClient method), [29](#)  
 updateSignal(ParameterList) (Java method), [54](#)  
 url (pyCDB.client.CDBSignal attribute), [30](#)

## V

VALID\_TYPES (Java field), [60](#)  
 validate\_variant() (in module pyCDB.client), [31](#)

## W

writeAxis(int, double, double) (Java method), [58](#)  
 writeData(Data, long[]) (Java method), [58](#)  
 writeData(String, String, Data) (Java method), [78](#)  
 writeData(String, String, Data, long[]) (Java method), [78](#)  
 writeData(String, String, int, long[], Object, long[]) (Java method), [79](#)  
 writeNotAttachedSignal(double, double) (Java method), [58](#)  
 writeSignal(double) (Java method), [58](#)  
 WrongDimensionsException (Java class), [83](#)  
 WrongDimensionsException(String) (Java constructor), [83](#)