
AutoAnnotate Documentation

Release 1.1

Mike Kucera

Jan 18, 2018

1	Getting Started (Installing)	1
1.1	Requirements	1
1.2	Installation	1
2	Creating an Annotation Set	3
2.1	Create Annotation Set Dialog	3
2.2	ClusterMaker2 Algorithms	5
3	Working With AutoAnnotate	7
3.1	Panels	7
3.2	Changing Display Options	8
3.3	Managing Annotation Sets	9
3.4	Managing Clusters in the Side Panel	11
3.5	Managing Clusters in the Network View	12
3.6	Collapsing and Expanding Clusters	12
3.7	Summary Network	14
3.8	Working with WordCloud	14
3.9	Saving and Restoring a Session with Annotation Sets	15
3.10	Additional Options	15
4	Bug Reports and Enhancement Requests	17

Getting Started (Installing)

1.1 Requirements

- Cytoscape minimum version 3.3
- AutoAnnotate version 1.1
- WordCloud minimum version 3.1.0
- clusterMaker2 minimum version 0.9.5

1.2 Installation

Install Cytoscape

- If you don't have Cytoscape please download and install the latest release from <http://www.cytoscape.org/download.php>.

Install AutoAnnotate (latest version is 1.1)

- Open Cytoscape
- In the main menu select Apps > App Manager
- In the App Manager select AutoAnnotate in the list of All Apps and click the Install button.
- The WordCloud app will also be installed automatically along with AutoAnnotate

Install clusterMaker2 (optional, but highly recommended)


- In the App Manager select clusterMaker2 in the list of All Apps and click the Install button.

Alternatively AutoAnnotate and clusterMaker2 can be installed from the App Store website.

- AutoAnnotate on the App Store: <http://apps.cytoscape.org/apps/autoannotate>
- clusterMaker2 on the App Store: <http://apps.cytoscape.org/apps/clustermaker2>

Creating an Annotation Set

- Start by selecting a network view.
- In the main menu select **Apps > AutoAnnotate > New Annotation Set...**
- The Create Annotation Set dialog will open.




2.1 Create Annotation Set Dialog

This dialog is used to select the options for creating an Annotation Set.

Note: AutoAnnotate can detect if the network was created by EnrichmentMap, and if so sensible defaults for EnrichmentMap networks will be pre-selected.

Note: You may manually edit the contents and labels of the clusters after the Annotation Set is created, and you may create as many Annotation Sets as you like, so its not critical to get everything perfect at this step.



- Cluster Options
 - Use clusterMaker App
 - * If this is selected then clusterMaker will be used to calculate the clusters.
 - Cluster algorithm:
 - * Used to select the clustering algorithm provided by clusterMaker that will be used to calculate the clusters.
 - * See the table below for a list of clusterMaker algorithms currently available for use in AutoAnnotate.
 - Edge weight column
 - * Some clusterMaker algorithms use edge weights. This box is used to select the Edge Column to use for the edge weights. If the selected algorithm does not use edge weights then this box will be disabled (greyed out).
 - User-defined clusters
 - * You may provide your own clusters or use another clustering app to calculate the clusters. In this case the clusters must already be calculated and cluster identifiers must be stored in a Node Column.
 - * Any Node Column of may be used to identify the clusters.
 - * If you select a column of List type then the clusters may overlap.
- Label Options
 - Label Column
 - * Select a Node Column of type String that will be used to calculate the cluster labels.
 - Label Algorithm

- * Select an algorithm that is used to calculate the label for each cluster.

- * Two label algorithms are currently available:

1. WordCloud: Adjacent Words

Uses the WordCloud app to calculate the labels.

Based on the algorithm that is used by the “Cloud Style: Clustered Standard” option in WordCloud. Words in the label are chosen from the most frequent words in the cluster, and words that are adjacent to the most frequent words. The higher the “*adjacent word bonus*” is, the more likely adjacent words will be in the label.

2. WordCloud: Largest Words

Uses the WordCloud app to calculate the labels.

Basically the same algorithm that is used by the “Cloud Style: No Clustering” option in WordCloud. Words in the label are chosen from the most frequent words in the cluster.

- Click **Create Annotations** to create the new Annotation Set.

WARNING: AutoAnnotate will manage all groups (compound nodes) for the current network view. It is not recommended to manually create groups on a network view that has Annotation Sets. If you would like to manually create groups at any time please duplicate the network view.

2.2 ClusterMaker2 Algorithms

The following clusterMaker2 algorithms are made available in AutoAnnotate:

Algorithm	Weighted?
Affinity Propagation Cluster	Yes
Cluster Fuzzifier	Yes
Community cluster (GLay)	No
ConnectedComponents Cluster	Yes
Fuzzy C-Means Cluster	Yes
MCL Cluster	Yes
SCPS Cluster	Yes


Weighted algorithms require an edge attribute to be used as the edge weight.

Working With AutoAnnotate

3.1 Panels

When an Annotation Set is created the AutoAnnotate panels will be shown.


Note: To show or hide the AutoAnnotate panels at any time go to the main menu and select **Apps > AutoAnnotate > (Show/Hide) AutoAnnotate**.



- In the **Control Panel** on the left side is the main AutoAnnotate panel.
 - This panel is used to manage Annotation Sets and Clusters.
- In the **Results Panel** on the right side is the Display Options panel.
 - This panel is used to modify the look of the annotations.

3.2 Changing Display Options

This AutoAnnotate Display panel shows display options for the currently selected Annotation Set.



Border width Slider that adjusts the width of the cluster annotation borders.

Opacity Slider that adjusts the opacity of the cluster annotation fill color. Slide to the left for completely transparent. Slide to the right for completely opaque.

Font Size

Scale font by cluster size When selected the size of the label will be relative to the number of nodes in the cluster. Larger clusters will have larger labels. When deselected all the clusters will have the same size labels. Slider is used to adjust the font size or font scale.


Shape There are two available shapes: Ellipse or Rectangle

Hide Clusters Select to hide all cluster shape annotations.


Hide Labels Select to hide all cluster labels.

3.3 Managing Annotation Sets

The main AutoAnnotate panel is used to manage Annotation Sets.



- At the top of the panel is a dropdown box that allows you to switch between all of the Annotation Sets that have been created for the current network view.
 - Select the (none) option to clear all annotations from the view.
 - Select an Annotation Set to show all the annotations and labels in that set.
 - When switching between Annotation Sets all the annotations in the view will be erased and redrawn.
- * Similarly when switching Annotation Sets any collapsed groups will be expanded.



Next to the dropdown box is the menu button. Click this button to show the Annotation Set Menu.

New Annotation Set... Shows the Create Annotation Set dialog.

Rename Renames the current annotation set

Delete... Deletes the current annotation set

Manage Annotation Sets... Opens a Dialog where you can change the order that Annotation Sets appear in the dropdown and delete multiple Annotation Sets at a time.

Collapse All Collapses all the clusters *Warning:* Collapsing all clusters in a large network can be very slow. Consider using the “Create Summary Network” option instead.

Expand all Expands all the clusters

Create Summary Network Generates a new network where each node represents a collapsed cluster in the original network. Very similar in functionality to “Collapse All”, but runs significantly faster.

Layout Clusters Runs a layout algorithm that will layout each cluster separately using the perfuse layout (evenly spaced in a grid). Note: this operation cannot be undone.

Redraw annotations Older versions of Cytoscape have a bug that cause annotations to occasionally be drawn in the wrong location. Redrawing the annotations can fix this problem.

Recalculate Labels Recalculates the labels for all of the clusters in the annotation set. Typically used when the WordCloud options or label algorithm changes. More on this later in this guide. (Note: to recalculate labels for a subset of clusters please select the clusters in the table and right-click to show the cluster context menu.)

Label Options Opens a dialog where you can select the label algorithm that is used when recalculating labels.

Settings... Opens the AutoAnnotate settings dialog.


3.4 Managing Clusters in the Side Panel

On the main AutoAnnotate panel you will find the Cluster Table.

- This table shows each cluster, along with its label, the number of nodes in the cluster and a checkbox icon that indicates if the cluster has been collapsed.
- Click on a table header to sort the table by that column.

Most of the cluster options are available using the context menu of the cluster table.

- Start by selecting one or more clusters in the cluster table.
- To select multiple clusters hold the Shift key to select a range or hold the Ctrl/Command key to select individual clusters.



Right click with the mouse to show the cluster table context menu.

Rename... Enabled when 1 cluster is selected. Manually set the label for the cluster.

Delete Delete the cluster annotations Note: This does not delete the nodes or edges in the cluster, just the annotations.

Merge Enabled when 2 or more clusters are selected. Merges the clusters into a single cluster. A new label for the merged cluster is automatically calculated.

Extract Clusters... Creates a new Annotation Set consisting of the selected clusters.


Collapse Collapses the selected clusters **Warning:** can be very slow when many clusters are selected.

Expand Expands the selected clusters

Generate Summary Network Generates a new network where each node represents a collapsed cluster that is selected in the original network. Very similar in functionality to “Collapse”, but runs significantly faster.

Recalculate Labels Recalculates the labels for the selected clusters.

3.5 Managing Clusters in the Network View




- Manually create a cluster
 - Select some nodes in the network view, right click to show the context menu, and select **Apps > AutoAnnotate > Create Cluster**.
 - * A new cluster is created.
 - * A label for the cluster is automatically calculated.
 - * The cluster is selected in the cluster table.
- Deleting a node in an existing cluster will remove the node from the cluster.
 - A new label is calculated.

3.6 Collapsing and Expanding Clusters

AutoAnnotate can summarize a network by collapsing some or all of the clusters automatically.

Warning: As of Cytoscape 3.4 this function can be very slow, depending on the number of clusters. Consider using the “Generate Summary Network” option instead because it runs significantly faster. See the section below this one for more details.



The Cytoscape “group nodes” feature is used to collapse and expand clusters. A group node is a compound node that contains other nodes.

When AutoAnnotate collapses a cluster it first creates a group node that contains all the nodes in the cluster and then the group node is collapsed. When the cluster is collapsed Cytoscape will create “meta-edges” between the group node and any other nodes it is connected to. The collapsed group nodes and the meta-edges provide a summary of the network.


The annotation for a cluster is not shown when the cluster is collapsed. The label shown for a collapsed cluster is generated by Cytoscape and depends on the Group Preferences (explained below).

When a cluster is expanded the group node is expanded and then deleted. The annotations are shown again.

It is not recommended to manually collapse or expand clusters using the Groups menu as this may confuse AutoAnnotate. When switching between Annotation Sets all collapsed clusters are automatically expanded.

Clusters can be collapsed or expanded from the Annotation Set Menu or the Cluster Table Context Menu. Before collapsing clusters make sure that node aggregation is enabled. This will tell Cytoscape to generate a label for the group node that aggregates the labels of the nodes in the cluster.

- In the main menu select **Edit > Preferences > Group Preferences**.
- In the Group Preferences dialog select the option for Enable Attribute Aggregation.



3.7 Summary Network

Using group nodes to collapse clusters can be very slow when there are many clusters. As an alternative a new option called **Generate Summary Network** has been added in AutoAnnotate 1.1. A new network is generated where each node represents a collapsed cluster in the original network, and each edge represents a meta-edge between clusters. The resulting network is very similar to the results you get from collapsing the clusters.

The main differences are:

- Nodes that are not contained in any cluster are ignored.
- There is no relationship maintained back to the original network. If the clusters in the original network change the summary network will not be updated. Create a new summary network at any time to reflect the current state of the clusters.
- A new visual style for the summary network is created based on the original visual style.

A summary network can be created from the Annotation Set Menu or the Cluster Table Context Menu.

Summary network uses the same attribute aggregation settings as collapsed group nodes. Available from **Edit > Preferences > Group Preferences** then expand the Default Aggregation Settings section.

3.8 Working with WordCloud

Labels for clusters are generated by the WordCloud app.


- You may change the label generation options from the WordCloud panel.
- In the main menu select **Apps > WordCloud > Show WordCloud**.
- In the WordCloud panel you can change the label options, for example using the Normalize slider.
- Back in the AutoAnnotate panel select **Recalculate Labels** to regenerate labels using the updated WordCloud options.

3.9 Saving and Restoring a Session with Annotation Sets

There are currently some limitations in Cytoscape with regards to saving and restoring annotations and groups in a session file. For this reason when you save your session all the clusters must be expanded first. This is planned to be fixed in future versions of Cytoscape and AutoAnnotate.

3.10 Additional Options

AutoAnnotate has various warning dialogs that can be dismissed permanently by selecting “don’t show this again”. If you would like to see these warnings again they can be re-enabled. In the main menu select **Edit > Preferences > Properties**, then select **org.baderlab.autoannotate** in the Properties Dialog.



Bug Reports and Enhancement Requests

Thank you for using AutoAnnotate. Your feedback is important and helps us to improve the App for everyone. If you have any comments, requests for new features or enhancements, or bugs to report, please enter them on the AutoAnnotate github issue tracker.

<https://github.com/BaderLab/AutoAnnotateApp/issues>